

Invisishield™ Drapes

Unique products to help improve safety and patient care



**Protect and
see clearly**

Overview

Medline's Invisishield specialty drapes are designed to make the clinician's job easier, while improving safety for both patients and staff.

Invisishield patient drapes

Medline's Invisishield drapes are clear so that clinicians can easily see theatre equipment and the incision area. The impervious, tear-resistant material provides exceptional protection and offers superb drapeability.

Clear, impervious and drapeable

Configurations for all major procedures and applications can be found in each section of our different Invisishield patient drapes.

- Transparent – Clinicians can easily view the area around incision sites.
- Impervious – Made of tear-resistant polyethylene (PE) and polyurethane films, Invisishield drapes provide superb protection from fluids.
- Drapeable – Invisishield drapes have just the right amount of flexibility and softness for superb drapeability.

Useful features

Invisishield drapes have handy features that help make the clinician's job easier. Our fluid pouches feature handy tabs on the adhesive liner that make removal easy. Capped, closable ports eliminate the need to cut or pierce the pouch for drainage, improving fluid control and reducing set-up times.

Table of contents

General accessories	3
Orthopaedic drapes	8
Ophthalmic drapes	10
Other drapes	13
Dermadine Plus and Dermadrape Incise Drapes	15
Equipment drapes	16
Stockinettes	26
Cross reference	27

General accessories

Towel drapes

Description

- Strong adhesive with no residue
- Non-linting, tear-resistant, release liner
- Strong, drapeable, low-glare matte-finish film

Use

Prepping, tourniquet and instrument cover, squaring, isolating sterile from non-sterile areas, fluid barrier, non-linting surface

Item number	Dimensions	Packaging
IDE1000	45 x 30 cm adhesive only: 45 x 5.1 cm	10 ea/bx; 40 ea/cs
IDE1010	45 x 61 cm adhesive only: 45 x 5.1 cm	10 ea/bx; 40 ea/cs



Easy-peel liner:
less force is needed to separate it
from the adhesive section.
The reinforced adhesive section
reduces the chance of wrinkling.

Isolation bags

Description

- Soft, strong, clear film with seamless bottom
- Cotton drawstrings for easy closure
- Cuffed opening to improve sterile delivery

Use

Isolate foot, hand, non-linting surface

Item number	Dimensions	Packaging
ICE1003	50 x 50 cm	10 ea/bx; 40 ea/cs
586560-1	60 x 65 cm	30 ea/bx; 180 ea/cs




Fluid pouches with drainage ports

Description

- Wide adhesive area for better adhesion
- Strong adhesive with no residue
- Non-linting, tear-resistant, paper release liner
- Mouldable strip to enhance pouch performance
- Micro screen for improved filtering system
- Suction port with cap for improved fluid control
- Heavy, low-glare matte-finish film

Use

CV, orthopaedic, neurology, urology, general



Extra-long, mouldable strip allows for customisation.

Item number	Dimensions	Packaging
IDE1016	50 x 50 cm Adhesive only: 50 x 15 cm	40 ea/cs
4638CE	61 x 51 cm Adhesive only: 51 x 13 cm	40 ea/cs
DYJPEFP020	51 x 56 cm Adhesive only: 46 x 5 cm	50 ea/cs


Instrument pouches

Description

- 2 compartments for cauteries or catheters
- New high-tack, aggressive adhesive for better adhesion to linen or non-woven materials
- Heavy, low-glare matte-finish film

Use

CV, orthopaedic, neurology



Item number	Dimensions	Packaging
IDE1018	17 x 29 cm Adhesive only (2): 17 x 3.8 cm	10 ea/bx; 40 ea/cs


Long instrument pouches

Description

- 3 long compartments for cauteries or catheters
- New high-tack, aggressive adhesive for better adhesion to linen or non-woven materials
- Heavy, low-glare matte-finish film

Use

CV, laparoscopy, gynaecology



Item number	Dimensions	Packaging
IDE1018L	25 x 47 cm Adhesive only (2): 25 x 5 cm	10 ea/bx; 40 ea/cs

Instrument pouches

Description




- 1 large compartment for instruments
- New high-tack adhesive tape
- Clear, strong polyethylene film

Use

CV, orthopaedic, neurology, general


Item number	Dimensions	Packaging
DYJPEISTPOM1	33 x 41 cm Adhesive only: 41 x 5 cm	100 ea/cs
DYJPEISTPOM2	33 x 45 cm Adhesive only: 41 x 5 cm	100 ea/cs
DYJPEISTPOM3	43 x 48 cm Adhesive (2): 43 x 5 cm	60 ea/cs
04A040CE	33 x 51 cm Adhesive only: 51 x 5 cm	180 ea/cs
96860CE	30.5 x 29 cm Adhesive (2): 30.5 x 5 cm	150 ea/cs
96890CE	43 x 49 cm Adhesive (2): 49 x 5 cm	250 ea/cs
9686CE*	27 x 29 cm Adhesive (2): 29 x 5 cm	150 ea/cs
9687CE*	36 x 86 cm Adhesive (2): 36 x 5 cm	120 ea/cs
9689CE*	43.5 x 49 cm Adhesive (2): 49 x 5 cm	150 ea/cs

*Made of non-transparent spunlace material



Set-up covers

Item number	Dimensions	Packaging
ICE2158	46 x 61 cm, adhesive	40 ea/cs
ICE2008	137 x 229 cm	20 ea/cs



Sterile sponge counter bags

Description

- Clear, drapeable polyethylene film
- 1 adhesive strip for attachment
- 3 sizes
- 3 compartment options



Sterile sponge counter bags allow for safe, effective and convenient counting, handling and disposal of contaminated swabs. They minimise the risk of pathogen inhalation and keep the surgical area clean and organised.

Item number	Dimensions	Packaging
IBE5242	52 x 42 cm, 40 compartments	50 ea/cs
IBE5637	56 x 37 cm, 6 compartments	50 ea/cs
IBE5652	56 x 52 cm, 20 compartments	50 ea/cs

Non-sterile sponge counter bags

Description

- Clear, drapeable polyethylene film
- Available with both a transparent and blue backing
- Can be hung from sponge counter rack NON50511, which attaches to most IV poles



Item number	Dimensions	Packaging
DYNJE1000	25 x 94 cm, 10 compartments	50 ea/cs
DYNJE1001	25 x 94 cm, 10 compartments	50 ea/cs
NON50511	Can hold up to 4 sets of bags	50 ea/cs

Non-sterile sponge counter bags can be purchased in blue or transparent for different customer preferences. The sponge counter bags are hung on a sponge counter rack in the operating theatre and are used after surgical procedures to quickly count lap sponges.

Orthopaedic drapes

Polyethylene U-drape



Description

- 'U' fenestration with wide adhesive area
- Advanced adhesive, providing excellent adhesion without residue
- Strong, drapeable, low-glare matte-finish film
- Split, paper liner for ease of application

Use

Prepping procedures: arthroscopies, hip, knee, shoulder, back, cranial, isolation and applications

Item number	Dimensions	Packaging
IDE1015	121 x 132 cm Adhesive only: 25 x 48 cm	10 ea/bx; 40 ea/cs
IDE1019	121 x 180 cm Adhesive only: 25 x 48 cm 'U' fenestration only: 15 x 48 cm	10 ea/bx; 40 ea/cs



These drapes include courtesy tabs for easy removal.

Shoulder pouches


Description

- Oval aperture with slit and adhesive areas
- Strong adhesive with no residue
- Heavy, low-glare matte-finish film for additional strength
- 2 mouldable strips and 1 foam strip to enhance pouch performance
- 2 suction ports with caps for improved fluid control

Use

Shoulder arthroscopy, implant applications

Item number	Dimensions	Packaging
IDE1067	91 x 73 cm Adhesive only: 25 x 48 cm Aperture: 19 x 25 cm	24 ea/cs



The capped, closable port eliminates the step of cutting or piercing for drainage and makes it easy to stop the flow at the end of the procedure.


U fluid collection pouches

Description

- 'U' aperture with adhesive area
- Strong adhesive with no residue
- 3 mouldable strips to enhance pouch performance
- Strong, drapeable, clear polyethylene film

Use

Shoulder arthroscopy, hip surgery, orthopaedic



Item number	Dimensions	Packaging
DYJPEFP090	76 x 112 cm Adhesive only: 25 x 51 cm Aperture: 15 x 46 cm	25 ea/cs


Patient isolation drapes

Description

- Less aggressive adhesive strip for easy removal
- Aggressive incise adhesive with no residue
- Clear, strong polyethylene film

Use

Hip pinning



Item number	Dimensions	Packaging
IDE1014	317 x 213 cm Adhesive only: 49 x 25 cm Adhesive strip only: 317 x 1.3 cm	5 ea/bx; 15 ea/cs

Patient isolation drapes with pouch


Description

- Strong incise adhesive with no residue
- Clear, strong polyethylene film
- 2 cord holders and instrument pouches for better organisation
- Large pouch with micro screen, suction port and cap for improved fluid control

Use

Hip pinning, nephroscopies, biliaries

Item number	Dimensions	Packaging
IDE1017	317 x 213 cm Adhesive only: 49 x 25 cm Adhesive strip only: 317 x 1.3 cm Pouch only: 51 x 51 cm	15 ea/cs
IDE1017W	317 x 254 cm Adhesive only: 85.1 x 27.9 cm Adhesive strip only: 317 x 1.3 cm Pouch only: 69.9 x 94 cm	15 ea/cs
5801CE	330 x 270 cm Adhesive only: 78 x 26 cm Adhesive strip only: 330 x 1.9 cm Pouch only: 25.5 x 28 cm	16 ea/cs



These drapes have a superb filter with deep placement for quick drainage and numerous small perforations.

Ophthalmic drapes

Aperture drapes



Description

- Circular aperture with adhesive area
- Non-linting, tear-resistant, paper release liner
- Thin, soft, drapeable, low-glare matte-finish film

Use

Cataract, corneal, refractive, lid/plastic, minor

Item number	Dimensions	Packaging
IDE1020	40 x 40 cm Adhesive only: 12 x 12 cm Aperture only: 6 cm circle	10 ea/bx; 40 ea/cs



Easy-peel liner: less force is needed to separate it from the adhesive section.


Ophthalmic drapes

Description

- Incise adhesive area optimised for skin
- Non-linting, tear-resistant, paper release liner
- 2 adhesive tabs for additional support
- Mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

Use

Cataract, vitreoretinal, corneal



Item number	Dimensions	Packaging
IDE1023	48 x 41 cm Adhesive only: 12 x 10 cm Adhesive tabs only: 5 x 5 cm	10 ea/bx; 40 ea/cs


Ophthalmic aperture pouches

Description

- Oval aperture with adhesive optimised for skin
- Non-linting, tear-resistant, paper release liner
- 2 adhesive tabs for additional support
- Mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

Use

Cataract, vitreoretinal, corneal



Item number	Dimensions	Packaging
IDE1024	48 x 41 cm Adhesive only: 6 x 3.5 cm Adhesive tabs only: 5 x 5 cm Aperture only: 6 x 3.5 cm	10 ea/bx; 40 ea/cs

Aperture drapes


Description

- Oval aperture with adhesive area
- Non-linting, tear-resistant, paper release liner
- Soft, drapeable, low-glare matte-finish film

Use

Cataract, vitreoretinal, corneal, lid/plastic, ENT

Item number	Dimensions	Packaging
IDE1030	83 x 129 cm Adhesive only: 12 x 11 cm Aperture only: 7 x 5 cm	10 ea/bx; 40 ea/cs



Ophthalmic incise drape with fluid pouch


Description

- Incise adhesive area
- Non-linting, tear-resistant, paper release liner
- Fluid pouch with mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

Use

Cataract, vitreoretinal

Item number	Dimensions	Packaging
IDE1061	121 x 129 cm Adhesive only: 12 x 11 cm Pouch only: 30 x 27 cm	10 ea/bx; 40 ea/cs



SMS ophthalmic incise drape with fluid pouch


Description

- Incise adhesive area
- Non-linting, tear-resistant, paper release liner
- Fluid pouch with mouldable strip for additional fluid control
- Soft, drapeable, repellent, low-linting SMS fabric

Use

Cataract, vitreoretinal

Item number	Dimensions	Packaging
IDE1062	122 x 129.5 cm Adhesive only: 13 x 11.5 cm Pouch only: 31 x 28 cm	10 ea/cs



Other drapes


Aperture pouches

Description

- Oval aperture with large adhesive area
- Strong adhesive with no residue
- Mouldable strip for additional fluid control
- Heavy, low-glare matte-finish film

Use

Vaginal, ophthalmic



Item number	Dimensions	Packaging
IDE1072	47 x 23 cm Adhesive only: 30 x 22 cm Aperture only: 11 x 9 cm	10 ea/bx; 40 ea/cs

Clear adhesive drape

Item number	Dimensions	Packaging
5810CE	150 x 300 cm	17 ea/cs


Urological drapes

Description

- Strong adhesive with no residue
- Latex-free finger cot
- Strong, low-glare matte-finish film
- Non-linting, tear-resistant, paper release liner

Use


TUR, urology, gynaecology



Item number	Dimensions	Packaging
IDE1071	60 x 91 cm Adhesive only: 60 x 5 cm Finger cot adhesive only: O.D. 10 x I.D. 3 cm	10 ea/bx; 40 ea/cs

Clear epidural drapes

Item number	Dimensions	Packaging
8461CE	60 x 100 cm, diameter 9 cm Adhesive fenestration	150 ea/cs
8465CE	80 x 140 cm, 6 x 10 cm Adhesive oval fenestration Pre-cut line for easy removal	150 ea/cs



Cesarean pouches


Description

- High-moisture, permeable incise film for extended adhesion
- Strong adhesive with no residue
- Large-cell foam to provide fluid control
- 360° fluid pouch for complete fluid containment
- 2 cord holders for line organisation
- Heavy, low-glare matte-finish film for additional pouch strength
- 2 suction ports with caps for improved fluid control

Use

C-section, trauma, general

Item number	Dimensions	Packaging
IDE2057	92 x 76 cm Adhesive only: 38 x 36 cm	10 ea/cs



Wound ring protectors


Description

- Strong, drapeable, low-glare matte-finish film
- 5 ring diameters to meet different surgical needs
- Application with adhesive tape

Use

General, laparotomy, cardiothoracic, colectomy

Item number	Dimensions	Packaging
IDE3707	90 x 90 cm, ring diameter: 7 cm	40 ea/cs
IDE3712	90 x 90 cm, ring diameter: 12 cm	40 ea/cs
IDE3718	90 x 90 cm, ring diameter: 18 cm	40 ea/cs
IDE3723	90 x 90 cm, ring diameter: 23 cm	40 ea/cs
IDE3728	90 x 90 cm, ring diameter: 28 cm	40 ea/cs


















The wound ring protectors have a flexible but rigid ring at the fenestration site to help keep the drape in the correct position. The frosted clear plastic drape around the fenestration helps maintain the sterile barrier.

Dermadine Plus and Dermadrape Incise Drapes

Description

- Thin, conformable, low-glare matte-finish film
- Strong adhesive with no residue
- Breathable, water-repellent and waterproof
- Comfortable and enveloping
- Long-lasting incise drapes
- Dermadine Plus: impregnated with an antibacterial solution
- Dermadrape: clear incise drapes



Dermadine item number		Dimensions	Packaging	Dermadrape item number	
42020P		20 x 20 cm Adhesive: 15 x 20 cm	20 ea/bx 120 ea/cs	32020	
n/a		35 x 25 cm Adhesive: 30 x 25 cm	20 ea/bx 120 ea/cs	33525	
44040P		40 x 40 cm Adhesive: 34 x 40 cm	15 ea/bx 90 ea/cs	34040	
n/a		50 x 48 cm Adhesive: 44 x 48 cm	15 ea/bx 90 ea/cs	35050	
47045P		70 x 45 cm Adhesive: 60 x 45 cm	15 ea/bx 90 ea/cs	37045	
47060P		70 x 60 cm Adhesive: 60 x 60 cm	10 ea/bx 60 ea/cs	n/a	
47080P		70 x 80 cm Adhesive: 60 x 80 cm	10 ea/bx 60 ea/cs	n/a	
48060P		80 x 60 cm Adhesive: 60 x 60 cm	10 ea/bx 60 ea/cs	38060	
48080P		80 x 80 cm Adhesive: 70 x 80 cm	10 ea/bx 60 ea/cs	38080	


Equipment drapes

Circular banded bags

Description

- Clear, drapeable polyethylene film
- Latex-free elastic bands
- Versatile, often used to cover intensifiers and other equipment

Item number	Dimensions	Stretched width	Packaging
ICE60020S	30 x 25 cm	55 cm	25 ea/cs
ICE60028S	45 x 35 cm	70 cm	25 ea/cs
ICE60030S	56 x 40 cm	80 cm	25 ea/cs
ICE60040S	58 x 50 cm	100 cm	25 ea/cs
ICE60050S	72 x 60 cm	120 cm	25 ea/cs
ICEBBG070	75 x 65 cm	130 cm	50 ea/cs
ICE60060S	140 x 70 cm	165 cm	10 ea/cs




Rectangular banded bags

Description

- Clear, drapeable polyethylene film
- Latex-free elastic bands
- Versatile, often used to cover intensifiers and other equipment

Item number	Dimensions	Packaging
ICE63036R	76 x 91 cm	25 ea/cs
ICE63628R	91 x 71 cm	25 ea/cs
ICE63636R	91 x 91 cm	25 ea/cs
ICEBBG095	95 x 85 cm	25 ea/cs
ICE63654S	95 x 130 cm	25 ea/cs
ICE64848S	122 x 122 cm	25 ea/cs
ICE65656S	142 x 142 cm	25 ea/cs



Banded bags for X-ray shields

Description

- Clear, drapeable polyethylene film
- Latex-free elastic bands
- A design suitable for most X-ray panels

Item number	Dimensions	Packaging
ICE90120	120 x 90 cm	25 ea/cs
ICE8090	90 x 80 cm	25 ea/cs



Medline's new banded bags help healthcare professionals cover protective shields (MAVIG), which are generally used in angiography procedures. These items provide a sterile barrier between the equipment and the patient undergoing the procedure.

C-arm set Philips

Description

- Clear banded bags for X-ray tank/intensifier
- Clear polyethylene drape to cover the C-arm


DYNJE4526 for BV29/300

DYNJE4527 for BV300

DYNJE4530 for BV25/26

DYNJE4569 for Pulsera/Endura

Item number	Dimensions	Packaging
DYNJE4526	Banded bags: 60 x 115 cm 70 x 130 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs
DYNJE4527	Banded bags: 60 x 115 cm rect.: 95 x 85 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs
DYNJE4530	Banded bags: 60 x 115 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs
DYNJE4569	Banded bags: 60 x 115 cm 70 x 140 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs




C-arm set Ziehm

Description

- Clear banded bags for X-ray tank/intensifier
- Clear polyethylene drape to cover the C-arm

For starter set SPTE-SS634

Item number	Dimensions	Packaging
DYNJE4531	Banded bags: 60 x 115 cm 70 x 130 cm C-arm cover: 35 x 185 cm	10 ea/bx; 20 ea/cs



C-arm set Siemens


Description

- Clear banded bags for X-ray tank/intensifier
- Clear polyethylene drape to cover the C-arm


DYNJE4532 for SIREMOBIL N4

- Clear polyethylene drape with clips for starter set SPTE-SS628

DYNJE4556 for ISO-C/Arcadis/Orbic, 2 parts



Item number	Dimensions	Packaging
DYNJE4532	Banded bags: 60 x 115 cm 60 x 115 cm C-arm cover: 35 x 183 cm	10 ea/bx; 20 ea/cs
DYNJE4547	Banded bags: 70 x 130 cm 60 x 115 cm C-arm cover: 45 x 183 cm	20 ea/cs
DYNJE4556	Top: 51 – 82 x 130 cm Bottom: 53 – 46 x 134 cm	10 ea/bx; 20 ea/cs




Fluoroscope drapes

Description

- 2 adhesive strips for attachment, with 2 adhesive tabs for proper fit
- Clear, drapeable polyethylene film

Item number	Dimensions	Packaging
ICE1012	91 x 111 cm Adhesive only (2): 91 x 1.2 cm Adhesive tabs only: 5 x 10 cm	10 ea/bx; 40 ea/cs



Full-size C-arm drapes

Description

- 3 adhesive polyethylene straps for securement
- Clear, drapeable polyethylene film

ICE4400R: has 2 latex-free elastic bands for additional securement

Item number	Dimensions	Packaging
ICE4400	104 x 188 cm	20 ea/cs
ICE4400R	104 x 188 cm	20 ea/cs
ICE4405	90 x 225 cm	75 ea/cs




Contoured C-arm drapes

Description

- 4 adhesive polyethylene straps for securement
- Clear, drapeable polyethylene film
- Contoured to fit full C-arm

Item number	Dimensions	Packaging
ICE4430	152 x 267 cm	20 ea/cs



Mini C-arm drapes


Description

- Held in place with an elastic band on the back
- 3 adhesive polyethylene straps for additional securement
- Clear, drapeable polyethylene film

ICE4410: fits Fluoroscan® mini C-arm models

ICE4420: fits Fluoroscan® and OEC6600/6800 mini C-arm models

Item number	Dimensions	Packaging
ICE4410	137 x 160 cm	20 ea/cs
ICE4420	137 x 216 cm	20 ea/cs




Full-size C-arm kits

Description

- Clear, drapeable polyethylene film
- Latex-free bands

Item number	Dimensions	Packaging
ICE4500	Intensifier drape (1): 91 x 91 cm C-Arm drape (1): 41 x 218 cm Footswitch drape (1): 46 x 41 cm Banded bag (1): 71 x 31 cm	20 ea/cs
ICE8130	C-Arm drape (1): 10/14 x 225 cm Banded bag (2): 75 x 90 cm	20 ea/cs
ICE8132	C-Arm drape (1): 10/14 x 225 cm Banded bag (1): 75 x 90 cm	20 ea/cs
ICE8230	C-Arm drape (1): 8/10 x 205 cm Banded bag (2): 75 x 75 cm	20 ea/cs
ICE8232	C-Arm drape (1): 8/10 x 205 cm Banded bag (1): 75 x 75 cm	20 ea/cs
ICE8330	C-Arm drape (1): 8/10 x 180 cm Banded bag (2): 75 x 75 cm	20 ea/cs
ICE8332	C-Arm drape (1): 8/10 x 180 cm Banded bag (1): 75 x 75 cm	20 ea/cs



Remote control covers

Description

- Clear polyethylene film bag
- Pinch-top closure

Item number	Dimensions	Packaging
DYNJE4635	10 x 28 cm	50 ea/bx; 200 ea/cs
DYNJE4641	18 x 35 cm	50 ea/bx; 200 ea/cs




Video camera covers

Description

- Adhesive polyethylene strap to hold it in place
- Adhesive strip at the base for additional securement
- Clear, drapeable polyethylene film
- Cable cover

Item number	Dimensions	Packaging
ICE4250	13 x 244 cm, reverse	25 ea/cs
ICE4300	23 x 244 cm, reverse	25 ea/cs
ICE4240	13 x 240 cm, telescopic with tab	25 ea/cs
ICE4260	15 x 240 cm, with ring	25 ea/cs
ICE4270	15 x 298 cm, concertina	25 ea/cs
ICE4200	18 x 244 cm, elastic end	25 ea/cs
ICE4100	8 x 200 cm, telescopic with tab	25 ea/cs
ICE4245	18 x 250 cm, telescopic with tab	25 ea/cs
ICE4280	18 x 250 cm, telescopic with closed system	25 ea/cs
ICE4290	14 x 250 cm, with ring and closed system	25 ea/cs



Microscope drapes

Description

- All fitted with clear, drapeable polyethylene film

ICE5010, ICE5040

- Fits 48 mm lens (Zeiss® OPMI)
- 3 adhesive PE straps for securement
- Fits 1 pair of oculars

ICE5200, ICE5210, ICE5030


- Fits 65 mm lens (Zeiss MD)
- 3 adhesive PE straps for securement
- Fits 3 pairs of oculars

ICE5220

- Fits 65 mm lens (Zeiss MD)
- 5 adhesive PE straps for securement
- Fits 3 pairs of oculars

ICE5300, ICE5330

- Fits 68 mm lens (Leica®)
- 3 adhesive PE straps for securement
- Fits 3 pairs of oculars



Item number	Dimensions	Packaging
-------------	------------	-----------

Fits 48 mm lens (Zeiss OPMI)

ICE5010	51 x 163 cm	20 ea/cs
---------	-------------	----------

ICE5040	51 x 267 cm	40 ea/cs
---------	-------------	----------

Fits 65 mm lens (Zeiss MD)

ICE5200	117 x 203 cm	10 ea/cs
---------	--------------	----------

ICE5030	117 x 305 cm	10 ea/cs
---------	--------------	----------

ICE5210	117 x 267 cm	10 ea/cs
---------	--------------	----------

ICE5220	137 x 381 cm	10 ea/cs
---------	--------------	----------

Fits 68 mm lens (Leica)

ICE5300	104 x 163 cm	10 ea/cs
---------	--------------	----------

ICE5320	137 x 267 cm	10 ea/cs
---------	--------------	----------

ICE5330	137 x 381 cm	10 ea/cs
---------	--------------	----------

Microscope handle covers

ICE5400	15 x 35 cm	20 ea/cs
---------	------------	----------



X-ray cassette covers

Description

- Clear, drapeable polyethylene film

ICE4010: release liner for facilitating sterile technique

Item number	Dimensions	Packaging
ICE4000	Large: 61 x 102 cm	20 ea/cs
ICE4010	Standard with liner: 61 x 64 cm	48 ea/cs
ICE4050	Standard: 61 x 64 cm	50 ea/cs
8720CEA	Large: 54 x 91 cm	48 ea/cs




Thermo cover

Description

- Clear polyethylene film
- Aluminium coated foil
- Temperature resistance up to 60° C

Item number	Dimensions	Packaging
DYNJE7000	68 x 86 cm	40 ea/cs



Probe cover kits

Description

- 1 probe cover: soft and drapeable polyurethane film
- 1 ultrasound gel packet: high-quality, sterile ultrasound gel for transducer
- 2 elastic bands for securement

Item number	Dimensions	Packaging
Kits with gel		
ICE5900	10 x 122 cm	20 ea/cs
ICE5910	15 x 61 cm	20 ea/cs
ICE5920	15 x 122 cm	20 ea/cs
ICE5930	15 x 244 cm	20 ea/cs

Kits without gel		
ICE18120	18 x 120 cm	25 ea/cs
ICE1830	18 x 30 cm	25 ea/cs
ICE1880	18 x 80 cm	25 ea/cs



High-quality polyurethane is soft and drapeable to conform to the transducer and has a smooth finish for easy application and glide. The slightly matte surface helps prevent glare.

Probe cover kits with adhesive

Description

- 1 probe cover: clear polyethylene film with adhesive tip
- 1 ultrasound gel packet; high-quality, sterile ultrasound gel for transducer
- 2 elastic bands and 2 adhesive tapes for additional securement

Item number	Dimensions	Packaging
Kits with gel		
US.20.180.040-G	18 x 40 cm	50 ea/bx; 200 ea/cs
US.20.180.080-G	18 x 80 cm	50 ea/bx; 200 ea/cs
US.20.180.120-G	18 x 120 cm	50 ea/bx; 200 ea/cs
US.20.180.160-G	18 x 160 cm	50 ea/bx; 200 ea/cs
Kits without gel		
US.20.180.040	18 x 40 cm	100 ea/bx; 200 ea/cs
US.20.180.080	18 x 80 cm	100 ea/bx; 200 ea/cs
US.20.180.120	18 x 120 cm	100 ea/bx; 200 ea/cs
US.20.180.160	18 x 160 cm	100 ea/bx; 200 ea/cs



Medline's SonoGuard Plus Sterile Ultrasound Probe Cover eliminates the need for contact gel use. The cover comes with adhesive ClearView film for attaching the ultrasound transducer head to the inside of the cover. This reduces setup time. And as a result of its non-slip properties, it is optional to use a second elastic band or adhesive strip for additional securement.




Cover for thyroid collar

Description

- Keeps thyroid collar sanitary and clean
- Comfortable spunlace fabric
- 2 adhesive tape strips for closing and securement
- Fits standard thyroid shields
- Easy to use and dispose of after use

Use

Medical and dental X-rays, orthopaedic, radiology, CT radiation, mammography



During procedures requiring fluoroscopy, an X-ray protective thyroid collar is worn to shield the wearer from the harmful effects of radiation. However, thyroid collars absorb sweat and begin to build up bacteria. Make-up can also rub off on the thyroid collar. Medline disposable thyroid collars were designed to wrap around and universally fit any thyroid collar shield. Our disposable thyroid collar cover will reduce or prevent the build-up of bacteria, make-up and other debris.

Item number	Dimensions	Packaging
SPTE-ICE9000	42 x 21.5 cm	25 ea/bx; 200 ea/cs

Stockinettes

Impervious stockinettes

Description

- Offers exceptional, high-level barrier properties and optimal fluid control, thanks to the use of an impervious thermoplastic elastomer material
- Gives the patient comfort by way of a soft, knitted jersey fabric on the inside.
- Available in 6 sizes

Use

Orthopaedic, arthroplasty, as well as hand, knee, ankle, foot, limb and vascular surgeries

Item number	Width 'A'	Width 'B'	Packaging
22520CE	15 cm	76 cm	20 ea/cs
22525CE	20 cm	97 cm	20 ea/cs
22530CE	23 cm	91 cm	20 ea/cs
22534CE	23 cm	122 cm	16 ea/cs
22535CE	30 cm	122 cm	12 ea/cs
22540CE	36 cm	122 cm	12 ea/cs



Impervious stretchable stockinettes


Description

- Allows for an ideal fit over the patient's anatomy, due to the use of a stretchable, bilaminate material
- Facilitates set-up with the help of 2 pulling ties
- Both lower and upper extremities can be covered, as it is available in 2 sizes

Use

Orthopaedic, arthroplasty, as well as hand, knee, ankle, foot, limb and vascular surgeries

Item number	Width 'A'	Width 'B'	Packaging
22620CE	17 cm	76 cm	18 ea/cs
22634CE	24 cm	118 cm	10 ea/cs



Invisishield cross reference

Page	Description	Product number Medline	Product number 3M®	Product number Microtek®
3	Towel drape	IDE1000	1000	ISO1000
4	Isolation bag	ICE1003	1003	ISO1003
3	Towel drape	IDE1010	1010	ISO1010
8	Patient isolation drape	IDE1014	1014	n/a
7	Polyethylene U drape	IDE1015	1015	ISO1015
4	Fluid pouch	IDE1016	1016	ISO1016
9	Patient isolation drape	IDE1017	1017	5960
9	Wide patient isolation drape	IDE1017W	n/a	6960
5	Instrument pouch	IDE1018	1018	ISO1018
5	Long instrument pouch	IDE1018L	1018L	n/a
7	Polyethylene U drape	IDE1019	1019	ISO1019
10	Aperture drape	IDE1020	1020	ISO1020
10	Ophthalmic incise drape	IDE1023	1023	n /a
11	Ophthalmic aperture pouch	IDE1024	1024	ISO1024
11	Aperture drape	IDE1030	1030	ISO1030
14	Incise drape	32020	1035	ISO1035
14	Incise drape	34040	1040	ISO1040
14	Incise drape	37045	1050	ISO1050
14	Incise drape	38080	1051	ISO1051
12	Ophthalmic incise drape with fluid pouch	IDE1061	1061	ISO1061
7	Shoulder pouch	IDE1067	1067	ISO1067
12	Urological drape	IDE1071	1071	ISO1071
13	Aperture pouch	IDE1072	1072	ISO1072
8	U fluid collection pouch	DYJPEFP090	n/a	n/a
6	Instrument pouch	DYJPEISTPOM1	n/a	n/a
6	Instrument pouch	DYJPEISTPOM2	n/a	n/a
6	Instrument pouch	DYJPEISTPOM3	n/a	n/a



All products are latex free and sterile, unless indicated in the description.



Adhesive



Matte film



Polyethylene
film

Dermadine Plus cross reference

Type of film	Description	Product number Medline	Product number 3M®
PU	Incise, sterile drape 15 x 20 cm; incise area 10 x 20 cm, Iodophor impregnated	42020P	6635
PU	Incise, sterile drape 44 x 35 cm; incise area 34 x 35 cm, Iodophor impregnated	44040P	6640
PU	Incise, sterile drape 66 x 60 cm; incise area 56 x 60 cm, Iodophor impregnated	47060P	6648
PU	Incise, sterile drape 66 x 45 cm; incise area 56 x 45 cm, Iodophor impregnated	47045P	6650
PU	Incise, sterile drape 66 x 85 cm; incise area 56 x 85 cm, Iodophor impregnated	47080P	6651
PU with applicator	Incise, sterile drape 30 x 20 cm; incise area 26 x 20 cm, Iodophor impregnated	42020P	6661EZ
PU with applicator	Incise, sterile drape 60 x 35 cm; incise area 35 x 35 cm, Iodophor impregnated	44040P	6640EZ
PU with applicator	Incise, sterile drape 90 x 60 cm; incise area 60 x 60 cm, Iodophor impregnated	47060P	6648EZ
PU with applicator	Incise, sterile drape 90 x 45 cm; incise area 60 x 45 cm, Iodophor impregnated	47045P	6650EZ
PU with applicator	Incise, sterile drape 90 x 85 cm; incise area 60 x 85 cm, Iodophor impregnated	48080P	6651EZ

Dermadrape cross reference

Type of film	Description	Product number Medline	Product number 3M®
PE	Incise, sterile drape 15 x 20 cm; incise area 10 x 20 cm	32020	1035
PE	Incise, sterile drape 29 x 25 cm; incise area 19 x 25 cm	32020	1037
PE	Incise, sterile drape 60 x 35 cm; incise area 35 x 35 cm	34040	1040
PE	Incise, sterile drape 82 x 45 cm; incise area 60 x 45 cm	37045	1050
PE	Incise, sterile drape 82 x 85 cm; incise area 60 x 85 cm	38080	1051
PU	Incise, sterile drape 20 x 20 cm; incise area 10 x 20 cm	32020	2035
PU	Incise, sterile drape 38 x 25 cm; incise area 28 x 25 cm	33525	2037
PU	Incise, sterile drape 38 x 41 cm; incise area 28 x 41 cm	34040	2040
PU	Incise, sterile drape 60 x 45 cm; incise area 50 x 45 cm	35050	2045
PU	Incise, sterile drape 90 x 60 cm; incise area 60 x 60 cm	38060	2050
PU	Incise, sterile drape 90 x 85 cm; incise area 60 x 85 cm	38080	2051



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com



Medline Ireland Ltd
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06
www.medline.eu/ie
ie-customerservice@medline.com


FOLLOW US  

Isolation bag 586560-1 on page 4 is a class I sterile medical device.  SENGEWALD Klinikprodukte GmbH
Adressstr. 2
D-83101 Rohrdorf-Thansau, Germany  0123

Dermadine Plus is a class III medical device, and Dermadrape is a class IIa sterile medical device.  TIASET di Lorenzo Carresi
SP Loc. Botriolo Km 5100
52026 Castelfranco-Piandisco (AR), Italy  0373 Istituto Superiore di Sanità

Wound ring protectors are class IIa sterile medical devices.  Medline International France SAS
44110 Châteaubriant, France  2797 (BSI)

C-arm covers ICE4405, ICE81xx, ICE82xx and ICE83xx on pages 19 and 20 and camera covers ICE4280 and ICE4290 are class I sterile medical devices.  Medbar Tıbbi Malzemeler Turizm San. ve Tic. A.Ş.
Gazimir-Izmir, Turkey  1984

Probe cover kits with gel ICE59xx on page 23 are procedure packs according to Article 12 of the MDD 93/42/EEC or Article 22 of the MDR 2017/745. Procedure pack producer:  Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France

Probe cover kits with adhesive on page 24 are class I sterile medical devices intended to be used by healthcare professionals.

Kits without gel:  Promecon GmbH
Kreuzerstraße 10
20095 Hamburg, Germany  0537

Kits with gel: Procedure pack according to Article 12 of the MDD 93/42/EEC  Probe cover: Promecon GmbH
Kreuzerstraße 10
20095 Hamburg, Germany   Get: Turkcaz Medikal
Akçaburgaz, Muhsin Yazıcıoğlu Cd. No:45/5
34522 Esenyurt/Istanbul, Turkey  2292

All other products listed are class I sterile medical devices, except the cover for thyroid collar on page 25.  Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France  2797 (BSI)

These products are intended to be used by healthcare professionals.
Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2022 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML164_EN 05/2022.

Procedure pack producer:
Promecon GmbH
Kreuzerstraße 10
20095 Hamburg, Germany