



PATIENT DISPOSABLE BASICS



Perioperative and Operating Theatre Supplies

Meeting your medical needs with unsurpassed quality and selection

Disposable OR Towels

Our Disposable OR Towels are manufactured to the highest specifications and are available in both green and blue. The portfolio offers X-ray and non-X-ray detectable towels to meet different customer needs.

- » Absorbent
- » 90% cotton, 10% viscose
- » Pre-washed



Item No.	NHS Code	Description	Packaging
28200004CEA	VJD973	4 ABSORBENT COTTON OR TOWELS GREEN 66 x 43 cm	4/pk, 20/cs
28700002CEA	VJD327	2 ABSORBENT COTTON OR TOWELS BLUE 66 x 43 cm	2/pk, 40/cs
28700004CEA	VEB128	4 ABSORBENT COTTON OR TOWELS BLUE 66 x 43 cm	4/pk, 20/cs
28700002XR	VEB127	2 ABSORBENT X-RAY COTTON OR TOWELS BLUE 66 x 43 cm	2/pk, 40/cs
28700004XR	VEB129	4 ABSORBENT X-RAY COTTON OR TOWELS BLUE 66 X 43 cm	4/pk, 20/cs

Sterile ABD Pads

Medline's absorbent ABD pads feature a soft non-woven outer layer that quickly wicks fluid to a cellulose centre. This thick layer of cellulose quickly absorbs and disperses fluids laterally to prevent pooling. The hydrophobic back impedes strikethrough, and all edges are sealed to keep moisture in the dressing and prevent linting.

- » Superabsorbent
- » Soft non-woven outer layer
- » Cellulose centre to prevent pooling
- » Hydrophobic back to impede strikethrough
- » Sealed edges to prevent leakage and linting

Item No.	NHS Code	Description	Packaging
ABD1010	EJA099	ABD PAD ST 10 x 10 cm	1/pch, 25/bx, 300/cs
ABD1020	EJA101	ABD PAD ST 10 x 20 cm	1/pch, 25/bx, 200/cs
ABD1525	n/a	ABD PAD ST 15 x 25 cm	1/pch, 15/bx, 180/cs
ABD2020	EJA102	ABD PAD ST 20 x 20 cm	1/pch, 15/bx, 120/cs
ABD2040	EJA103	ABD PAD ST 20 x 40 cm	1/pch, 10/bx, 80/cs



Cotton Applicators

Cotton applicators are devices used to apply a substance, e.g. medication, to the ear/nose/throat area or for wound care treatments.

- » Tightly wound heads for durability and confidence in sample selection
- » Choice of sterile and non-sterile

Item No.	Description	Packaging
830000	SMALL NON-STERILE COTTON APPLICATOR	100/pk, 10/bx, 20/cs
830010	SMALL STERILE COTTON APPLICATOR	2/pk, 100/pch, 15/cs
831510	BIG STERILE COTTON APPLICATOR	2/pk, 50/pch, 15/cs
831578	BIG NON-STERILE COTTON APPLICATOR	50/pk, 20/cs



Lap Sponges

Medline offers a broad range of quality lap sponges in a variety of sizes for maximum absorbency. All sponges are delinted and X-ray detectable to promote operating theatre safety.

- » Pre-washed, delinted
- » 100% cotton
- » Highly absorbent
- » X-ray detectable

Sterile Lap Sponges, White, Pre-Washed

Item No.	Description	Packaging
LW2P089004	LAP SPONGE WHITE ST 8 x 90 cm 4PLY	2/pk, 50/cs
LW2P203004	LAP SPONGE WHITE ST 20 x 30 cm 4PLY	2/pk, 90/cs
LW2P454504	LAP SPONGE WHITE ST 45 x 45 cm 4PLY	2/pk, 60/cs
LW2P454506	LAP SPONGE WHITE ST 45 x 45 cm 6PLY	2/pk, 40/cs
LW5P454504	LAP SPONGE WHITE ST 45 x 45 cm 4PLY	5/pk, 20/cs
LW5P454506	LAP SPONGE WHITE ST 45 x 45 cm 6PLY	5/pk, 20/cs



Sterile Lap Sponges, Green, Pre-Washed

Item No.	Description	Packaging
LG2P089004	LAP SPONGE GREEN ST 8 x 90 cm 4PLY	2/pk, 50/cs
LG2P454506	LAP SPONGE GREEN ST 45 x 45 cm 6PLY	2/pk, 40/cs
LG5P454504	LAP SPONGE GREEN ST 45 x 45 cm 4PLY	5/pk, 20/cs
LG5P454506	LAP SPONGE GREEN ST 45 x 45 cm 6PLY	5/pk, 20/cs



Lap Sponge Counter

Medline's lap sponge counter system helps promote a safer operating theatre by helping you track all of your sponges following a surgical procedure.

- » Easy and accurate counting bags
- » Assists and organises the counting process in the operating theatre
- » Advanced sponge counter rack for holding 4 sets of bags at once

Item No.	Description	Packaging
DYNJE1000	SPONGE COUNTER BAGS, CLEAR	1/pch, 50/bx, 250/cs
DYNJE1001	SPONGE COUNTER BAGS, BLUE	1/pch, 50/bx, 250/cs
NON50511	SPONGE COUNTER BAG RACK	1/bx, 1/cs



Does not come with an IV pole

HELP MAKE YOUR OPERATING THEATRE SAFE, ONE TOWEL AND LAP SPONGE AT A TIME.

When bacteria attaches to lint particles from drapes, gowns and towels, they can contaminate instruments, devices, gloves or even an open surgical wound. This can result in infections and other complications leading to poor outcomes. To help prevent this from occurring, Medline provides a portfolio with high-quality OR towels and lap sponges.

You can always count on Medline to deliver the highest quality products at the best value, where and when you need them.

Innovation

Medline's heritage of quality and innovation dates back to 1910. From the introduction of coloured surgeon's gowns back in the 1940s to our new OR towels, Medline is at the forefront of enhancing safety for patients and staff.

Quality

As a manufacturer, we have complete control over the entire production process, including the raw materials and manufacturing methods. We have more than 550 quality inspectors worldwide who ensure that we maintain the strictest manufacturing and quality standards, which meet or exceed different local regulations.

Value

Medline delivers directly to end users, which eliminates the middleman, providing you with the best prices and highest value.

Selection

Medline offers a variety of OR towels and lap sponges to meet the needs of different customers.



Contact your account manager or call
+44 844 334 52 37 to order our products today.



**ALWAYS
ON.**

Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238

www.medline.eu/uk
uk-customerservice@medline.com

Lap Sponge is a class IIa sterile medical device.



Sterile ABD Pad is a class I sterile medical device.



OR Towel is a class I sterile medical device.



Cotton Applicator is a class I medical device.



These products are intended to be used by healthcare professionals. Before use, please consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct any errors that may have occurred.

© 2018 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. ML419-EN06/MG 12/2018.