


Operating Theatre Supplies

Selection of OR Towels, Lap Sponges, Gauzes and ABD Pads



**Basic
disposables
for patient
care**





Disposable OR Towels

Medline's Disposable OR Towels are manufactured to the highest specifications and are available in both green and blue. The portfolio offers X-ray and non-X-ray detectable towels to meet different customer needs.

- Absorbent
- 100% cotton
- Pre-washed



Item number	NHS Code	Description	Packaging
28200004CEA	VJD973	OR Towels; absorbent, cotton, green - 61 x 43 cm	4/pk, 20/cs
28700002CEA	VJD327	OR Towels; absorbent, cotton, blue - 61 x 43 cm	2/pk, 40/cs
28700004CEA	VEB128	OR Towels; absorbent, cotton, blue - 61 x 43 cm	4/pk, 20/cs
28700002XR	VEB127	OR Towels; absorbent, x-ray, cotton, blue - 61 x 43 cm	2/pk, 40/cs
28700004XR	VEB129	OR Towels; absorbent, x-ray, cotton, blue - 61 x 43 cm	4/pk, 20/cs



Lap Sponges

Medline's quality Lap Sponges offer maximum absorbency and come in a variety of sizes. All sponges are delinted and X-ray detectable to promote operating theatre safety. These sponges are:

- Pre-washed, delinted
- 100% cotton
- Highly absorbent
- X-ray detectable

Sterile Lap Sponges, White, Pre-Washed

Item number	Description	Packaging
LW2P089004	Lap Sponge; white, 4-ply, sterile - 8 x 90 cm	2/pk, 50/cs
LW2P089006	Lap Sponge; white, 6-ply, sterile - 8 x 90 cm	2/pk, 50/cs
LW2P203004	Lap Sponge; white, 4-ply, sterile - 20 x 30 cm	2/pk, 90/cs
LW5P203004	Lap sponge; white, 4-ply, sterile - 20 x 30 cm	5/pk, 45/cs
LW2P454504	Lap Sponge; white, 4-ply, sterile - 45 x 45 cm	2/pk, 60/cs
LW5P454504	Lap sponge; white, 4-ply, sterile - 45 x 45 cm	5/pk, 20/cs
SS5P454504A	Lap sponge; white, 4-ply, sterile - 45 x 45 cm*	5/pk, 20/cs
LW2P454506	Lap Sponge; white, 6-ply, sterile - 45 x 45 cm	2/pk, 40/cs
LW5P454506	Lap sponge; white, 6-ply, sterile - 45 x 45 cm	5/pk, 20/cs
LW2P506004	Lap sponge; white, 4-ply, sterile - 50 x 60 cm	2/pk, t45/cs

X-Ray detectable for better tracking

100% cotton, pre-washed, soft and absorbent



Sterile Lap Sponges, Green, Pre-Washed

Item number	Description	Packaging
LG2P089004	Lap Sponge; green, 4-ply, sterile - 8 x 90 cm	2/pk, 50/cs
LG2P089006	Lap sponge; green, 6-ply, sterile - 8 x 90 cm	2/pk, 50/cs
LG5P203004	Lap sponge; green, 4-ply, sterile - 20 x 30 cm	5/pk, 30/cs
LG2P454504	Lap sponge; green, 4-ply, sterile - 45 x 45 cm	2/pk, 60/cs
LG5P454504	Lap Sponge; green, 4-ply, sterile - 45 x 45 cm	5/pk, 20/cs
LG2P454506	Lap sponge; green, 6-ply, sterile - 45 x 45 cm	2/pk, 40/cs
LG5P454506	Lap Sponge; green, 6-ply, sterile - 45 x 45 cm	5/pk, 20/cs
LG2P506004	Lap sponge; green, 4-ply, sterile - 50 x 60 cm	2/pk, 45/cs
LG2P506006	Lap sponge; green, 6-ply, sterile - 50 x 60 cm	2/pk, 30/cs



Also available in green to conceal blood stains

Loop sewed into the sponge for easy removal and better tracking



Non-Sterile Sponge Counter Bags

Medline's lap sponge counter system helps promote a safer operating theatre by helping you track all of your sponges, following a surgical procedure. Our non-sterile counter bags offer:

- Ease of use and accuracy
- Assists and organises the counting process in the operating theatre
- Advanced sponge counter rack for holding 4 sets of bags at once

Item number	Description	Packaging
DYNJE1000	Sponge Counter Bag; clear	1/pch, 50/bx, 250/cs
DYNJE1001	Sponge Counter Bag; blue	1/pch, 50/bx, 250/cs
NON50511	Sponge Counter Bag Rack	1/bx, 1/cs




Does not come with an IV pole

Sterile Sponge Counter Bags

Medline's sterile Sponge Counter Bags allow for safe, effective, convenient counting, handling and disposal of contaminated swabs. They minimise the risk of pathogen inhalation and keep the surgical area clean and organised. Our sterile counter bags feature:

- Clear, drapeable polyethylene film
- 1 adhesive strip for attachment
- 3 sizes
- 3 options for number of compartments

Item number	Description	Packaging
IBE5242	52 x 42 cm, 40 compartments	50 ea/cs
IBE5637	56 x 37 cm, 6 compartments	50 ea/cs
IBE5652	56 x 52 cm, 20 compartments	50 ea/cs



Cotton Applicators

Medline's Cotton Applicators are devices used to apply a substance, for example, medication, to the ear/nose/throat area or for wound care treatments. With our applicators, you have:

- Tightly wound heads for durability and confidence in sample selection
- The choice of sterile and non-sterile

Item number	Description	Packaging
830000A	Cotton Applicator; small, non-sterile	100/pk, 10/bx, 20/cs
830010A	Cotton Applicator; small, sterile	2/pk, 100/pch, 15/cs
831510A	Cotton Applicator; big, sterile	2/pk, 50/pch, 15/cs
831578A	Cotton Applicator; big, non-sterile	50/pk, 20/cs



Gauzes

Medline offers a variety of gauze swabs, available in multiple sizes and plies. Gauze is ideal for wound dressings, wound packing and general wound care.

Woven Gauzes

Item number	Size	Ply	Sterility	XR	Packaging
GS10P102012XR	10 x 20 cm	12	Sterile	Yes	10 pcs/pack 60 packs/case
GS10P757512XR	7.5 x 7.5 cm	12	Sterile	Yes	10 pcs/pack 60 packs/case
GS2P101012	10 x 10 cm	12	Sterile	No	2 pcs/pack 800 packs/case
GS5P050508	5 x 5 cm	8	Sterile	No	5 pcs/pack 540 packs/case
GSNS101008	10 x 10 cm	8	Non-sterile	No	100 pcs/pack 10 packs/case



Non-woven Gauzes

Item number	Size	Ply	Sterility	XR	Packaging
NWNS30101004	10 x 10 cm	4	Non-sterile	No	100 pcs/pack 10 packs/case

Sterile ABD Pads

Medline's absorbent ABD Pads feature a soft non-woven outer layer that quickly wicks fluid to a cellulose centre. This thick layer of cellulose quickly absorbs and disperses fluids laterally to prevent pooling. The hydrophobic back impedes strikethrough, and all edges are sealed to keep moisture in the dressing and prevent linting. General benefits and features of these pads include:

- Superabsorbent
- Soft non-woven outer layer
- Cellulose centre to prevent pooling
- Hydrophobic back to impede strikethrough
- Sealed edges to prevent leakage and linting

Item number	NHS Code	Description	Packaging
ABD1010	EJA099	ABD Pad; sterile - 10 x 10 cm	1/pch, 25/bx, 300/cs
ABD1020	EJA101	ABD Pad; sterile - 10 x 20 cm	1/pch, 25/bx, 200/cs
ABD1525	n/a	ABD Pad; sterile - 15 x 25 cm	1/pch, 15/bx, 180/cs
ABD2020	EJA102	ABD Pad; sterile - 20 x 20 cm	1/pch, 15/bx, 120/cs
ABD2040	EJA103	ABD Pad; sterile - 20 x 40 cm	1/pch, 10/bx, 80/cs



Help make your operating theatre safe, one towel and lap sponge at a time.

When bacteria attach to lint particles from drapes, gowns and towels, they can contaminate instruments, devices, gloves or even an open surgical wound. This can result in infections and other complications leading to poor outcomes. To help prevent this from occurring, Medline provides a portfolio with high-quality OR towels and lap sponges.

You can always count on Medline to deliver the highest quality products, at the best value, wherever and whenever you need them.

Innovation

Medline's heritage of quality and innovation dates back to 1910. From the introduction of coloured gowns for surgeons back in the 1940s, to our new OR towels, Medline is at the forefront of enhancing safety for patients and staff.

Quality

As a manufacturer, we have complete control over the entire production process, including raw material and manufacturing methods. We have more than 550 quality inspectors worldwide, who ensure that we maintain the strictest manufacturing and quality standards, which meet or exceed different local regulations.

Value

Medline delivers directly to end users, which eliminates the middleman, providing you with the best prices and highest value.

Selection

Medline offers a variety of OR towels and lap sponges to meet the needs of different customers.



For more information on our products, please contact your Medline account manager or visit us at: medline.eu/uk



Medline Industries Ltd.
 3rd Floor
 Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom
 Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Industries Ltd.
 1 Grant's Row
 Mount Street Lower
 Dublin 2
 Ireland
 Tel.: +353 1 691 73 06
www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

- | | | |
|--|--|--|
| Sterile Sponge Counter Bags are class I sterile medical devices. |  Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France |  2797 (BSI) |
| Lap Sponge is a class IIa sterile medical device. |  Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France |  2797 (BSI) |
| Sterile ABD Pad is a class I sterile medical device. |  Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France |  2797 (BSI) |
| OR Towel is a class IIa sterile medical device. |  Shanghai International Holding Corp. GmbH (Europe), 20527 Hamburg, Germany |  0123 (TUV SUD) |
| Cotton Applicator is a class I non-sterile medical device. |  CMC Medical Devices & Drugs S.L., C/ Maracaibo Norte, CP29006, Málaga, Spain |  |
| Woven Gauzes are class IIa sterile medical devices. |  Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France |  2797 (BSI) |
| Non-woven Gauzes are class I non-sterile medical devices. |  Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France |  |

These products are intended to be used by healthcare professionals. Before use, please consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.
 © 2024 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML419_EN 03/2024