

# Operating Theatre Supplies

Selection of OR Towels, Lap Sponges and ABD Pads



**Basic  
disposables  
for patient  
care**



Disposable OR Towels

Medline’s Disposable OR Towels are manufactured to the highest specifications and are available in both green and blue. The portfolio offers X-ray and non-X-ray detectable towels to meet different customer needs. These towels are:

- Absorbent
- 100% cotton
- Pre-washed



Item number	NHS Code	Description	Packaging
28200004CEA	VJD973	OR Towels; absorbent, cotton, green - 61 x 43 cm	4/pk, 20/cs
28700002CEA	VJD327	OR Towels; absorbent, cotton, blue - 61 x 43 cm	2/pk, 40/cs
28700004CEA	VEB128	OR Towels; absorbent, cotton, blue - 61 x 43 cm	4/pk, 20/cs
28700002XR	VEB127	OR Towels; absorbent, x-ray, cotton, blue - 61 x 43 cm	2/pk, 40/cs
28700004XR	VEB129	OR Towels; absorbent, x-ray, cotton, blue - 61 x 43 cm	4/pk, 20/cs

Sterile ABD Pads

Medline’s absorbent ABD Pads feature a soft non-woven outer layer that quickly wicks fluid to a cellulose centre. This thick layer of cellulose quickly absorbs and disperses fluids laterally to prevent pooling. The hydrophobic back impedes strikethrough, and all edges are sealed to keep moisture in the dressing and prevent linting. General benefits and features of these pads include:

- Superabsorbent
- Soft non-woven outer layer
- Cellulose centre to prevent pooling
- Hydrophobic back to impede strikethrough
- Sealed edges to prevent leakage and linting



Item number	NHS Code	Description	Packaging
ABD1010	EJA099	ABD Pad; sterile - 10 x 10 cm	1/pch, 25/bx, 300/cs
ABD1020	EJA101	ABD Pad; sterile - 10 x 20 cm	1/pch, 25/bx, 200/cs
ABD1525	n/a	ABD Pad; sterile - 15 x 25 cm	1/pch, 15/bx, 180/cs
ABD2020	EJA102	ABD Pad; sterile - 20 x 20 cm	1/pch, 15/bx, 120/cs
ABD2040	EJA103	ABD Pad; sterile - 20 x 40 cm	1/pch, 10/bx, 80/cs

Cotton Applicators

Medline’s Cotton Applicators are devices used to apply a substance, for example, medication, to the ear/nose/throat area or for wound care treatments. With our applicators, you have:

- Tightly wound heads for durability and confidence in sample selection
- The choice of sterile and non-sterile



Item number	Description	Packaging
830000	Cotton Applicator; small, non-sterile	100/pk, 10/bx, 20/cs
830010	Cotton Applicator; small, sterile	2/pk, 100/pch, 15/cs
831510	Cotton Applicator; big, sterile	2/pk, 50/pch, 15/cs
831578	Cotton Applicator; big, non-sterile	50/pk, 20/cs

Lap Sponges

Medline’s quality Lap Sponges offer maximum absorbency and come in a variety of sizes. All sponges are delinted and X-ray detectable to promote operating theatre safety. These sponges are:

- Pre-washed, delinted
- 100% cotton
- Highly absorbent
- X-ray detectable



Sterile Lap Sponges, White, Pre-Washed

Item number	Description	Packaging
LW2P089004	Lap Sponge; white, 4-ply, sterile - 8 x 90 cm	2/pk, 50/cs
LW2P203004	Lap Sponge; white, 4-ply, sterile - 20 x 30 cm	2/pk, 90/cs
LW2P454504	Lap Sponge; white, 4-ply, sterile - 45 x 45 cm	2/pk, 60/cs
LW2P454506	Lap Sponge; white, 6-ply, sterile - 45 x 45 cm	2/pk, 40/cs
LW5P454504	Lap Sponge; white, 4-ply, sterile - 45 x 45 cm	5/pk, 20/cs
LW5P454506	Lap Sponge; white, 6-ply, sterile - 45 x 45 cm	5/pk, 20/cs

Sterile Lap Sponges, Green, Pre-Washed

Item number	Description	Packaging
LG2P089004	Lap Sponge; green, 4-ply, sterile - 8 x 90 cm	2/pk, 50/cs
LG2P454506	Lap Sponge; green, 6-ply, sterile - 45 x 45 cm	2/pk, 40/cs
LG5P454504	Lap Sponge; green, 4-ply, sterile - 45 x 45 cm	5/pk, 20/cs
LG5P454506	Lap Sponge; green, 6-ply, sterile - 45 x 45 cm	5/pk, 20/cs

Non-sterile sponge counter bags

Medline’s lap sponge counter system helps promote a safer operating theatre by helping you track all of your sponges, following a surgical procedure. Our non-sterile counter bags offer:

- Ease of use and accuracy
- Assists and organises the counting process in the operating theatre
- Advanced sponge counter rack for holding 4 sets of bags at once

Item number	Description	Packaging
DYNJE1000	Sponge Counter Bag; clear	1/pch, 50/bx, 250/cs
DYNJE1001	Sponge Counter Bag; blue	1/pch, 50/bx, 250/cs
NON50511	Sponge Counter Bag Rack	1/bx, 1/cs




Sterile sponge counter bags

Medline’s sterile Sponge Counter Bags allow for safe, effective, convenient counting, handling and disposal of contaminated swabs. They minimise the risk of pathogen inhalation and keep the surgical area clean and organised. Our sterile counter bags feature:

- Clear, drapeable polyethylene film
- 1 adhesive strip for attachment
- 3 sizes
- 3 options for number of compartments

Item number	Description	Packaging
IBE5242	52 x 42 cm, 40 compartments	50 ea/cs
IBE5637	56 x 37 cm, 6 compartments	50 ea/cs
IBE5652	56 x 52 cm, 20 compartments	50 ea/cs





# Help make your operating theatre safe, one towel and lap sponge at a time.

When bacteria attach to lint particles from drapes, gowns and towels, they can contaminate instruments, devices, gloves or even an open surgical wound. This can result in infections and other complications leading to poor outcomes. To help prevent this from occurring, Medline provides a portfolio with high-quality OR towels and lap sponges.

You can always count on Medline to deliver the highest quality products, at the best value, wherever and whenever you need them.

## Innovation

Medline's heritage of quality and innovation dates back to 1910. From the introduction of coloured gowns for surgeons back in the 1940s, to our new OR towels, Medline is at the forefront of enhancing safety for patients and staff.

## Quality

As a manufacturer, we have complete control over the entire production process, including raw material and manufacturing methods. We have more than 550 quality inspectors worldwide, who ensure that we maintain the strictest manufacturing and quality standards, which meet or exceed different local regulations.

## Value

Medline delivers directly to end users, which eliminates the middleman, providing you with the best prices and highest value.

## Selection

Medline offers a variety of OR towels and lap sponges to meet the needs of different customers.



For more information on our products, please contact your Medline account manager or visit us at: [uk.medline.eu](http://uk.medline.eu)



### Medline Industries Ltd

3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
Fax: +44 844 334 5238

[uk.medline.eu](http://uk.medline.eu) | [uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

FOLLOW US  

Lap Sponge is a class IIa sterile medical device.

 Medline International France SAS  
1 rue Charles Lindbergh  
44100 Châteaubouillant, France  2797 (BSI)

Sterile ABD Pad is a class I sterile medical device.

 Medline International France SAS  
1 rue Charles Lindbergh  
44100 Châteaubouillant, France  2797 (BSI)

OR Towel is a class IIa sterile medical device.

 Shanghai International Holding Corp. GmbH (Europe)  
Löffelstraße 80  
20037 Hamburg, Germany  0123 (TUV SUD)

Cotton Applicator is a class I medical device.

 Applimed SA  
Zone Industrielle, route de Pra de Plan 1  
CH 1618 Châblé St Denis, Switzerland  1639 SCS Belgium NV

These products are intended to be used by healthcare professionals. Before use, please consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2020 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML419\_EN 10/2020