

A close-up photograph of a laparoscopic surgical procedure. A surgeon's hand, wearing a white sterile glove, is holding a laparoscopic grasper. The grasper is positioned near a patient's abdomen, which is covered with blue sterile drapes. Several other laparoscopic ports are visible, with instruments inserted into the abdominal cavity. The background is slightly blurred, showing more of the surgical field and equipment.

**For safe and
uninterrupted
laparoscopy**

Essential Laparoscopic Suite





As a large-scale provider of operating theatre-related products, laparoscopy is a core part of Medline's expertise. Our suite of fundamental laparoscopic products enables you to perform minimally invasive procedures, with flexibility and consistency. It is a comprehensive solution that facilitates the achievement of continuity in the operating theatre while focusing on safe surgery.


MED-LAP ORIS Insufflation Tubing Set

For proper CO₂ insufflation

The ORIS Insufflation Tubing Set completes Medline's laparoscopic portfolio. By filtering and insufflating clean CO₂, it allows clinicians to effectively establish pneumoperitoneum prior to and during laparoscopic procedures.

Medline offers a line of items that are fully compatible for use in conjunction with different insufflators.

Image	ORIS Insufflation Tubing Set	Compatibility range
	ORISISO with ISO 15 M/22 F connector	Karl Storz Richard Wolf Stryker ConMed
	ORIS15 with 15-mm male connector	Karl Storz Richard Wolf
	ORISLL with Leur-lock connector	Olympus
	ORISPUSH with flexible female connector	Stryker



Benefits and features:


- Offers flexible tubing and a highly secure filter
- Filters with 99.9% BFE and VFE
- Fully compatible range

Item code	Description	Length	Packaging
ORISISO	Operating theatre insufflation set with ISO 15 M/22 F connector	3 m	25 ea/cs
ORISS15	Operating insufflation set with 15-mm connector	3 m	25 ea/cs
ORISLL	Operating theatre insufflation set with Luer-lock connector	3 m	25 ea/cs
ORISPUSH	Operating theatre insufflation set with female connector	3 m	25 ea/cs

MED-LAP Veress Needle

For safe establishment of the pneumoperitoneum

The insufflation needle is equipped with a retracting, spring-loaded blunt stylet and is designed for simple and safe establishment of pneumoperitoneum during laparoscopic procedures.



Benefits and features:


- Non-abrasive lint-free applicator sponge
- Ready-to-use kit: bottle and sponge
- Screw-cap bottle with mild alcohol solution

Item code	Description	Length	Gauge	Insufflation eyes	Packaging
ORVN12	Veress insufflation needle	120 mm	14	2	12 ea/cs
ORVN15	Veress insufflation needle	150 mm	14	2	12 ea/cs

MED-LAP Anti-Fog Solution

For preventing foggy laparoscopic lenses

The Anti-Fog Solution prevents laparoscopic lenses from misting, thereby providing an unobstructed view during laparoscopic procedures.



Benefits and features:


- Non-abrasive lint-free applicator sponge
- Ready-to-use kit: bottle and sponge
- Screw-cap bottle with mild alcohol solution

Item code	Description	Packaging
NONFB100	6 ml	20 ea/cs

MED-LAP Suction-irrigation Set

For keeping the working field clean and clear of liquids and smoke

This suction-irrigation system enables easy and efficient aspiration or irrigation, in order to keep the abdominal cavity clear of blood, liquids and smoke during laparoscopic procedures. It, therefore, ensures safe surgical intervention.



Item code	Description	Suction-irrigation cannula			Packaging
		Outer diameter	Length	Lateral holes	
OR1200S	Suction-irrigation set with single spike	5 mm	33 cm	8	12/cs
OR1201D	Suction-irrigation set with double spike	5 mm	33 cm	8	12/cs
OR1245S	Suction-irrigation set with single spike	5 mm	45 cm	8	12/cs
ORC45	Suction-irrigation cannula	5 mm	45 cm	8	10/cs
ONG299601CST12	Suction-irrigation set with single spike	10 mm	33 cm	8	5/cs



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US

ONGxxx set is a class IIa sterile medical device. 0459 (G-MED) Rocamed France
ZI de Signes, Allée de Stockholm
83870 Signes, France

ORISxxx are class IIa sterile medical devices. 2460 (DNV Product Assurance AS) OBELIS S.A.
Bd. Général Wahnis, 53
1030 Brussels, Belgium

The Veress needles are class IIa sterile medical devices. 1639 (SGS Belgium NV) Shanghai International Holding Corp. GmbH (Europe)
Erfelstraße 80
20537 Hamburg, Germany

The anti-fog solution is a class IIa sterile medical device. 2797 (BSI, The Netherlands) Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands

OR12xxx and ORC45 are class IIa sterile medical devices. 0197 (TÜV Rheinland) Luxus Lebenswelt GmbH
Kochstraße 1
47877 Willich, Germany

These products are intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2023 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1048_EN 03/2023