



**Protecting  
illumination  
devices in  
operating  
theatres**

Medline offers sterile single-use light handle covers to help prevent the time-consuming and expensive task of reprocessing reusable light handles in the operating theatre. Our light handle covers are easy to assemble and flat-layered packaging is used to help reduce waste.

Designed in different styles to meet your needs, here are some tips to help you choose from amongst our universal, flexible or clip-on light handle covers.

## For use with reusable light handles

### Reusable Light Handle



### Universal Light Handle Cover

Long grips that adapt easily to the light handle



## For use with universal light handle adaptors

### Adapter with Flange



### Flexible Light Handle Cover

Full coverage of the adapter flange to help lower the risk of contamination



### Adapter with Ring



### Clip-On Light Handle Cover



Internal grips that clip firmly onto the handle adapter to help ensure secure placement



Item number	Description	Packaging	Pack component number
16558A	Universal light handle cover	200/case	9053EU-BULK
FLHC1	Flexible light handle cover (1 eaches)	400/case	9100A-EU
FLHC2	Flexible light handle cover (2 eaches)	200/case	N/A
AR-LHCB01	Semi-rigid light handle cover (1 eaches)	100/case	N/A
AR-LHCB02	Semi-rigid light handle cover (2 eaches)	60/case	N/A

**Looking for universal light handle adaptors?** Please contact your Medline account manager or visit us at: [www.medline.eu/uk](http://www.medline.eu/uk)

Light handle covers are class I sterile medical devices intended to be used by healthcare professionals.

FLHC1, FLHC2 and 16558A are certified by BSI.  Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Chateaubriant, France  2797

AR-LHCB01 and AR-LHCB02 are certified by Presafe Denmark A/S.  ArcRoyal, Virginia Road, Kells, Co. Meath, Ireland  0543

We reserve the right to correct errors that may occur within this brochure.

© 2023 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1344\_EN 02/2023