



SKIN HEALTH SOLUTION

MARATHON®

No-Sting, Cyanoacrylate Skin Protectant

Product Overview

Marathon no-sting cyanoacrylate skin protectant is a non-cytotoxic and fast-drying liquid barrier film for the protection of damaged or intact skin. It is applied as a liquid and dries within a minute, adhering to the contours of the skin to form a transparent flexible film. Marathon no-sting cyanoacrylate skin protectant will wear off naturally as the skin regenerates.



Item Details

Product	Item No.	Size	Order Quantity	NHS Supply Chain Code
Marathon Skin Protectant	MSC093001	0.5 ml	5 ea/bx	ELY 560

Indications

Marathon no-sting cyanoacrylate skin protectant is intended to protect intact or damaged skin from the effects of moisture, friction, shear and adhesive stripping.

Contact Details

Customer Service
Tel.: +44 844 334 5237
Email: uk-customerservice@medline.com

Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Fax: +44 844 334 5238
www.medline.eu/uk

FEATURES

- Protects skin:**
- Exposed to irritation from moisture such as urine, faeces, digestive juices, perspiration and wound drainage
 - Exposed to friction and shear from bedding, clothing and shoes
 - Against irritation caused by adhesive products