

# MED-RIGID

Disposable rigid suction canisters



**An accurate  
and practical  
suction system**



# MED-RIGID Suction Canisters

## Discover strong, durable construction and ease of use.

This is a disposable suction system based on a canister/holder combination for the collection of bodily fluids, liquids and aspiration waste. The canister holder is reusable and the rigid canister is for single-patient use. The MED-RIGID system features single-use disposable hard-locking lid attachments and a clear transparent canister with ability to stand alone. This device encourages the proper disposal of infectious liquid medical waste and enhances worker safety.

Special brackets and stands are available for convenient mounting configurations and mobility.



## MED-RIGID Canister Holders

Reusable Canister Holders are used to support the MED-RIGID canister in any slide mount, such as the wall support or mobile trolley. It is built with a pre-attached red vacuum tube for connection to MED-RIGID canister and an ON/OFF switch for vacuum control.

Item number	Description	Packaging
OR159	MED-RIGID Canister Holder, for 1.2-l canister	20/cs
OR139	MED-RIGID Canister Holder, for 2 or 3-l canister	20/cs



## MED-RIGID Canisters

MED-RIGID Canisters are for single-patient use. They are available in 1.2, 2-l, 3-l and even 12-l volumes, each featuring a hard, disposable body and lid.

Item number	Description	Packaging
OR212	MED-RIGID Canister, 1.2-l	40/cs
OR220	MED-RIGID Canister, 2-l	40/cs
OR230	MED-RIGID Canister, 3-l	40/cs
ORLVC	MED-RIGID Canister, 12-l For this specific canister, please use its dedicated metallic trolley (Item number: OR12LTR)	4/cs



## MED-RIGID Canisters with Suction Tubing

MED-RIGID Canisters are available in kits, including a non-sterile ORNEX tubing. ORNEX is a low-memory flexible tubing with ribbed external surface, which increase its resistance.

Item number	Description	Packaging
OR312	MED-RIGID Canister 1.2-l + 5 mm x 1.8-m tubing	30/cs
OR319	MED-RIGID Canister 1.2-l + 6 mm x 3.0-m tubing	30/cs
OR320	MED-RIGID Canister 2-l + 5 mm x 1.8-m tubing	30/cs



## MED-RIGID Tandem Tubing

This is used to facilitate the connection of several canisters and automatic transfer of fluids from one canister to another.

Item number	Description	Packaging
OR542	Tandem Tubing for MED-RIGID Canister	50/cs



# Hardware and Dedicated Equipment

## Disposable specimen traps

The disposable specimen traps are used to collect and retain samples for analysis. They can be affixed to the accessory port of the MED-RIGID canisters and connected 7-mm or 10-mm suction tubing.

Item number	Description	Packaging
OR123	Short specimen sock (9 cm), with male connector (7-mm diameter)	100/cs
OR125	Standard specimen sock (23 cm), with male connector (7-mm diameter)	100/cs
OR124	Large specimen sock (23 cm), with male connector (10-mm diameter)	100/cs
OR84	Specimen cup (10-cm diameter)	20/cs



## Mounting brackets

Item number	Description	Packaging
DYNDSCWL	Plastic wall mount	20/cs
ORCLAMP	Plastic rail clamp	20/cs
FKD0036	Converter for Medline suction canister, compatible with VacSax holder	10/cs
ORCONVS	Converter for Medline suction canister, compatible with Serres rail clamp	20/cs



## Trolleys

These products enable advanced and customised configurations.

Item number	Description	Packaging
DYNDSCC	Metallic canister trolley (adjustable height between 61 and 106 cm, can hold 4 canisters)	1/cs
ORBASE	Canister trolley (50-cm fixed height, can hold 2 canisters)	1/cs
ORSTAND	Canister holder with adhesive base	10/cs
ORSTANDR	Canister holder with roll stand	2/cs



## Reusable tubing

Our reusable manifold tubing enables customised setup by connecting a multiple-canister network to the vacuum source.

Item number	Description	Packaging
DYND2MNF	Manifold tubing for 2 canisters	10/cs
DYND3MNF	Manifold tubing for 3 canisters	10/cs
DYND4MNF	Manifold tubing for 4 canisters	10/cs



**Medline Industries Ltd.**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237

[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

**Medline Ireland Ltd.**  
1 Grant's Row  
Mount Street Lower  
Dublin 2  
Ireland  
Tel.: +353 1 691 73 06

[www.medline.eu/ie](http://www.medline.eu/ie)  
[ie-customerservice@medline.com](mailto:ie-customerservice@medline.com)

FOLLOW US  

ORCLAMP, ORCONVS and ORSTAND are class I non-sterile medical devices intended to be used by healthcare professionals.

 Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France



MED-RIGID Canister Holder, MED-RIGID Canister, MED-RIGID Canister with Suction Tubing, MED-RIGIDD Tandem Tubing, Disposable specimen traps, DZNDSCWL, DYNDSCC, ORBASE, DYNDxMNF are class I non-sterile medical devices intended to be used by healthcare professionals.



Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France



Before use, please consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093, USA. ML546\_EN 11/2023