

# MediGuard Special Nitrile

Non-sterile powder-free nitrile exam gloves



## Benefits and features

- Comfortable fit
- Reliable grip
- Acceptance quality level (AQL) for EN455-1: 1.0
- Force at break (FAB - median value): 7.0 newtons
- PPE-certified Type B
- Tested against chemotherapy drugs

## Recommended for use in:



Labor and Delivery



Environmental Services



Pharmacy/Oncology



Noncritical Care Unit



Phlebotomy Care Unit



Critical Care Unit




+ Textured fingertips

+ Colour: Violet blue


+ Length: min. 240 mm

**Thickness**  
Finger: 0.10 mm  
Palm: 0.08 mm  
Cuff: 0.07 mm

## Innovative and sustainable packaging technology




Flat-layered packaging



Colour-coded sizing



Recyclable materials



## Consistency is key.

We are committed to quality and helping you deliver optimal healthcare, as seen in our:



## We're committed to corporate social responsibility



### Social impact

Our primary manufacturing facility has created and/or supported more than 3,000 jobs directly



### Environmental impact

ISO 14001-certified to assess environmental risks, measure energy savings and help reduce carbon emissions



### WRAP Certification

Promoting labour and human rights in Medline's supply chain in partnership with Worldwide Accredited Production (WRAP)

## Ordering information

| Item number | Size    | Length | Width  | Packaging           |
|-------------|---------|--------|--------|---------------------|
| MGSNXS      | X-small | 245 mm | 74 mm  | 200 ea/bx, 10 bx/cs |
| MGSNS       | Small   | 249 mm | 84 mm  | 200 ea/bx, 10 bx/cs |
| MGSNM       | Medium  | 240 mm | 97 mm  | 200 ea/bx, 10 bx/cs |
| MGSNL       | Large   | 251 mm | 106 mm | 200 ea/bx, 10 bx/cs |
| MGSNXL      | X-large | 252 mm | 113 mm | 180 ea/bx, 10 bx/cs |

## Regulations and quality standards

- Class 1 non-sterile medical device per Medical Device Regulation (EU) 2017/745 and the UK Medical Device Regulation
- Complies with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended
- Compliant with standards EN ISO 21420; EN ISO 374-1; EN ISO 374-2; EN 16523-1; EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Compliant with Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food



Explore how Medline's gloves range, including MediGuard Special Nitrile, makes a difference around the world. Contact your local Medline representative or visit our website: [www.medline.eu/uk](http://www.medline.eu/uk)

This product is a class I non-sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1343\_EN 06/2023