

Medispo Essential

Interventional latex surgical glove



Comfort, dexterity and protection: The new mould has been designed using insight gained from 20 years of customer feedback, in order to better fit clinicians' hands. This new mould, therefore, offers high-level comfort and dexterity.

You may use Medispo Essential as your only glove, or as an overglove when double-gloving, to facilitate perforation detection and help protect you and your patients.

Protection: With computer control and predictive quality analytics enabled by artificial intelligence (AI), each glove undergoes rigorous quality control testing before its release. In addition, each glove is 100% air-leaked and visually inspected.

Skin-friendly formulation: Unique chlorine-free formulation which is tested according to ISO EN 10993.

Easy donning: Skin-friendly polymer coating allows for easy donning, even with damp hands, designed to save you valuable time.

Consistent calibrated grip: The outer surface of Medispo Essential gloves is specially treated to provide a consistent grip. The grip is monitored and controlled during production, and each lot must pass stringent internal grip controls before being released.

Social responsibility: Medispo Essential is manufactured in a WRAP-certified facility.

Worldwide Responsible Accredited Production (WRAP) is an independent, objective, non-profit organisation dedicated to promoting safe, lawful, humane and ethical manufacturing around the world.



Packaging



Decreasing our environmental footprint by using paper-paper peel pouch.



Primary material

- Made from natural rubber latex
 - Contains 50 ug/g or less water-extractable protein
 - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established.
 - **Caution:** This product contains natural rubber latex, which may cause allergic reactions, including anaphylactic response.
- Powder-free to help prevent powder-related complications
- A versatile, sterile, powder-free latex option

Specifications

- Freedom from holes, per EN 455-1: Acceptable quality level (AQL) of 0.65
- Finger thickness: 0.20 ±0.04 mm
- Colour: cream
- Texture: Microtextured
- Cuff design: Tapered and beaded to reduce roll-down

Ordering information

Item code	Size	Packaging
MLES55	5 ½	50 pairs/box, 400 pairs/case
MLES60	6	50 pairs/box, 400 pairs/case
MLES65	6 ½	50 pairs/box, 400 pairs/case
MLES70	7	50 pairs/box, 400 pairs/case
MLES75	7 ½	50 pairs/box, 400 pairs/case
MLES80	8	50 pairs/box, 400 pairs/case
MLES85	8 ½	50 pairs/box, 400 pairs/case
MLES90	9	50 pairs/box, 400 pairs/case

Recommended for use in

- Sterile procedures and general surgery

Regulations and quality standards

- Factory EN ISO 13485-compliant
- Class IIa sterile medical device; compliant with Medical Device Regulation (MDR) (EU) 2017-745
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



Learn more. Contact your Medline account manager or visit: www.medline.eu/uk



These gloves are class IIa sterile medical devices intended for use by healthcare professionals. Before use, consult the instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1530_EN 11/2024



HBM Medical
Unit 7, The Courtyard, Fonthill Business Park,
Fonthill Road, Dublin 22, D22 XA07, Ireland
www.hbmmedical.com
info@hbmmedical.com



2797
MD Class IIa
2777
PPE III Category