


Nasal Oxygen Catheter

For efficient delivery of supplemental oxygen



Designed with cost-effectiveness and versatility in mind, Medline's Nasal Oxygen Catheter is suitable for use across all hospital wards. It facilitates the administration of supplemental oxygen or enhanced oxygen flow to patients in need of respiratory support.

This device is offered in various diameters to cater to various patient requirements.




PVC straight conical step connector for secure connection to an oxygen source

Ordering Information

Reference	OD*	Length	Connector	Packaging
820608	08	40	Straight conical	100/bx
820610	10	40	Straight conical	100/bx
820612	12	40	Straight conical	100/bx
820614	14	40	Straight conical	100/bx

*Outer diameter



Detail A: Polyether foam for secure and comfortable placement



Detail B: Barbed fitting; compatible with various oxygen delivery systems

Learn more. Contact your Medline representative or visit our website www.medline.eu/uk

Medline Nasal Oxygen Catheter is a medical devices of class IIa sterile, intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.



Shanghai International Holding Corporate GmbH (Europe)
Elffesstraße 90
20537 Hamburg, Germany

We reserve the right to correct errors that may occur within this document.

© 2023 Medline Industries LP. is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1416_EN 08/2023