

ORIS Insufflation Tubing Sets

For proper CO₂ insufflation







Benefits and features

The ORIS insufflation tubing sets complete Medline's laparoscopic portfolio.

- Filter and insufflate clean CO₂, allowing clinicians to effectively establish pneumoperitoneum prior to and during laparoscopic procedures
- Offer flexible tubing and a highly secure filter
- Single-use devices


Medline offers a line of items that are fully compatible for use in conjunction with different insufflators.

Image	ORIS Insufflation Tubing Set	Compatibility range
	ORISISO with ISO 15 M/22 F connector	Karl Storz Richard Wolf Stryker ConMed
	ORIS15 with 15-mm male connector	Karl Storz Richard Wolf
	ORISLL with Leur-lock connector	Olympus
	ORISPUSH with flexible female connector	Stryker

Ordering information

Medline's offering provides you with high-quality solutions and best-value laparoscopic essentials.

Item number	Description	Length	Packaging
ORISISO	OR Insufflation Set with ISO 15 M/22 F connector	3 m	25 ea/cs
ORIS15	OR Insufflation Set with 15-mm connector	3 m	25 ea/cs
ORISLL	OR Insufflation Set with Luer-lock connector	3 m	25 ea/cs
ORISPUSH	OR Insufflation Set with female connector	3 m	25 ea/cs



ORISLL
Luer-lock connector



ORISPUSH
Female connector



ORISISO
ISO 15 male/22 female connector

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk

This product is a class IIa sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



LOTUS NL B.V.
Koningin Julianaplein 10, 1e verdieping
2595AA, the Hague, the Netherlands



2460 (DNV Product Assurance AS)

We reserve the right to correct errors that may occur within this brochure.

© 2022 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. ML1315_EN 09/2022.