



OR STANDARD Powdered Latex Surgical Glove

- Made from natural rubber latex
 - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established
 - **Caution:** This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses
- A versatile, sterile, powder latex option
- Space-saving folded packaging design

Item Number	Size	Packaging
MSG4155	5 ½	50 pr/bx, 4 bx/cs
MSG4160	6	50 pr/bx, 4 bx/cs
MSG4165	6 ½	50 pr/bx, 4 bx/cs
MSG4170	7	50 pr/bx, 4 bx/cs
MSG4175	7 ½	50 pr/bx, 4 bx/cs
MSG4180	8	50 pr/bx, 4 bx/cs
MSG4185	8 ½	50 pr/bx, 4 bx/cs
MSG4190	9	50 pr/bx, 4 bx/cs

Recommended for use in:
General surgery

Regulations and quality standards

- Certified to EN ISO 13485
- Medical Device Directive 93/42/EEC
- Certified to "Personal Protective Equipment" council directive 89/686/EEC Category III

SPECIFICATIONS:

- Freedom from holes per EN455: AQL 0.65
- Finger thickness: 0.20 mm
- Length: 285 mm
- Colour: white
- Grip: lightly textured
- Former (mould) design: anatomical, curved to reduce hand fatigue
- Cuff design: tapered and beaded to reduce roll-down

To learn more about our products, please contact your Medline representative or visit us at www.medline.com/uk

