



## OR STANDARD PF Interventional Latex Surgical Glove

- Made from natural rubber latex
  - Contains 10µg/g or less of water extractable protein
  - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established
  - **Caution:** This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses
- Powder-free to avoid powder-related complications
- A versatile, sterile, powder-free latex option
- Space-saving folded packaging design

Item Number	Size	Packaging
MSG5255	5 ½	50 pr/bx, 4 bx/cs
MSG5260	6	50 pr/bx, 4 bx/cs
MSG5265	6 ½	50 pr/bx, 4 bx/cs
MSG5270	7	50 pr/bx, 4 bx/cs
MSG5275	7 ½	50 pr/bx, 4 bx/cs
MSG5280	8	50 pr/bx, 4 bx/cs
MSG5285	8 ½	50 pr/bx, 4 bx/cs
MSG5290	9	50 pr/bx, 4 bx/cs

### Recommended for use in:

General surgery

### Regulations and quality standards

- Certified to EN ISO 13485
- Medical Device Directive 93/42/EEC
- Class IIa sterile Medical Device
- Certified to "Personal Protective Equipment" council directive 89/686/EEC Category III

### SPECIFICATIONS:

- Freedom from holes per EN455: AQL 0.65
- Finger thickness: 0.22 mm
- Length: 291 mm
- Colour: Cream
- Grip: Lightly textured
- Former (Mould) Design: Anatomical, curved to reduce hand fatigue
- Cuff Design: Tapered and beaded to reduce roll-down

To learn more about our products, please contact your Medline representative or visit us at [www.medline.com/uk](http://www.medline.com/uk)

