



**Oral care
solutions**

MEDLINE

MOUTH RINSE

- de Mundreinigung
- fr Bain de bouche
- es Enjuague
- nl Mondwater
- sv Munsköljningsmedel
- it Collutorio
- fi Suuvesi
- pt Enxagüe bucal
- da Mundskyllemiddel
- el Στοματικό διάλυμα
- es Ústni voda
- no Munnvann
- pl plyn do płukania ust

5 FLOZ (145 ml) REF MDS096815E

Individual items to implement step-by-step protocols



DenTips® oral swab untreated

Swabs are used to remove mucus build-up from the patient's mouth. They may be used for oral cleaning in combination with mouthwash, moisturiser or other solutions.

MDS096202EU	Swab individually wrapped	500 ea/case (500 swabs/case)
MDS096206EU	Swabs wrapped in a bag of 20	50 ea/case (1,000 swabs/case)

CE Class I non-sterile products*



DenTips oral swab treated with dentifrice

With a sponge treated with dentifrice, DenTips stimulate oral tissue while lifting debris and mucus build-up from the patient's mouth. The swab can be used alone or in combination with toothpaste or moisturiser. In order to activate the toothpaste, use patient saliva or a small trickle of water.

MDS096502EU	Swab individually wrapped	500 ea/case (500 swabs/case)
MDS096504EU	Swabs wrapped in a bag of 10	100 ea/case (1,000 swabs/case)

CE Class I non-sterile products*



Suction swab treated with dentifrice

This swab can be connected to a standard suction line and controlled by using the thumb port. It includes a sponge treated with dentifrice, and is used for oral cleaning when a toothbrush is not needed. The suction function aids in the removal of secretions, and the angled head assists you in reaching the back of the mouth for a more thorough cleaning.

MDS096525EU	Swab individually wrapped	100 ea/case (100 swabs/case)
-------------	---------------------------	---------------------------------

CE Class I non-sterile products*



Suction toothbrush untreated

This toothbrush can be connected to a standard suction line and controlled by using the thumb port. The super-soft toothbrush with a sponge is used for mechanical debridement of the teeth and gums. The suction function aids in the removal of secretions.

MDS096575EU	Toothbrush individually wrapped	100 ea/case (100 toothbrushes/case)
-------------	---------------------------------	--

CE Class I non-sterile products*



Toothbrushes

The adult and mini toothbrushes with exceptionally soft bristles are gentle on the teeth and gums and provide effective oral cleaning.

MDS096082	Adult toothbrush	72 ea/case
TB001	Mini toothbrush	500 ea/case 100 ea/box

Mouth Moisturiser

Applying Medline's Mouth Moisturiser with DenTips to a patient's lips and oral tissues protects them from drying out.

MDS096803E	Mouth Moisturiser	100 ea/case
------------	-------------------	-------------



Single-tray kits designed for convenience and improved processes

The kit contains every item that you need.



The cup is provided as the perfect place to pour the prescribed solution to perform the oral care.



Individual Oral Care Kit standard

The kit contains:

- 2 DenTips green oral swabs treated with dentifrice
- 1 packet of Medline Mouth Rinse (15 ml)
- 1 packet of Medline Mouth Moisturiser gel (3 g).

MDS096013MEU	VAPprevent Oral Care Kit standard	100 ea/case (100 kits/case)
--------------	-----------------------------------	--------------------------------

CE Class I non-sterile products*



Individual Oral Care Kit with suction swabs

The kit contains:

- 2 suction swabs treated with dentifrice
- 1 DenTips oral swab.

Also available with Mouth Rinse and Mouth Moisturiser

MDS09651300EU	Solution-free VAPprevent Kit with suction swabs	100 ea/case (100 kits/case)
MDS096513MEU	VAPprevent Kit with suction swabs and solution	100 ea/case (100 kits/case)

CE Class I non-sterile products*



Individual Oral Care Kit with suction toothbrush

The kit contains:

- 1 suction toothbrush untreated
- 1 DenTips blue oral swab untreated.

Also available with Mouth Rinse and Mouth Moisturiser

MDS09657200EU	Solution-free VAPprevent Kit with suction toothbrush	100 ea/case (100 kits/case)
MDS096572MEU	VAPprevent Kit with suction toothbrush	100 ea/case (100 kits/case)

CE Class I non-sterile products*

Medline Mouth Rinse ingredients: aqua, glycerine, poloxamer 407, flavour, betaine, potassium sorbate, hydroxyethyl cellulose, sucralose, citric acid
Medline Mouth Moisturiser ingredients: glycerin, sucralose, carbomer, sodium hydroxide, flavor, sorbic acid, hydroxyethyl cellulose, disodium EDTA
The products listed above are intended to be used by healthcare professionals.
Before use, consult instructions and precautions on the corresponding labelling.*

Accessories



Flexible suction tubing connector

The flexible tubing connector with two female ends makes it unnecessary to cut tubing to fit suction tools.

VAPCONNECT	Suction tubing connector	25 ea/box (1,000 ea/case)
------------	--------------------------	------------------------------

CE This product is a class I non-sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.**

24-Hour Kits drive the highest levels of compliance.

Medline's 24-Hour Kits offer all necessary components for a full day of oral care on one patient, and are conveniently packaged in a dispenser box. The two different types of trays are colour-coded to help healthcare workers quickly know which one to use, and also to let others easily check whether proper care has been performed. Furthermore, the dispenser box has three convenient storage solutions thanks to the handy plastic/paper strap, allowing you to hang the kit on the wall or ventilator, or you can place it on the bedside table.

Code	Description	Toothbrush tray	Swab tray
Q400	Q4 24-Hour Oral Care Kit with moisturiser	2x 1 suction toothbrush, 1 suction swab, 1 DenTips, 1 moisturiser	4x 1 suction swab, 1 DenTips, 1 moisturiser
Q800	Q8 24-Hour Oral Care Kit with moisturiser	2x 1 suction toothbrush, 1 suction swab, 2 DenTips, 1 moisturiser	1x 1 suction swab, 2 DenTips, 1 moisturiser


CE Class I non-sterile medical device



Medline Industries Ltd
 3rd Floor
 Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom
 Tel.: +44 844 334 5237
 Fax: +44 844 334 5238
uk.medline.eu | uk-customerservice@medline.com

FOLLOW US  

*  Medline International France SAS
 5, rue Charles Lindbergh
 44110 Châteaubriant, France

 Medline Industries Inc.
 Three Lakes Drive
 Northfield, IL 60093, USA

**  Shanghai International Holding
 Corp. GmbH (Europe)
 Eiffelstraße 80
 20537 Hamburg, Germany

 Conod Medical Co., Ltd.
 No. 38 Hongfeng Road, Baimao Industrial Park, Guli Town
 215532 Changshu City, Jiangsu Province, P.R. China

All VAPrevent products are manufactured by Medline Industries, Inc. These products are medical devices. Before use, please refer to the instructions for any warnings or precautions.

We reserve the right to correct errors that may occur within this brochure.
 © 2020 Medline and DenTips are registered trademarks of Medline Industries, Inc. VAPrevent is a trademark of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060, USA. ML427_EN 12/2020.