



Surgical Gloves Ordering Guide

Latex surgical gloves

Standard portfolio

	Uses	Texture	Thickness	Colour	Sizes	Item numbers	Packaging
Signature Latex Performance and protection	General surgery	Smooth grip	0.22 mm	Tan	5.5 – 9	MSG5455 – MSG5490	50/bx, 200/cs
Signature Latex Micro A higher level of tactile sensitivity	General surgery, neurosurgery, vascular, plastic, ophthalmology	Smooth grip	0.20 mm	Tan	5.5 – 9	MSG5855 – MSG5890	50/bx, 200/cs
Signature Latex Essential Superior donning and flexibility	General surgery, top layer when double gloving	Smooth grip	0.23 mm	Cream/translucent	5.5 – 9	MSG5955 – MSG5990	50/bx, 200/cs
OR Standard PF Latex powder-free	General surgery	Lightly textured grip	0.21 mm	Cream/translucent	5.5 – 9	MSG5255 – MSG5290	50/bx, 200/cs

Underglove





Signature Latex Green Dark green underglove with aloe	Cardiac, trauma, orthopaedics, high-risk procedures	Smooth grip	0.19 mm	Dark green	5.5 – 9	MSG5555 – MSG5590	50/bx, 200/cs
---	---	-------------	---------	------------	---------	-------------------	---------------

Powdered latex



OR Standard Latex powder	General surgery	Lightly textured grip	0.21 mm	White	5.5 – 9	MSG4155 – MSG4190	50/bx, 200/cs
------------------------------------	-----------------	-----------------------	---------	-------	---------	-------------------	---------------

Latex-free surgical gloves

Standard portfolio

	Uses	Material	Texture	Thickness	Colour	Sizes	Item numbers	Packaging
Neolon® 2G Standard polychloroprene 	General surgery, ideal for staff with skin sensitivities and double gloving	Polychloroprene	Lightly textured grip	0.22 mm	Brown	5.5 – 9	MSG6055 – MSG6090	50/bx, 200/cs
SensiCare® PI Evolution PI blend 	General surgery	Polyisoprene	Smooth grip	0.21 mm	Cream/translucent	5.5 – 9	MSG9555 – MSG9590	50/bx, 200/cs
SensiCare® PI The performance of latex without the proteins 	General surgery	Polyisoprene	Smooth grip	0.23 mm	Cream/translucent	5.5 – 9	MSG9055 – MSG9090	50/bx, 200/cs
SensiCare® PI Micro Enhanced tactile sensitivity 	General surgery, neurosurgery, vascular, plastic, ophthalmology	Polyisoprene	Smooth grip	0.20 mm	Cream/translucent	5.5 – 9	MSG9655 – MSG9690	50/bx, 200/cs

Undergloves

DermAssure™ Green Made without chemical accelerators 	General surgery, ideal for skin sensitivities and use with chemicals	Polychloroprene	Smooth grip	0.21 mm	Dark green	5.5 – 9	MSG6555I – MSG6590I	50/bx, 200/cs
SensiCare® PI Green with Aloe Specialty dark green underglove 	Cardiac, trauma, orthopaedics, high-risk procedures	Polyisoprene	Smooth grip	0.21 mm	Dark green	5.5 – 9	MSG9255 – MSG9290	50/bx, 200/cs



Medline Industries Ltd
 3rd Floor
 Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom
 Tel.: +44 844 334 5237
 Fax: +44 844 334 5238

www.medline.eu/uk | uk-customerservice@medline.com

FOLLOW US  

These products are class IIa medical devices. Before use, please consult instructions and precautions on the corresponding labelling.



Medline International France SAS
 5, rue Charles Lindbergh
 44110 Châteaubriant, France



We reserve the right to correct errors that may occur in this brochure.

© 2020 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML225-EN 03/2020.