



# SURGICAL GLOVES ORDERING GUIDE

## LATEX SURGICAL GLOVES

### STANDARD PORTFOLIO

	Uses	Texture	Thickness	Colour	Sizes	Item Numbers	Packaging
<b>Signature Latex</b> Performance and protection	General surgery	Smooth grip	0.23 mm	Tan	5.5 – 9	MSG5455 – MSG5490	50/bx, 200/cs
<b>Signature Latex Micro</b> A higher level of tactile sensitivity	General surgery, neurosurgery, vascular, plastic, ophthalmology	Smooth grip	≥ 0.18 mm	Tan	5.5 – 9	MSG5855 – MSG5890	50/bx, 200/cs
<b>Signature Latex Essential</b> Superior donning and flexibility	General surgery, top layer when double gloving	Smooth grip	0.21 mm	Cream/translucent	5.5 – 9	MSG5955 – MSG5990	50/bx, 200/cs
<b>OR Standard PF</b> Latex powder-free	General surgery	Lightly textured grip	0.22 mm	Cream/translucent	5.5 – 9	MSG5255 – MSG5290	50/bx, 200/cs

### SPECIALTY

<b>Signature Latex Grip</b> Extra texture for a better grip	General surgery, orthopaedics, neurosurgery, ophthalmology	Textured grip	0.24 mm	Translucent	5.5 – 9	MSG2955 – MSG2990	50/bx, 200/cs
<b>Signature Latex Ortho</b> Extra thickness, extra protection	General orthopaedics, trauma, high-risk procedures	Lightly textured grip	0.31 mm	Brown	6.0 – 9	MSG5660 – MSG5690	50/bx, 200/cs
<b>Signature Latex Green</b> Dark green underglove with aloe	Cardiac, trauma, orthopaedics, high-risk procedures	Smooth grip	0.20 mm	Dark green	5.5 – 9	MSG5555 – MSG5590	50/bx, 200/cs

### POWDERED LATEX

<b>OR Standard</b> Latex powder	General surgery	Lightly textured grip	0.20 mm	White	5.5 – 9	MSG4155 – MSG4190	50/bx, 200/cs
------------------------------------	-----------------	-----------------------	---------	-------	---------	-------------------	---------------

# LATEX-FREE SURGICAL GLOVES

## STANDARD PORTFOLIO

	Uses	Material	Texture	Thickness	Colour	Sizes	Item Numbers	Packaging
<b>Neolon® 2G</b> Standard polychloroprene	General surgery, ideal for staff with skin sensitivities and double gloving	Polychloroprene	Lightly textured grip	0.23 mm	Brown	5.5 – 9	MSG6055 – MSG6090	50/bx, 200/cs
<b>SensiCare® PI Evolution</b> PI blend	General surgery	Polyisoprene	Smooth grip	0.23 mm	Cream/translucent	5.5 – 9	MSG9555 – MSG9590	50/bx, 200/cs
<b>SensiCare® PI</b> The performance of latex without the proteins	General surgery	Polyisoprene	Smooth grip	0.23 mm	Cream/translucent	5.5 – 9	MSG9055 – MSG9090	50/bx, 200/cs
<b>SensiCare® PI Micro</b> Enhanced tactile sensitivity	General surgery, neurosurgery, vascular, plastic, ophthalmology	Polyisoprene	Smooth grip	0.20 mm	Cream/translucent	5.5 – 9	MSG9655 – MSG9690	50/bx, 200/cs

## SPECIALTY

<b>DermAssure™ Green</b> Made without chemical accelerators	General surgery, ideal for skin sensitivities and use with chemicals	Polychloroprene	Smooth grip	0.20 mm	Dark green	5.5 – 9	MSG6555I – MSG6590I	50/bx, 200/cs
<b>SensiCare® PI Ortho</b> Extra thickness, extra protection	General orthopaedics, trauma, high-risk procedures	Polyisoprene	Smooth grip	0.31 mm	Cream/translucent	6.0 – 9	MSG9460 – MSG9490	50/bx, 200/cs
<b>SensiCare® PI Green with Aloe</b> Specialty dark green underglove	Cardiac, trauma, orthopaedics, high-risk procedures	Polyisoprene	Smooth grip	0.20 mm	Dark green	5.5 – 9	MSG9255 – MSG9290	50/bx, 200/cs
<b>SensiCare® Shield</b> Radiation attenuation	Fluoroscopic procedures to shield hands from scattered radiation	Polyisoprene and tungsten	Lightly textured grip	0.30 mm	Black	6.5 – 9	MSG3965 – MSG3990	5/bx, 20/cs

These products are class IIa medical devices. Before use, consult instructions for use and precautions on the corresponding labelling.

CE 2797 (BSI)  Medline Industries France SAS  
11 rue de la République  
44100 Nantes Cedex 03, France



**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
Fax: +44 844 334 5238  
[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

**ALWAYS  
ON.**