












# SURGICAL GLOVE ORDERING GUIDE

## LATEX SURGICAL GLOVES

### STANDARD PORTFOLIO

	Uses	Texture	Thickness	Color	Sizes	Ordering Info	Packaging
<b>Signature Latex</b> Performance + Protection 	General surgery	Smooth grip	0.23mm	Tan	5.5-9	MSG5455-MSG5490	50/bx, 200/cs
<b>Signature Latex Micro</b> A Higher Level of Tactile Sensitivity 	General surgery, neuro, vascular, plastic, eye surgery	Smooth grip	0.20mm	Tan	5.5-9	MSG5855-MSG5890	50/bx, 200/cs
<b>Signature Latex Essential</b> Superior Donning & Flexibility 	General surgery, top layer when double gloving	Smooth grip	0.21mm	Cream/translucent	5.5-9	MSG5955-MSG5990	50/bx, 200/cs
<b>Triumph® Classic</b> A Smart Alternative to Powdered Latex 	General surgery	Smooth grip	0.22mm	White	5.5-9	MSG5055-MSG5090	50/bx, 200/cs
<b>OR Standard PF</b> Latex Powder-Free 	General surgery	Lightly textured grip	0.21mm	Cream/translucent	5.5-9	MSG5255-MSG5290	50/bx, 200/cs

### SPECIALTY

<b>Signature Latex Grip</b> Extra Texture For A Better Grip 	General surgery, orthopedics, neurosurgery, ophthalmology	Textured	0.24mm	Translucent	5.5-9	MSG2955-MSG2990	50/bx, 200/cs
<b>Signature Latex Ortho</b> Extra Thickness, Extra Protection 	General orthopedics, trauma, & high-risk procedures	Lightly textured grip	0.31mm	Brown	6.0-9	MSG5660-MSG5690	50/bx, 200/cs
<b>Triumph® Green</b> Dark Green Underglove with Aloe 	Cardiac, trauma, ortho, & high-risk procedures	Smooth grip	0.22mm	Dark Green	5.5-9	MSG2555-MSG2590	50/bx, 200/cs
<b>Signature Latex Green</b> Dark Green Underglove with Aloe 	Cardiac, trauma, ortho, & high-risk procedures	Smooth grip	0.22mm	Dark Green	5.5-9	MSG5555-MSG5590	50/bx, 200/cs

### POWDERED LATEX

<b>OR Standard</b> Latex Powder 	General surgery	Lightly textured grip	0.20mm	White	5.5-9	MSG4155-MSG4190	50/bx, 200/cs
--	-----------------	-----------------------	--------	-------	-------	-----------------	---------------

# LATEX-FREE SURGICAL GLOVES

## STANDARD PORTFOLIO

	Uses	Material	Texture	Thickness	Color	Sizes	Ordering Info	Packaging
<b>Neolon® 2G</b> Standard Polychloroprene	General surgery; ideal for skin sensitivities, chemicals	Polychloroprene	Lightly Textured	0.23mm	Brown	5.5-9	MSG6055-MSG6090	50/bx, 200/cs
<b>SensiCare® PI Evolution</b> PI Blend	General surgery	Polyisoprene	Smooth grip	0.23mm	Cream/translucent	5.5-9	MSG9555-MSG9590	50/bx, 200/cs
<b>SensiCare® PI</b> The Performance of Latex without the Proteins	General surgery	Polyisoprene	Smooth grip	0.23mm	Cream/translucent	5.5-9	MSG9055-MSG9090	50/bx, 200/cs
<b>SensiCare® PI Micro</b> Enhanced Tactile Sensitivity	General surgery, neuro, vascular, plastic, & eye surgery	Polyisoprene	Smooth grip	0.20mm	Cream/translucent	5.5-9	MSG9655-MSG9690	50/bx, 200/cs

## SPECIALTY

<b>DermAssure™ Green</b> Not Made with Chemical Accelerators	General surgery; ideal for staff with skin sensitivities and double gloving	Polychloroprene	Smooth grip	0.20mm	Dark Green	5.5-9	MSG6555I-MSG6590I	50/bx, 200/cs
<b>SensiCare® PI Ortho</b> Extra Thickness, Extra Protection	General orthopedics, trauma, & high-risk procedures	Polyisoprene	Smooth grip	0.31mm	Cream/translucent	6.0-9	MSG9460-MSG9490	50/bx, 200/cs
<b>SensiCare® PI Green with Aloe</b> Specialty Dark Green Underglove	Cardiac, trauma, ortho, & high-risk procedures	Polyisoprene	Smooth grip	0.20mm	Dark Green	5.5-9	MSG9255-MSG9290	50/bx, 200/cs
<b>SensiCare® Shield</b> Radiation Attenuation	Fluoroscopic procedures to shield hands from scattered radiation	Polyisoprene and Tungsten	Lightly Textured	0.30mm	Black	6.5-9	MSG3965-MSG3990	5/bx, 20/cs

This product is a class IIa medical device. Before use, consult instructions for use and precautions on the corresponding labeling.

CE 0086 (BSI)  Medline International France SAS  
11 rue de la République  
91000 Châtenay-Malabry, France



**Medline Industries Ltd**  
No 5 Booths Park  
Chelford Road  
Knutsford, Cheshire WA16 8GS  
Tel: +44 844 334 5237  
Fax: +44 844 334 5238  
[www.medline.com/uk](http://www.medline.com/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)