

Surgical Gloves Ordering Guide



Latex-free surgical gloves



	SensiCare® Neoprene Standard polychloroprene	SensiCare® PI Evolution PI blend	SensiCare® PI The performance of latex without the proteins	SensiCare® PI Micro Enhanced tactile sensitivity	DermAssure™ Green Made without chemical accelerators	SensiCare® PI Green with Aloe Specialty dark green underglove
Item numbers	MSG6355 – MSG6390	MSG9555 – MSG9590	MSG9055 – MSG9090	MSG9655 – MSG9690	MSG6555 – MSG6590	MSG9255 – MSG9290
Sizes	5.5 – 9	5.5 – 9	5.5 – 9	5.5 – 9	5.5 – 9	5.5 – 9
Material	Polychloroprene	Polyisoprene	Polyisoprene	Polyisoprene	Polychloroprene	Polyisoprene
Thickness	0.22 mm	0.21 mm	0.23 mm	0.20 mm	0.21 mm	0.21 mm
Texture	Lightly textured grip	Smooth grip	Smooth grip	Smooth grip	Smooth grip	Smooth grip
Colour	Brown	Cream/translucent	Cream/translucent	Cream/translucent	Dark green	Dark green
Packaging	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs
Uses	General surgery, ideal for staff with skin sensitivities and double gloving	General surgery	General surgery	General surgery, neurosurgery, vascular surgery, plastic surgery, ophthalmology	General surgery, ideal for skin sensitivities and use with chemicals	Cardiac, trauma, orthopaedics, high-risk procedures

Latex surgical gloves



	Signature Latex Performance and protection	Signature Latex Micro A higher level of tactile sensitivity	Signature Latex Essential Superior donning and flexibility	Signature Latex OR Powder-free	Signature Latex Green Dark green underglove with aloe
Item numbers	MSG5455 – MSG5490	MSG5855 – MSG5890	MSG5955 – MSG5990	MSG7255 – MSG7295	MSG5555 – MSG5590
Sizes	5.5 – 9	5.5 – 9	5.5 – 9	5.5 – 9.5	5.5 – 9
Thickness	0.22 mm	0.20 mm	0.23 mm	0.21 mm	0.19 mm
Texture	Smooth grip	Smooth grip	Smooth grip	Lightly textured grip	Smooth grip
Colour	Tan	Tan	Cream/translucent	Cream/translucent	Dark green
Packaging	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs	50/bx, 200/cs
Uses	General surgery	General surgery, neurosurgery, vascular surgery, plastic surgery, ophthalmology	General surgery, top layer when double gloving	General surgery	Cardiac, trauma, orthopaedics, high-risk procedures

Make a
material
difference







Make a material difference



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
uk.medline.eu | uk-customerservice@medline.com

FOLLOW US  

MSG72xx gloves are class IIa medical devices certified by DEKRA Certification B.V. (notified body number 0344).



HBM Medical
Coliemore House, Coliemore Road, Dalkey, Co.
Dublin, Ireland

CE 0344 (DEKRA Certification B.V.)

These products are class IIa medical devices. Before use, please consult instructions and precautions on the corresponding labelling.



Medline International France SAS
5, rue Charles Lindbergh
44110 Châteaubouillant, France

CE 2797 (BSI)

We reserve the right to correct errors that may occur in this brochure.

© 2021 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML225_EN 04/2021.