

QuickSuite® Operating Room Table Kits

All-in-one kits for your operating theatre



**Improve
efficiency.
Promote
safety.**



Reduce operating theatre turnover time. Promote patient safety.

Medline's QuickSuite OR Table Kits contain an assortment of disposable components to be used before a procedure to set up the operating theatre table and after the procedure to rapidly clear it for the next patient. QuickSuite lets you conquer some of your largest theatre challenges by helping you to:

Boost efficiency

QuickSuite OR Table Kits contain everything you need to quickly and efficiently set up the operating theatre table for the next procedure.

Promote safety

QuickSuite items are single-use and should be discarded after the procedure, reducing the risk of cross-contamination. With QuickSuite, you get a collection of products, including our Sahara Absorbent OR Table Sheet, that work together to keep your patient dry and comfortable. The impervious layer of our Sahara sheet protects the theatre table pads since it prevents liquids from touching it.

Lower costs


The different QuickSuite kits have everything you need for setting up the operating theatre table, which helps you eliminate laundering costs, as well as replacement costs for linens and operating theatre table pads.

Reduce moisture: Sahara Absorbent OR Table Sheet.

Excessive moisture is one of the major causes of skin maceration and breaking. Tests show that the Sahara Absorbent OR Table Sheet absorbs 2,056 per cent of its mass or over 3.7 litres/m² of fluid.⁽¹⁾ This helps reduce the risk of skin breakdown through moisture management.

BENEFITS

- Core polymer technology rapidly wicks moisture away, maintaining skin integrity.
- There are 3 different sizes for different usage: 102 x 230 cm, 102 x 152 cm and 100 x 310 cm.
- Sahara is quilted for a smooth surface.
- Disposable design creates a low risk for harbouring bacteria.
- Sahara is low linting.
- Sahara Grippy is ideal for surgical procedures that require a tilted position.



⁽¹⁾Data on file

Choose your stock kits.

Medline offers a range of components and kits to satisfy all your needs.

Article No.	Sahara absorbent table sheet			Draw/lift sheet		Armboard covers 72.5 x 33.5 cm	Fitted armboard cover 57 x 31 cm	Head positioner cover		Patient body strap 83.5 x 10 cm	Packaging (ea/cs)
	102 x 152 cm	102 x 230 cm	100 x 310 cm	102 x 152 cm	152 x 198 cm			38 x 39 cm	45 x 55 cm		
DYND4060SBE	x										20
DYND4090SBE		x									20
QSS01 ⁽²⁾			x								20
DYNDD54060E				x							50
DYNDD56078E					x						48
QSACOVER						x					100
DYKQSHPC1CEA								x			300
QSHPC									x		300
ORS01										x	24
DYKQSUITEL2CEA		x		x							24
DYKQSUITEL2CEASB ⁽²⁾		x		x							24
DYKQSUITEL3CEA		x		x		x					16
DYKQSUITEL4CEA		x		x		x		x			16
DYKQS01	x			x		x					30
DYKQSCO1 ⁽³⁾		x		x			x			x	16

⁽²⁾These items include a slippery back sheet. All other Sahara sheets include a non-slip back sheet

⁽³⁾This item includes 3 disposable underpads.

Choose your compatible accessories.

In addition to our disposable components and kits, Medline offers some other components for patient comfort and modesty during examinations and surgeries, in the perioperative and post-operative periods.



BENEFITS

- Low linting
- Breathable
- Strong and tear-resistant
- Gentle and non-irritating
- Single-use, minimising cross-contamination between patients



Article No.	Description	Size	Packaging (ea/cs)
DORB01 ⁽⁴⁾	Disposable Patient Blanket	110 x 210 cm	30
ORS01	Disposable Body Strap	83.5 x 10 cm	24

⁽⁴⁾This product offers the additional benefit of being able to be heated in the warming cabinet. In addition, it comes individually wrapped in a PE bag.

Advance your operating theatre.



Promote safety.

Lower the risk of cross-contamination with the Sahara Absorbent Theatre Table Sheet.



Improve efficiency.

Reduce the time it takes to turn over your operating theatre with easy, customisable kits.



Reduce costs.

Eliminate waste with custom-made kits and avoid the associated costs of reusable linens.

For more information on our products, please contact your account manager or visit us at: www.medline.eu/uk



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk | uk-customerservice@medline.com

FOLLOW US  

These products (except DORB01) are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

 Medline International France SAS
S, rue Charles Lindbergh
44110 Châteaubriant, France



We reserve the right to correct any errors that may occur in this brochure.

© 2020 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML129-EN 01/2020.