

QuickSuite® and Floor Safety

Moisture management and floor safety product catalogue



**Improve
efficiency.
Promote
safety.**



Medline's moisture management solution

As a global provider of comprehensive moisture management solutions, Medline fulfils the needs of facilities like yours by offering a wide range of high-quality and cost-effective products that fully meet clinical requirements.

Choosing a moisture management system is not only an important step in providing the best and most efficient care possible, but it is also fundamental in ensuring the protection of equipment.

QuickSuite's Sahara operating table sheets are designed to manage moisture and optimise the operating theatre workflow. They help reduce clean-up time by absorbing and retaining fluids during procedures and keeping the patient and equipment dry.

With a variety of QuickSuite accessories, customisable QuickSuite kits and Liqui-Loc™ floor mats, the possibilities are endless for perfectly tailoring your moisture management solution to your needs.

Table of contents

QuickSuite	4
Sahara operating table sheets	4
QuickSuite accessories	8
QuickSuite operating table kits	12
Floor safety	14
Floor safety products	14


A photograph of an operating theatre table. The table is covered with a white sterile drape. A blue sterile drape is draped over the side of the table, and a blue sterile drape is draped over the end of the table. A blue sterile drape is also draped over the arm of the table, which is secured with a white strap. The table is positioned on a blue floor. The background is a plain white wall.

**Designed to
optimise the
operating theatre
workflow**

Sahara operating table sheets

Boost efficiency. Promote safety. Lower costs.

Medline's Sahara operating table sheets have a wicking function paired with an impervious backing, which promotes effective moisture management and the optimisation of the operating table's lifespan.



- Quilted top layer made of non-woven polypropylene to ensure comfort and reduce linting
- Middle absorbent layer reinforced by core polymer technology that rapidly wicks moisture away, maintaining skin integrity
- Impervious bottom layer that keeps all fluids away from the operating table

1**Boost efficiency**

Allows for quick and efficient setup and cleaning of the operating table

2**Increase comfort**

Keeps the patient dry and comfortable

3**Lower costs**

Eliminates laundering costs and replacement costs for linens and operating table pads

Sahara operating table sheets


Choose the solution that suits you best.

Select the Sahara operating table sheet that best suits your needs in terms of absorbency and patient positioning. With a variety of product options, you can choose the most appropriately sized surgical table cover to help you protect the operating table.

Sahara operating table sheets are available in a variety of configurations to meet your needs for different procedures, and come with two types of backings:


Grippy backing

Textured polyethylene that increases the friction between the sheet and the operating table's surface



Standard backing

Smooth polyethylene that allows for increased manoeuvrability



Sahara operating table short sheets

Designed for lithotomy positioning

Sahara super absorbent operating table sheet (short)

Grippy

This shorter sheet has been designed for gynaecological and urological procedures.

102 x 152 cm | Reference: DYND4060SBE



Sahara operating table long sheets

Designed for supine positioning

Sahara super absorbent operating table sheet (long)

Grippy

This sheet is ideal for low- to medium-fluid procedures lasting less than three hours.

102 x 229 cm | Reference: DYND4090SBE



Sahara fitted operating table sheet

Grippy

This fitted sheet has been designed for a better grip on the table, stretcher or gurney.

102 x 229 cm | Reference: DYNDE4090FTD



Sahara protective operating table sheet

Grippy

This non-absorbent sheet can be used for short, low-fluid surgical procedures.

102 x 229 cm | Reference: DYND4090NA



Sahara Essential operating table sheet

Grippy

This eco-friendly sheet provides an ideal solution for less invasive surgical procedures with low to moderate fluid usage.

80 x 220 cm | Reference: SL80220E



Sahara economy operating table sheet

Standard

This sheet is ideal for low- to medium-fluid procedures lasting less than three hours.

102 x 229 cm | Reference: DYND4090E



Item number	Description	Backing	Size	Packaging
DYND4090SBE	Sahara super absorbent operating table sheet (long)	Grippy	102 x 229 cm	20/cs
DYNDE4090FTD	Sahara fitted operating table sheet	Grippy	102 x 229 cm	20/cs
DYND4090NA	Sahara protective operating table sheet	Grippy	102 x 229 cm	40/cs
DYND4060SBE	Sahara super absorbent operating table sheet (short)	Grippy	102 x 152 cm	20/cs
QSS01	Sahara angiography operating table sheet	Standard	102 x 310 cm	40/cs
SL80220E	Sahara Essential operating table sheet	Grippy	80 x 220 cm	10/cs
DYND4090E	Sahara economy operating table sheet	Standard	102 x 229 cm	20/cs

QuickSuite accessories

A variety of products to complete your moisture management solution

In addition to the Sahara operating table sheets, Medline offers other components to complete your moisture management solution and ensure patient comfort and modesty during examinations and surgeries, in the perioperative and post-operative periods.

Operating table protection

Head positioning cover

This cover is available in two sizes and helps protect the operating table headboard from being damaged.

38 x 39 cm | Reference: DYKQSHPC1CEA



45 x 55 cm | Reference: QSHPC



Arm board cover with strap

This cover helps protect operating table arm boards from being damaged, and is available in a standard and fitted version.

Standard | 72.5 x 33.5 cm | Reference: QSACOVER



Fitted | 75 x 39 cm | Reference: QSACEL01



Patient positioning

Patient body strap

This solution ensures safe patient positioning while they are on the operating table.

10 x 83.5 cm | Reference: ORS01



Draw/lift sheets

These lift sheets are available in a variety of weight options, and with or without handles. They are optimal for repositioning and transferring patients to and from the operating table.

≥ 224 kg | Without handles | 152 x 198 cm | Reference: DYNDDS6078E



≥ 200 kg | Without handles | 102 x 152 cm | Reference: DYNDDS4060E



≥ 220 kg | With handles (can be used for overhead lifting) | 90 x 175 cm (including handles) | Reference: QSL501



Item number	Description	Size	Packaging
DYKQSHPC1CEA	Head positioning cover	38 x 39 cm	300/cs
QSHPC	Head positioning cover	45 x 55 cm	300/cs
QSACOVER	Arm board cover with strap	72.5 x 33.5 cm	100/cs
QSACEL01	Fitted arm board cover with strap	75 x 39 cm	100/cs
ORS01	Patient body strap	10 x 83.5 cm	24/cs
DYNDDS6078E	Draw/lift sheet	152 x 198 cm	100/cs
DYNDDS4060E	Draw/lift sheet	102 x 152 cm	100/cs
QSL501	Draw/lift sheet with handles	75 x 39 cm	100/cs

Patient warmth

Disposable patient blankets

Disposable patient blanket, white

This blanket can be used for patient comfort and modesty during examinations and surgeries.

110 x 210 cm | Reference: DORB01



Disposable patient blanket, navy

This lightweight blanket made of recycled fluff helps regulate the patient's temperature and mask stains.

110 x 210 cm | Reference: DORB03



Outer layer

- Polypropylene for flame resistance




Inner layer

- High-quality polyester fluff for breathability

Outer layer

- Polypropylene for flame resistance
- Dark blue colour to mask biological fluids and stains, improving patient dignity



Inner layer

- Recycled high-quality polyester fluff for breathability

Item number	Description	Size	Packaging
DORB01	Disposable patient blanket, white	110 x 210 cm	30/cs
DORB03	Disposable patient blanket, navy	110 x 210 cm	30/cs

Did you know?

Our navy disposable patient blanket is made of 60% recycled polyester.



**Designed to
provide warmth,
comfort and
dignity**

QuickSuite operating table kits

Reduce waste and boost efficiency with our kits.

Medline's QuickSuite operating table kits contain an assortment of disposable components to be used before a procedure to set up the operating theatre table, and afterwards to rapidly clear it for the next patient.

Reference: DYKQSUITEL2CEA, 2 components:

1 x Sahara super absorbent operating table sheet (long)
1 x draw/lift sheet

Grippy



Reference: DYKQSUITEL2CEASB, 2 components:

1 x Sahara economy operating table sheet
1 x draw/lift sheet

Standard



Reference: DYKQSUITEL3CEA, 3 components:

1 x Sahara super absorbent operating table sheet (long)
1 x draw/lift sheet
2 x arm board covers with strap

Grippy



Reference: DYKQSUITEL4CEA, 4 components:

1 x Sahara super absorbent operating table sheet (long)
1 x draw/lift sheet
2 x arm board covers with strap
1 x head positioning cover

Grippy



Gynaecology/urology kits

Reference: DYKQS01, 3 components:

1 x Sahara super absorbent operating table sheet (short)
1 x draw/lift sheet
2 x arm board covers with strap

Grippy



Reference: DYKQSC02, 3 components:

1 x Sahara super absorbent operating table sheet (short)
1 x draw/lift sheet
2 x underpads

Grippy



Item number	Sahara sheet (1x)		Draw/lift sheet (1x)	Arm board covers with strap (2x)	Head positioner cover (1x)	Underpads (3x)	Packaging
	102 x 229 cm	102 x 152 cm	102 x 152 cm	72.5 x 33.5 cm	38 x 39 cm	60 x 90 cm	
DYKQSUITEL2CEA	x		x				24/cs
DYKQSUITEL2CEASB	x		x				24/cs
DYKQSUITEL3CEA	x		x	x			16/cs
DYKQSUITEL4CEA	x		x	x	x		16/cs
DYKQS01		x	x	x			30/cs
DYKQSC02		x	x			x	16/cs

QuickSuite custom kits


Create your own kit.

You can create your own customised QuickSuite kit by choosing one type of Sahara operating table sheet and combining it with as many QuickSuite accessories as you would like. Are you interested in receiving more information or a quote for your tailored solution? Contact us!

Sahara operating table sheet (choose 1)



Sahara super absorbent operating table sheet (long)
DYND4090SBE



Sahara super absorbent operating table sheet (short)
DYND4060SBE




Sahara fitted operating table sheet
DYNDE4090FTD



Sahara protective operating table sheet
DYND4090NA




QuickSuite accessories (choose as many as you like)



Head positioning cover
DYKQSHPC1CEA (38 x 39 cm)
QSHPC (45 x 55 cm)



Arm board covers with strap
QSACOVER (slip-on)
QSACEL01 (fitted)



Body straps
ORS01



Lift sheets
DYND4060E (short)
DYND4078E (long)
QSL01 (with handles)

A photograph showing the lower half of a person wearing blue medical scrubs and blue slip-on shoes. They are standing on a white, wrinkled cloth that is laid out on a light-colored floor. To the right, the edge of a medical table with a blue drape is visible. The scene is brightly lit, suggesting an operating theatre environment.

**Designed to
optimise the
operating theatre
workflow**

Floor safety products

Working in wet environments means being exposed to hazards, including slips and falls. Wet floors can compromise care and add additional challenges for healthcare facilities and staff.

Risk environments:



Operating theatre
and lab stations



Examination
rooms



Scrub sink
areas




Maternity
wards



Central sterilisation
areas

Liqui-Loc™ floor mat

- **Impervious anti-slip back**
- **Drip-safe**, locking in fluids to avoid cross-contamination and complicated handling
- **Easy disposal** for hygiene and safety
- For:
 - Standing and placing tools on
 - Placement in staff-patient proximity areas (operating table, examination room, scrub sinks)
 - Use where large volumes of fluid are generated (gynaecology, urology, arthroscopy)



Item number	Description	Size	Absorbency	Anti-slip	Packaging
ORLLPBS	Liqui-Loc - small	71.1 x 101.6 cm	2.2 l	Yes	13/cs
ORLLPBM	Liqui-Loc - medium	71.1 x 142.2 cm	3.0 l	Yes	10/cs
ORLLPBL	Liqui-Loc - large	71.1 x 182.8 cm	3.9 l	Yes	8/cs

Double-sided floor pad

- **Absorbs on both sides** to contain a higher volume of fluids
- **Is drip-safe**, locking in fluids to avoid cross-contamination and ease disposal
- **Speeds up the cleaning process** and absorbs fluid spills on the floor
- **Improves the overall level of hygiene**, which promotes cost and time savings
- **Do not stand on**



Item number	Description	Size	Absorbency	Anti-slip	Packaging
ORDSFMAT	Double-sided floor pad	75 x 38 cm	7 l	No	50/cs

Advance your operating theatre.



Promote safety

Lower the risk of cross-contamination with the Sahara absorbent theatre table sheet.



Improve efficiency

Reduce the time it takes to turn over your operating theatre with easy, customisable kits.



Reduce costs

Eliminate waste with custom-made kits and avoid the associated costs of reusable linens.

For more information on our products, please contact your account manager or visit us at: www.medline.eu/uk



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06
www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

These products (except DORB01, DPRB03, ORLLPBL, ORLLPBM, ORLLPBS) are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct any errors that may occur in this brochure.

© 2025 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML129_EN 04/2025

 Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France

