

QuickSuite® and Ultrasorbs®

Catalogue of solutions



**Operating
theatre
efficiency
and protection**

QuickSuite helps your operating theatre run smarter, smoother and safer.

Every minute counts.

Hospital operating theatres are a critical but costly resource in the healthcare process. The efficiency and utilisation of operating theatres impact hospitals' finances, clinical effectiveness and general patient outcomes.¹

Keep your procedures on track with QuickSuite and Ultrasorbs for operating theatres for a faster, more efficient turnover.

There's a real need for protection.

The incidence of pressure ulcers in surgical patients continues to increase.³ Apart from pressure, moisture, temperature, friction and shear are four of the main extrinsic factors associated with pressure ulcer development. Keeping patients' skin dry is essential to help prevent pressure ulcers and moisture-associated skin damage (MASD).

Choose from a selection of disposable linens, kits, or individual components designed to wick moisture away from the skin.

Prevent to save.

The role of surface contamination in the transmission of nosocomial pathogens is increasingly being recognised. The inanimate environment in operating theatres (e.g. tables, floors, and equipment surfaces) has been considered a potential source of pathogens that may cause surgical site infections (SSIs).⁵

Our disposable solutions help reduce the risk of cross-contamination and maximise barrier protection during theatre turnover.



Increase efficiency

The cost of inefficiency, due to turnaround delays, is estimated to be **£24.77** per minute.²



Improve safety

Approximately **23%** of hospital-acquired pressure ulcers develop during surgeries lasting longer than 3 hours.⁴



Contain costs

The cost of SSIs to the National Health Service is an additional **£10,523** per patient.⁶

All levels of protection

Elevate patient safety and operational efficiency in the operating theatre with Medline's Sahara and Ultrasorbs absorbent table covers.

Choose the most optimal solution for your surgical procedures:

- **Sahara OR Table Sheets** for standard procedures with a low to medium level of fluids
- **Ultrasorbs Breathable Superabsorbent OR Table Sheet** for high-fluid procedures (cardiac, general/thoracic, orthopaedic); offers enhanced moisture management and maintains an optimal microenvironment for the skin

Both solutions have an impervious design that protects operating tables from fluids, thereby reducing risk of harbouring bacteria.

Sahara Protective OR Table Sheet

This sheet offers effective protection for the operating table during short low-fluid procedures. It has a durable impervious back sheet with a soft non-woven top sheet.



Sahara Superabsorbent OR Table Sheet

This superabsorbent sheet helps to elevate patient comfort and safety during standard procedures in the operating theatre by simultaneously protecting patients' skin and the operating table from fluids.

It is also available in a fitted configuration that guarantees fast and easy application.



Ultrasorbs Breathable Drypads

This ultra-soft, multi-layer breathable disposable sheet offers a specially designed solution for rapid and effective moisture management. Its superabsorbent polymer core wicks moisture away from patients' skin and is dry to the touch within minutes.

Recommended for patients at a high risk of skin breakdown, the sheet is available with a reinforced back sheet for patient turning and repositioning (150 to 170 kg weight capacity). A sterile version is available.




Choose the most suitable solution for your operating theatre

A comparison of Sahara and Ultrasorbs

Sahara

- Core polymer technology to rapidly wick moisture away; maintains skin integrity*
- Quilted for a smooth surface; promotes patient comfort*
- Available in 3 different sizes for different procedures:
 - » 102 x 230 cm
 - » 102 x 152 cm
 - » 100 x 310 cm
- Sahara non-slippery back film for enhanced grip to the operating table surface
- Low risk of harbouring bacteria because it is disposable

(*) Except Sahara Protective



Ultrasorbs

- Air-permeable, superabsorbent core technology to rapidly wick moisture away; reduces body heat build-up and maintains an optimal microenvironment for skin integrity
- Soft top sheet for enhanced patient comfort to help reduce friction and shear
- Available in sterile versions and with a reinforced back sheet; allows you to turn, reposition and transfer patients weighing up to 170 kg
- Enhanced absorbency
- Low risk of harbouring bacteria because it is disposable

Below is the ordering information for the different configurations of the above-mentioned products.

Product	Item code	Size	Absorbency level	Low linting	Air-permeable	Non-slippery back sheet	Lifting capacity	Sterile	Packaging
Sahara®	DYND4090NA	102 x 230 cm	Low	X		X			40
	DYND4060SBE	102 x 152 cm	High	X		X			20
	DYND4090SBE	102 x 230 cm	High	X		X			20
	DYNDE4090FTD	102 x 230 cm	High	X		X			20
	QSS01 (Angio)	100 x 310 cm	High	X					20
Ultrasorbs®	MSCE4090PUPS	101 x 228 cm	Very high	X	X				30
	MSCE4090PPSTA	101 x 228 cm	Very high	X	X			X	25
	USAP4090ES	101 x 228 cm	Very high	X	X		170 kg		25
Extrasorbs®	EXTRASRBE150	80 x 210 cm	Very high	X	X		150 kg		30



Also looking to protect your floors? Ask your Medline account manager about Liqui-Loc™ floor pad solutions.

Barrier protection

QuickSuite linen kits shield your operating table pads from bodily fluids and bacteria, and reduce replacement costs.



Pick the perfect kit ...

Time is money—and safety is everything. Speed up your set-up and take-down processes with a QuickSuite complete theatre turnover kit.

QuickSuite Kit

Each set is single-use, so every case starts clean. And everything is at your fingertips in one convenient kit. You can order a QuickSuite kit with different components for various surgical procedures and depending on your needs. The ability to customise your kit helps reduce waste.

Items	Sahara Absorbent OR Table Sheet		Draw/lift sheet	Armboard cover	Headboard cover	Fitted armboard cover	Disposable body strap	Underpads (2 eaches)	Packaging
	102 x 152 cm	102 x 230 cm	102 x 152 cm	72.5 x 33.5 cm	38 x 39 cm	75 x 39 cm	83.5 x 10 cm	60 x 90 cm	
DYKQSUITEL2CEA		X	X						24
DYKQSUITEL2CEASB*		X	X						24
DYKQSUITEL3CEA		X	X	X					16
DYKQSUITEL4CEA		X	X	X	X				16
DYKQS01	X		X	X					30
DYKQSC01		X	X			X	X		16
DYKQSC02	X		X					X	16

... Or order individual items


Shop ‘off-the-shelf’ to buy only those disposable items you need separately.

Sahara OR Table Sheet

The Sahara sheet contains core polymer technology that rapidly wicks moisture away and helps maintain skin integrity. Its quilted pattern helps moisture to spread evenly throughout the sheet and, thanks to sealed edges, the liquid is prevented from leaking. In addition, the impervious layer and disposable design create low risk for harbouring bacteria.

Item code	Description	Size	Packaging
DYND4060SBE	Sahara Superabsorbent OR Table Sheet	102 x 152 cm	20
DYND4090SBE	Sahara Superabsorbent OR Table Sheet	102 x 230 cm	20
QSS01(*)	Sahara Absorbent Angiography Table Sheet	100 x 310 cm	20
DYND4090NA	Sahara Protective OR Table Sheet	102 x 230 cm	40
DYNDE4090FTD	Sahara Absorbent Fitted Sheet	102 x 230 cm	20

(*)These items include a standard back sheet. All other Sahara sheets include a non-slip back sheet



Individual components

These products are available as single items or collectively in a package in Medline’s QuickSuite OR Table Kit. Every kit contains an assortment of disposable components that may be used to easily and efficiently set up and clean the operating table.

Draw/Lift Sheets

Medline’s soft disposable lift sheets are optimal for repositioning patients and transferring them to and from the operating table. Configurations with handles are also available; designed for overhead lifting systems.

Item code	Description	Size	Packaging	Weight capacity
DYND4060E	Draw/lift sheet; white	102 x 152 cm	20	Up to 200 kg
DYND4078E	Draw/lift sheet; white	152 x 198 cm	20	Up to 224 kg
QSL501	Draw/lift sheet with handles	90 x 175 cm	20	Up to 220 kg



Headboard, Armboard and Leg Stirrups Covers

Designed to protect the operating table headboard, armboards and leg stirrups during surgeries. The impervious polyethylene layer keeps the equipment clean and dry, reducing the potential for cross-contamination and extending the period of usability.

Item code	Description	Size	Packaging
QSACOVER	Armboard cover with straps	72.5 x 33.5 cm	100
QSACEL01	Fitted armboard cover with straps	152 x 198 cm	100
QSLCEL01	Fitted stirrups cover with straps	102 x 52 cm	100
DYKQSHPC1CEA	Operating table head positioner cover	38 x 39 cm	300



Disposable Patient Blanket

Medline’s Disposable Patient Blanket provides warmth, comfort and dignity to patients during examinations and surgery, in the perioperative and postoperative periods.

Item code	Description	Size	Packaging
DORB01	Disposable Patient Blanket	110 x 210 cm	30
DORB02	Disposable Patient Blanket (bulk packaging)	110 x 210 cm	26



Patient Body Straps

These are soft and strong foam straps for securing the patient to the operating table or trolleys and limiting unwanted movements and injuries during surgeries.


Item code	Description	Size	Packaging
ORS01	Disposable body strap; hook-and-loop closure, x2	10 x 83.5 cm	24



Floor Safety Solutions

Our complete floor safety portfolio addresses the hazard of falls, and the need for quick fluid removal and to improve safety in the work environment.

Item code	Description	Size	Absorbency level	Anti-slip	Packaging
ORLLPBS	Liqui-Loc; small	71.1 x 101.6 cm	2.2 ltr	Yes	13/cs
ORLLPBM	Liqui-Loc; medium	71.1 x 142.2 cm	3.0 ltr	Yes	10/cs
ORLLPBL	Liqui-Loc; large	71.1 x 182.8 cm	3.9 ltr	Yes	8/cs
ORDSFMAT	Double-sided floor pad	75 x 38 cm	7.0 ltr	No	50/cs



Looking to protect your floor as well?

Support for success

We understand that introducing new protocols or products can be a significant investment in time and effort. We can help. Our team of experts offers you clinical and product support to help ensure successful implementation.

Our team of experts can provide you with:

- Product-related services and training
- Specialised information on product usage and best practices
- A comprehensive skin management programme.

**Find out more about our clinical and educational support services.
Contact your Medline representative today.**

For more information, please contact your Medline account manager or visit our website: uk.medline.eu




Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
uk.medline.eu
uk-customerservice@medline.com

FOLLOW US  

1. Meredith, J. O. et al. (2011) "Are we operating effectively? A lean analysis of operating theatre changeovers," Operations management research, 4(3-4), pp. 89-98.
2. Ang, W. W. et al. (2016) "The cost of trauma operating theatre inefficiency," Annals of medicine and surgery (2012), 7, pp. 24-29.
3. Chen, Hong-Lin, Chen X., Wu J. (2012). The Incidence of Pressure Ulcers in Surgical Patients of the Last 5 Years: A Systematic Review. WOUNDS, 24 (9), 234-241.
4. Primiano, M., Friend, M., McClure, C., Nardi, S., Fix, L., Schafer, M., Savochka, K., McNett, M. (2011). Pressure ulcer prevalence and risk factors during prolonged surgical procedures. AORN journal, 94(6), 555-66.
5. Donskey, C. J. (2013) "Does improving surface cleaning and disinfection reduce health care-associated infections?," American journal of infection control, 41(5 Suppl), pp. S12-9.
6. McDonald, S. and Bott, A. (2020) "Surgical site infections," Surgery, 38(3), pp. 150-154.

QuickSuite kits, Sahara sheets, draw/lift sheets, headboard/armboard, leg stirrup covers and patient body straps are class I medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Ultrisorbs and Extrasorbs are class I medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



We reserve the right to correct errors that may occur within this brochure.

© 2020 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML1087_EN 12/2020