

# Medline's Safety Lancet

A solution for capillary blood sampling



**Quick and  
simple to  
use**



# How to Use Medline's Safety Lancet



## Advantages

- All needles are gamma-sterilised and protected by a sterility cap
- Fully enclosed needle housing ensures the needle is hidden before and after use to prevent sharps injuries
- Self-destructing mechanism prevents reuse of the device
- Silicone-coated, ultra-sharp needles, ideally positioned during skin penetration
- Various, colour-coded product versions to meet specific blood sample applications and address the variety of patient skin types
- Wide range of designs provides simple, comfortable and secure puncture
- Easy to handle, intuitive activation in only two steps, with no pre-loading, for both push button and contact activations
- Wide and long safety cap for easy removal; simply twist and pull to remove the needle cover
- Precise design for convenient finger positioning and ideal blood sample collection



Reference	Description	Needle	Penetration depth	Colour	Packaging
SLN2118	Normal flow	Needle; 21 G	1.8 mm	Green	200/cs
SLN1818	High flow	Needle; 18 G	1.8 mm	Orange	200/cs
SLN2816	Micro flow	Needle; 28 G	1.6 mm	Blue	200/cs
SLBLADE12	Pediatric	Blade	1.2 mm	Purple	200/cs

**For more information** on our products, please contact your Medline account manager or visit us at: [www.medline.eu/uk](http://www.medline.eu/uk)

These products are class IIa medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Shanghai International Holding Corp.GmbH (Europe)  
Eiffelstraße 90  
20537 Hamburg, Germany



All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML421\_EN 01/2023