

SensiCare[®] Extended Cuff



Longer cuff for better protection.

Choose our 280-mm non-sterile powder-free nitrile exam gloves.

Benefits and features

- Extended cuff intended to protect you longer
- Designed to offer a comfortable fit and reliable grip
- Acceptable quality level (AQL) for EN455-1: 1.0
- Force at break (FAB) median value: 12.4 newtons
- Personal protective equipment (PPE)-certified: Type B
- Tested against chemotherapy drugs*


Ordering information

Item number	Size	Packaging
MDSEXTXS	X-small	100 ea/bx, 10 bx/cs
MDSEXTS	Small	100 ea/bx, 10 bx/cs
MDSEXTM	Medium	100 ea/bx, 10 bx/cs
MDSEXTL	Large	100 ea/bx, 10 bx/cs
MDSEXTXL	X-large	90 ea/bx, 10 bx/cs




Sustainable packaging technology

MicroStop™ film



Flat-layered packing



Colour-coded sizing



100% recyclable materials



Recommended for use in:



Labour and delivery



Sterile processing



Laboratory



Pharmacy/
oncology



Non-critical
care unit



Critical care
unit

Regulations and quality standards

- Class I non-sterile medical device, per Medical Device Regulation (EU) 2017/745 and the UK Medical Device Regulation
- Compliant with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended
- Compliant with standards EN ISO 21420; EN ISO 374-1; EN ISO 374-2; EN 16523-1; EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Compliant with Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food



CE 2797 BSI Group The Netherlands B.V.
1066 EP, Amsterdam, Netherlands



Explore how Medline's gloves range, including SensiCare Extended Cuff, makes a difference around the world. Contact your local Medline representative or visit our website: www.medline.eu/uk

*Not for use with Carmustine or Thiotepa.

This product is a class I non-sterile medical device intended for use by healthcare professionals. Before use, consult the instructions and precautions on the corresponding labeling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML153_EN 06/2024

Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France