

SensiCare® Neoprene

Sterile powder-free synthetic polychloroprene surgical gloves



- Made from second-generation synthetic polychloroprene
 - **Note:** This product is not made from natural rubber latex
 - Softer and more flexible than standard synthetic polychloroprene gloves
- Powder-free and coated with a synthetic polymer to assist with damp hand donning
- Tested for chemotherapy agent permeation, in accordance with ASTM D 6978-05
- Cost-effective, latex-free alternative to synthetic polyisoprene
- Space-saving folded packaging design
- Ideal for general surgery

Specifications

- Freedom from holes per EN455-1: acceptable quality level (AQL) of 0.65
- Finger thickness: 0.19 mm
- Glove length (per EN 455-2, size 7.5): 293 mm
- Colour: light brown
- Texture: smooth
- Former (mould) design: anatomical
- Cuff design: beaded, tapered and textured to reduce roll down

Ordering information

Item code	Size	Packaging
MSG6355	5 ½	50 pairs/box, 4 boxes/case
MSG6360	6	50 pairs/box, 4 boxes/case
MSG6365	6 ½	50 pairs/box, 4 boxes/case
MSG6370	7	50 pairs/box, 4 boxes/case
MSG6375	7 ½	50 pairs/box, 4 boxes/case
MSG6380	8	50 pairs/box, 4 boxes/case
MSG6385	8 ½	50 pairs/box, 4 boxes/case
MSG6390	9	50 pairs/box, 4 boxes/case

Recommended for use in

- General surgery

Regulations and quality standards

- ISO 13485-compliant
- Class IIa medical device; compliant with Medical Device Directive 93/42/EEC
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



These gloves are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML214_EN 11/2023

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk