

# SensiCare® PI Evolution



Made from a synthetic polyisoprene blend rather than natural rubber latex

SensiCare PI Evolution is the newest addition to Medline's revolutionary line-up of surgical gloves without natural rubber latex. The unique polyisoprene blend ensures protection for surgeons and staff who require an alternative to natural rubber latex.

As a result, SensiCare surgical gloves are softer, more elastic, and more comfortable than other latex alternatives, to satisfy clinical needs and support safety initiatives.

## Specifications

- Synthetic polyisoprene blend (not made with natural rubber latex)
- Synthetic polymer coating
- Finger thickness: 0.20 mm
- Length: 285 mm
- Colour: cream
- Grip: smooth

## Ordering information

Item code	Size	Packaging
MSG9555	5 ½	50 pairs/box, 200 pairs/case
MSG9560	6	50 pairs/box, 200 pairs/case
MSG9565	6 ½	50 pairs/box, 200 pairs/case
MSG9570	7	50 pairs/box, 200 pairs/case
MSG9575	7 ½	50 pairs/box, 200 pairs/case
MSG9580	8	50 pairs/box, 200 pairs/case
MSG9585	8 ½	50 pairs/box, 200 pairs/case
MSG9590	9	50 pairs/box, 200 pairs/case

## Recommended for use in

- General surgery

## Regulations and quality standards

- EN ISO 13485 compliant
- Class IIa sterile medical device; compliant with Medical Device Regulation (EU) 2017/745
- Category III personal protective equipment; compliant with PPE Regulation (EU) 2016/425



This glove is a class IIa sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML432\_EN 01/2021

Learn more. Contact your Medline account manager or visit: [uk.medline.eu](http://uk.medline.eu)