

SensiCare® PI Evolution



Made from a synthetic polyisoprene, blend rather than natural rubber latex

SensiCare® PI Evolution is the newest addition to Medline's revolutionary line-up of surgical gloves, made without natural rubber latex. The unique polyisoprene blend is designed to offer protection for surgeons and staff who require an alternative to natural rubber latex.

As a result, SensiCare® surgical gloves are designed to be softer, more elastic, and more comfortable than other latex alternatives, to help satisfy clinical needs and support safety initiatives.

Specifications

- Freedom from holes, per EN455-1: acceptable quality level (AQL) of 0.65
- Synthetic polyisoprene blend (not made with natural rubber latex)
- Synthetic polymer coating
- Finger thickness: 0.20 mm
- Glove length (per EN 455-2, size 7.5): 283 mm
- Colour: cream
- Texture: smooth

Ordering information

Item code	Size	Packaging
MSG9555	5 ½	50 pairs/box, 200 pairs/case
MSG9560	6	50 pairs/box, 200 pairs/case
MSG9565	6 ½	50 pairs/box, 200 pairs/case
MSG9570	7	50 pairs/box, 200 pairs/case
MSG9575	7 ½	50 pairs/box, 200 pairs/case
MSG9580	8	50 pairs/box, 200 pairs/case
MSG9585	8 ½	50 pairs/box, 200 pairs/case
MSG9590	9	50 pairs/box, 200 pairs/case

Recommended for use in

- General surgery

Regulations and quality standards

- EN ISO 13485-compliant
- Class IIa sterile medical device; compliant with Medical Device Directive 93/42/EEC
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France

This glove is a class IIa sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML432_EN 11/2023

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk