



# SENSICARE® PI GREEN WITH ALOE VERA

## STERILE POLYISOPRENE POWDER-FREE SURGICAL GLOVES

- Unique synthetic Isolex polyisoprene formulation
  - Chemical structure identical to natural rubber latex, without the proteins
  - Offers healthcare workers an alternative solution that performs like latex in terms of comfort, fit and feel
- Powder-free helps to reduce powder-related skin irritation and post-operative complications associated with granulomas
- Superior tactile sensitivity
- Aloetouch® Moisturising Technology
  - The world's first therapeutic medical glove technology – clinically proven to moisturise skin
  - Contains 100% pure freeze-dried aloe vera gel on the inner surface
  - Assists with donning in addition to moisturising hands
- Dark green colour allows for rapid detection of perforations in the outer glove
- Ideal for double gloving in cardiac, trauma, ortho and high-risk procedures

Item Number	Size	Packaging
MSG9255	5 ½	50 Pairs/Box, 200 Pairs/Case
MSG9260	6	50 Pairs/Box, 200 Pairs/Case
MSG9265	6 ½	50 Pairs/Box, 200 Pairs/Case
MSG9270	7	50 Pairs/Box, 200 Pairs/Case
MSG9275	7 ½	50 Pairs/Box, 200 Pairs/Case
MSG9280	8	50 Pairs/Box, 200 Pairs/Case
MSG9285	8 ½	50 Pairs/Box, 200 Pairs/Case
MSG9290	9	50 Pairs/Box, 200 Pairs/Case

### Recommended for use in:

Ideal for double gloving in cardiac, trauma, ortho and high-risk procedures

### Regulations and quality standards

- Certified to EN ISO 13485
- Medical Device Directive 93/42/EEC - Class IIA
- Personal Protective Equipment Directive 89/686/EEC - Category III

### SPECIFICATIONS:

- Freedom from Holes per EN 455: AQL 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- No residual of accelerators
- Cuff design: tapered and beaded to reduce roll-down
- Finger tip thickness: 0.20 mm
- Glove length (per EN 445-2, size 7.5): 305 mm
- Colour: green
- Grip: smooth



To learn more about our products, please contact your Medline representative or visit us at [www.medline.eu/uk](http://www.medline.eu/uk)