

SensiCare® PI Green with Aloe Vera

Sterile polyisoprene powder-free surgical gloves



- Unique synthetic Isolex polyisoprene formulation
 - Chemical structure identical to natural rubber latex, without proteins
 - Offers healthcare workers an alternative solution that performs like latex, in terms of comfort, fit and feel
- Powder-free to help reduce powder-related skin irritation and post-operative complications associated with granulomas
- Superior tactile sensitivity
- **AloeTouch® Moisturising Technology**
 - The world's first therapeutic medical glove technology—clinically proven to moisturise skin
 - Contains 100% pure freeze-dried aloe vera gel on the inner surface
 - Assists with donning, in addition to moisturising hands
- Dark-green colour allows for rapid detection of perforations in the outer glove

Specifications

- Freedom from holes, per EN455-1: acceptable quality level (AQL) of 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- No residue from accelerators
- Cuff design: tapered and beaded to help reduce roll-down
- Finger tip thickness: 0.21 mm
- Glove length (per EN 455-2, size 7.5): 301 mm
- Colour: dark green
- Texture: smooth

Ordering information


Item code	Size	Packaging
MSG9255	5 ½	50 pairs/box, 200 pairs/case
MSG9260	6	50 pairs/box, 200 pairs/case
MSG9265	6 ½	50 pairs/box, 200 pairs/case
MSG9270	7	50 pairs/box, 200 pairs/case
MSG9275	7 ½	50 pairs/box, 200 pairs/case
MSG9280	8	50 pairs/box, 200 pairs/case
MSG9285	8 ½	50 pairs/box, 200 pairs/case
MSG9290	9	50 pairs/box, 200 pairs/case

Recommended for use in

- Ideal for double-gloving in cardiac, trauma, orthopaedic and high-risk procedures

Regulations and quality standards

- Factory EN ISO 13485-compliant
- Class IIa sterile medical device; compliant with Medical Device Directive 93/42/EEC
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk

These gloves are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML570_EN 11/2023