

SensiCare® Pro Nitrile Gloves

Heavy-duty nitrile option with exceptional grip and 400-mm cuffs

A close-up photograph of a person's hands wearing blue nitrile gloves. The person is using surgical forceps to handle a surgical instrument. The gloves are long, extending up the forearm, and have a textured surface for grip. The background is a solid blue color.

**Longer gloves
for extra
protection and
enhanced
performance**

Extra protection and better grip

Medline's SensiCare Pro Nitrile Gloves provide exceptional protection, fit and performance. They are designed specially for staff who use harsh chemicals for general tasks and procedures, or clean, decontaminate and/or disinfect surgical instruments.

- All-around protection and performance
- Thick gloves with longer 400-mm cuffs; good performance in puncture resistance tests*; chemo-tested** and fully textured for a better grip

Ordering information

Item number	Size	Packaging
PRO400XS	X-small	50 ea/bx, 10 bx/cs
PRO400S	Small	50 ea/bx, 10 bx/cs
PRO400M	Medium	50 ea/bx, 10 bx/cs
PRO400L	Large	50 ea/bx, 10 bx/cs
PRO400XL	X-large	50 ea/bx, 10 bx/cs

Chemical code letter	Chemical	Permeation level of protection	Degradation (%)
K	Sodium hydroxide (40%)	4	-18.1
P	Hydrogen peroxide (30%)	4	+1.9
T	Formaldehyde (37%)	4	+5.8

EN ISO 374-5:2016	Level
Protection against bacteria and fungi	Pass
Protection against viruses	Pass



Ideal for specialty areas, including:



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238

uk.medline.eu
uk-customerservice@medline.com

FOLLOW US

*Test data is available upon request. **Do not use with Carmustine or Thiotepea.

SensiCare Pro gloves are class I non-sterile medical devices, intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. We reserve the right to correct errors that may occur within this document.

Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant, France