

SensiCare® Pro

Heavy-duty nitrile option with exceptional grip and 400-mm cuffs



**Extra-long
gloves for
enhanced
protection and
performance**

Extra protection and better grip

SensiCare® Pro nitrile exam gloves are designed to provide exceptional protection, fit, and performance. They are tailored for staff who use chemicals for general tasks and procedures, or clean, decontaminate and/or disinfect surgical instruments.

- All-around performance and protection for the hand and forearm
- Fully textured thick gloves with longer 400-mm cuffs for better grip
- Strong performance results in puncture resistance tests*
- Force at break (FAB) median value: 12.4 newtons
- Tested against chemotherapy drugs**

Ordering information


Item number	Size	Packaging
PRO400XS	X-small	50 ea/bx, 10 bx/cs
PRO400S	Small	50 ea/bx, 10 bx/cs
PRO400M	Medium	50 ea/bx, 10 bx/cs
PRO400L	Large	50 ea/bx, 10 bx/cs
PRO400XL	X-large	50 ea/bx, 10 bx/cs

Sustainable packaging technology

Flat-layered packaging

Colour-coded sizing

100% recyclable materials



Regulations and quality standards

- Class 1 non-sterile medical device, per Medical Device Regulation (EU) 2017/745 and the UK Medical Device Regulation
- Compliant with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended
- Compliant with standards EN ISO 21420; EN ISO 374-1; EN ISO 374-2; EN 16523-1; EN ISO 374-4 and EN ISO 374-5, EN 388
- Tested for viral penetration, in accordance with EN ISO 374-5



Medline Industries Ltd.
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com


Medline Ireland Ltd.
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  



Ideal for specialised areas, including:



*Test data is available upon request. **Do not use with Carmustine or Thiotepa.

SensiCare Pro gloves is a class I non-sterile medical device intended for use by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML577_EN 09/2024

Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant, France

CE 2797 BSI Group The Netherlands B.V.
1066 EP, Amsterdam, Netherlands