

SensiCare® Free

Accelerator-free non-sterile powder-free nitrile exam gloves



While healthcare professionals are well-informed about latex allergies, they may overlook allergies triggered by chemical accelerators in glove manufacturing (Type IV allergy; allergic contact dermatitis).

SensiCare Free is designed to help address this concern. Over 66.67 per cent of patients with dermatitis cases were symptom-free after switching to an accelerator-free glove.¹

Switch to Medline's SensiCare Free gloves to help protect and prevent the related allergies and negative impact that chemical accelerators may have on your healthcare workers.

Benefits and features

- Free from chemical accelerators
- Designed to offer a comfortable fit and reliable grip
- Acceptable quality level (AQL) for EN455-1: 1.0
- Force at break (FAB) median value: 8.2 newtons
- Personal protective equipment (PPE)-certified: Type B
- Tested against chemotherapy drugs*

Ordering information


Item number	Size	Packaging
MDSFREE-XS	X-Small	200 ea/bx, 10 bx/cs
MDSFREE-S	Small	200 ea/bx, 10 bx/cs
MDSFREE-M	Medium	200 ea/bx, 10 bx/cs
MDSFREE-L	Large	200 ea/bx, 10 bx/cs
MDSFREE-XL	X-Large	180 ea/bx, 10 bx/cs




Thickness
Finger: 0.09 mm
Palm: 0.07 mm
Cuff: 0.05 mm

Sustainable packaging technology

MicroStop™ film



Flat-layered packing



Colour-coded sizing



100% recyclable materials




Recommended for use in:



Regulations and quality standards

- Class I non-sterile medical device, per Medical Device Regulation (EU) 2017/745 and the UK Medical Device Regulation
- Compliant with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended
- Compliant with standards EN ISO 21420; EN ISO 374-1; EN ISO 374-2; EN 16523-1; EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Compliant with Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food



CE 2797 BSI Group The Netherlands B.V.
1066 EP, Amsterdam, Netherlands



Explore how Medline's gloves range, including SensiCare Free, makes a difference around the world. Contact your local Medline representative or visit our website: www.medline.eu/uk

¹Gunasegaran, J. et al. (2024) 'Review on prevalence, risk factors, and research advancements on the use of medical gloves concerning hand dermatitis among health care workers', Safety and Health at Work [Preprint]. doi:10.1016/j.shaw.2024.02.005

*Not for use with Carmustine or Thiotepa

This product is a class I non-sterile medical device intended for use by healthcare professionals. Before use, consult the instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML521_EN 06/2024

Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France