

SensiCare® Ice

Non-sterile powder-free nitrile exam gloves



**Protection
at every
touch point**



SensiCare® Ice

Non-sterile powder-free nitrile exam gloves

Benefits and features

- Designed to offer a comfortable fit and reliable grip
- Acceptance quality level (AQL) for EN455-1: 1.0
- Force at break (FAB) median value: 8.1 newtons
- Personal protective equipment (PPE)-certified
- Tested against chemotherapy drugs


Recommended for use in:



Innovative and sustainable packaging technology


MicroStop™ film

Covers the box opening to control the number of gloves dispensed, helping to reduce waste



Flat-layer packaging

Our proprietary robotic smart packing system stacks gloves in a single flat layer within each box, helping to reduce waste and improve efficiency and satisfaction



Colour-coded sizing

Our convenient colour-coded sizing system cues you to quickly and easily choose the right size, helping to reduce waste and inefficiency



Recyclable materials

Our boxes and cases are made of 70 – 100% recycled materials, all of which are 100% recyclable



Medline Exam Gloves

Protection at every touch point

Our commitment to you

Our gloves protect your front line, your bottom line, and your supply chain via customised solutions for all your exam glove needs.

25 billion
eaches sold
annually
worldwide

Consistency is key.

We're committed to quality, as seen in our:

5-point quality process

Solid infrastructure

Rigorous product testing

bsi.
Certificate of Registration
QUALITY MANAGEMENT SYSTEM - ISO 13485:2003 & EN ISO 13485:2012
This is to certify that:
Medline International France SAS
S, rue Charles Lindbergh
Créteil-Montant
93110
France
Holds Certificate Number: MD 595395
and operates a Quality Management System which complies with the requirements of ISO 13485:2003 & EN ISO 13485:2012 for the following scope:
Design, development, manufacture and distribution of sterile surgical drapes, gowns, equipment covers, procedure packs, surgical masks, surgical personnel covering, surgical equipment, medical gloves, skin barrier film and silicone and foam wound care dressings.
For and on behalf of BSI: Frank Lee, EMEA Compliance & Risk Director
Original Registration Date: 16/07/2013 Effective Date: 01/10/2023
Latest Revision Date: 15/03/2017 Expiry Date: 28/02/26

bsi.
Certificate of Registration
ENVIRONMENTAL MANAGEMENT SYSTEM - ISO 14001:2015
This is to certify that:
Medline Industries, Inc.
HQ Functions
3 Lakes Dr.
Northfield
Illinois
60093
USA
Holds Certificate No: EMS 578261
and operates an Environmental Management System which complies with the requirements of ISO 14001:2015 for the following scope:
The environmental management system for the control of risks associated with injection molding, kit assembly, wet wipes, distribution of adult disposable incontinence products and testing lab of disposable medical devices and associated equipment.
For and on behalf of BSI: Carlos Pianga, SVS System Certification and Compliance
Original Registration Date: 2012-05-16 Effective Date: 2015-05-16
Latest Revision Date: 2017-07-10 Expiry Date: 2018-05-15
Page: 1 of 2
...making excellence a habit!

We are committed to environmental, social and governance (ESG), taking a people and planet approach.



Audit and Ethical Sourcing

We advocate for labour and human rights in Medline's supply chain through our Ethical Sourcing Programme and Worldwide Accredited Production (WRAP) certification.




Environment and Climate Change

Our glove manufacturing is ISO 14001-certified, allowing us to assess environmental risks and measure energy savings, and help reduce carbon emissions.

Read the **2023 Medline Global ESG Report** to learn more about our sustainable initiatives.

Ordering information

Item number	Size	Packaging
486800	X-small	200 ea/bx, 10 bx/cs
486801	Small	200 ea/bx, 10 bx/cs
486802	Medium	200 ea/bx, 10 bx/cs
486803	Large	200 ea/bx, 10 bx/cs
486804	X-large	180 ea/bx, 10 bx/cs



Regulations and Quality Standards

- Class 1 non-sterile medical device per Regulation (EU) 2017/745 and UK Medical Device Regulation
- Complies with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended;
- Complies with standards EN ISO 21420, EN ISO 374-1, EN ISO 374-2, EN 16523-1, EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food



2797 BSI Group The Netherlands B.V.
1066 EP, Amsterdam, Netherlands

Explore how Medline gloves range, including SensiCare® Ice, makes a difference around the world. Contact your local Medline representative or visit our website: www.medline.eu/uk



Medline Industries Ltd.
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd.
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US

This product is a class I non-sterile medical device intended for use by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant, France

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1333_EN 08/2024