

SensiCare® Ice

Non-sterile powder-free nitrile exam gloves



**Protection
at every
touch point**



SensiCare Ice

Non-sterile powder-free nitrile exam gloves

Features and benefits

Acceptable quality limit (AQL): 1.0

Force at break (FAB - median value): 8.1 newtons

- Comfortable fit
- Reliable grip
- PPE-certified
- Chemo-tested

Colour: Violet blue

Textured fingertips

Length: 249 mm

Thickness

Finger: 0.09 mm

Palm: 0.07 mm

Cuff: 0.05 mm

Recommended for use in:



Labour and Delivery



Environmental services



Pharmacy/Oncology



Noncritical Care Unit



Phlebotomy Care Unit




Critical Care Unit



Innovative and sustainable packaging technology


MicroStop™ film

Covers the box opening to control the number of gloves dispensed, helping to reduce waste



Flat-layer packaging

Our proprietary robotic smart packing system stacks gloves in a single flat layer within each box, helping to reduce waste and improve efficiency and satisfaction



Colour-coded sizing

Our convenient colour-coded sizing system cues you to quickly and easily choose the right size, helping to reduce waste and inefficiency



Recyclable materials

Our boxes and cases are made of 100% recycled materials, all of which can be recycled further



Medline Exam Gloves

Protection at every touch point

Our commitment to you

Our gloves protect your front line, your bottom line, and your supply chain via customised solutions for all your exam glove needs.

25 billion
pairs sold
annually
worldwide

Consistency is key.

We're committed to quality, as seen in our:

**5-point
quality
process**

**Solid
infra-
structure**

**Rigorous
product
testing**

bsi.
Certificate of Registration
QUALITY MANAGEMENT SYSTEM - ISO 13485:2003 & EN ISO 13485:2012

This is to certify that:
Medline International France SAS
5, rue Charles Lindbergh
Clichy-sous-Bois
91410
France

Hold Certificate Number: **MD 595395**
and operates a Quality Management System which complies with the requirements of ISO 13485:2003 & EN ISO 13485:2012 for the following scope:
Design, development, manufacture and distribution of sterile surgical drapes, gowns, equipment covers, procedure pods, surgical masks, surgical personnel covering, scrub equipment, medical gloves, skin barrier film and silicone and foam wound care dressings

For and on behalf of BSI:
Frank Lee, BSI Compliance & Risk Director
Original Registration Date: 16/07/2013
Latest Review Date: 15/02/2017
Effective Date: 01/10/2013
Expiry Date: 26/02/2020

bsi.
Certificate of Registration
ENVIRONMENTAL MANAGEMENT SYSTEM - ISO 14001:2015

This is to certify that:
Medline Industries, Inc.
HQ Functions
3 Lakes Dr.
Northfield
Illinois
60093
USA

Hold Certificate No: **EMS 578261**
and operates an Environmental Management System which complies with the requirements of ISO 14001:2015 for the following scope:
The environmental management system for the control of risks associated with injection molding, M assembly, test setup, distribution of adult disposable incontinence products and testing lab of disposable medical devices and associated equipment.

For and on behalf of BSI:
Carlos Piango, BSI System Certification and Compliance
Original Registration Date: 23/02/05-16
Latest Review Date: 20/07/18
Effective Date: 2015-05-16
Expiry Date: 2018-05-15

We're committed to corporate social responsibility



Social impact

More than 3,000 jobs directly created or supported by our primary manufacturing facility



Environmental impact

ISO 14001 certification to assess environmental risk, measure energy savings and help reduce carbon emissions



WRAP Certification

Promoting labour and human rights in Medline's supply chain in partnership with Worldwide Accredited Production (WRAP)

Ordering information

Item number	Size	Length	Width	Packaging
486800	X-small	246 mm	78 mm	200 ea/bx, 10 bx/cs
486801	Small	250 mm	84 mm	200 ea/bx, 10 bx/cs
486802	Medium	249 mm	96 mm	200 ea/bx, 10 bx/cs
486803	Large	246 mm	115 mm	200 ea/bx, 10 bx/cs
486804	X-large	248 mm	114 mm	180 ea/bx, 10 bx/cs



Regulations and Quality Standards

- Class 1 non-sterile medical device per Regulation (EU) 2017/745 and UK Medical Device Regulation
- Complies with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended;
- Complies with standards EN 420, EN ISO 374-1, EN ISO 374-2, EN 16523-1, EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food



Explore how Medline gloves make a difference around the world. Contact your local Medline representative or visit our website: www.medline.eu/uk



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US

This product is a class I non-sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1333_EN 01/2023