



## SENSICARE® PI MICRO

### STERILE POLYISOPRENE POWDER-FREE SURGICAL GLOVES

- Unique synthetic Isolex polyisoprene formulation
  - Chemical structure identical to natural rubber latex, without the proteins
  - Offers healthcare workers an alternative solution that performs like latex in terms of comfort, fit and feel
- Powder-free helps to reduce powder-related skin irritation and post-operative complications associated with granulomas
- Thinner than Senticare® PI Surgical Gloves for an enhanced tactile response
- Excellent choice as an outer glove when double-gloving

Item Number	Size	Packaging
MSG9655	5 ½	50 Pairs/Box, 200 Pairs/Case
MSG9660	6	50 Pairs/Box, 200 Pairs/Case
MSG9665	6 ½	50 Pairs/Box, 200 Pairs/Case
MSG9670	7	50 Pairs/Box, 200 Pairs/Case
MSG9675	7 ½	50 Pairs/Box, 200 Pairs/Case
MSG9680	8	50 Pairs/Box, 200 Pairs/Case
MSG9685	8 ½	50 Pairs/Box, 200 Pairs/Case
MSG9690	9	50 Pairs/Box, 200 Pairs/Case

#### Recommended for use in:

General, neuro, vascular, plastic & eye surgery

#### Regulations and quality standards

- Certified to EN ISO 13485
- Medical Device Directive 93/42/EEC - Class IIA
- Personal Protective Equipment Directive 89/686/EEC - Class III

#### SPECIFICATIONS:

- Freedom from Holes per EN 455: AQL 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- No residual of accelerators
- Cuff design: tapered and beaded to reduce roll-down
- Finger tip: 0.20 mm
- Glove length (per EN 445-2, size 7.5): 299 mm
- Colour: cream
- Grip: smooth



To learn more about our products, please contact your Medline representative or visit us at [www.medline.eu/uk](http://www.medline.eu/uk)