

Signature Latex Green with Aloe Vera

Sterile and powder-free latex surgical gloves



- Natural rubber latex formulation
 - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established
 - **Caution:** This product contains natural rubber latex, which may cause allergic reactions, including anaphylactic response
- **Aloetouch® Moisturising Technology**
 - The world's first therapeutic medical glove technology
 - Clinically proven to moisturise skin
 - Contains 100% pure freeze-dried aloe vera gel on the inner surface
 - Assists with donning, in addition to moisturising hands
- Powder-free and coated with a synthetic polymer to assist with damp-hand donning
- Dark-green colour allows for rapid detection of perforations in the outer glove

Specifications

- Freedom from holes per EN 455-1: acceptance quality limit (AQL) 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- No residue from accelerators
- Cuff design: tapered and beaded to reduce roll-down
- Finger tip: 0.19 mm
- Glove length: (per EN 455-2, size 7.5): 291 mm
- Colour: dark green
- Texture: smooth

Ordering information

Item code	Size	Packaging
MSG5555	5 ½	50 pairs/box, 200 pairs/case
MSG5560	6	50 pairs/box, 200 pairs/case
MSG5565	6 ½	50 pairs/box, 200 pairs/case
MSG5570	7	50 pairs/box, 200 pairs/case
MSG5575	7 ½	50 pairs/box, 200 pairs/case
MSG5580	8	50 pairs/box, 200 pairs/case
MSG5585	8 ½	50 pairs/box, 200 pairs/case
MSG5590	9	50 pairs/box, 200 pairs/case

Recommended for use in

- Ideal for double-gloving in cardiac, trauma, ortho and high-risk procedures

Regulations and quality standards

- Factory EN ISO 13485-compliant
- Class IIa medical device; compliant with Medical Device Directive 93/42/EEC
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



This glove is a class IIa medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML575_EN 06/2023

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk