

# Signature Latex

## Sterile latex powder-free surgical gloves



- Made from natural rubber latex
  - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established
  - **Caution:** This product contains natural rubber latex, which may cause allergic reactions, including anaphylactic responses
- Powder-free feature helps reduce powder-related skin irritation and post-operative complications associated with granulomas
- Exclusive Medline SH8 glove mould with an independent thumb design allows for an anatomical fit and more natural movement in the fingers, thumbs and palms
- Beaded, tapered and textured cuff design reduces roll down

### Specifications

- Freedom from holes per EN455-1: AQL 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- No residual of accelerators
- Finger tip thickness: 0.22 mm
- Length: 294 mm
- Colour: light brown
- Grip: smooth

### Ordering information

Item code	Size	Packaging
MSG5455	5 ½	50 pairs/box, 200 pairs/case
MSG5460	6	50 pairs/box, 200 pairs/case
MSG5465	6 ½	50 pairs/box, 200 pairs/case
MSG5470	7	50 pairs/box, 200 pairs/case
MSG5475	7 ½	50 pairs/box, 200 pairs/case
MSG5480	8	50 pairs/box, 200 pairs/case
MSG5485	8 ½	50 pairs/box, 200 pairs/case
MSG5490	9	50 pairs/box, 200 pairs/case

### Recommended for use in

- General surgery

### Regulations and quality standards

- EN ISO 13485 compliant
- Class IIa sterile medical device; compliant with Medical Device Regulation (EU) 2017/745
- Category III personal protective equipment; compliant with PPE Regulation (EU) 2016/425



This glove is a class IIa sterile medical device intended to be used by healthcare professionals.

Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML571\_EN 01/2021

**Learn more.** Contact your Medline account manager or visit: [uk.medline.eu](https://uk.medline.eu)