

Signature Latex Micro

Sterile latex powder-free surgical gloves



- Natural rubber latex formulation
 - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established
 - **Caution:** This product contains natural rubber latex, which may cause allergic reactions, including anaphylactic responses
- Powder-free feature helps reduce powder-related skin irritation and post-operative complications associated with granulomas
- Exclusive Medline SH8 glove mould with an independent thumb design allows for an anatomical fit and more natural movement in the fingers, thumbs and palms
- Thinner than Signature Latex Essential surgical gloves for enhanced tactile sensitivity
- Excellent choice as outer gloves when double-gloving

Specifications

- Freedom from holes per EN455-1: AQL 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- No residual of accelerators
- Cuff design: beaded, tapered and textured to reduce roll down
- Finger tip thickness: 0.20 mm
- Length: 292 mm
- Colour: Light brown
- Grip: smooth

Ordering information

Item code	Size	Packaging
MSG5855	5 ½	50 pairs/box, 200 pairs/case
MSG5860	6	50 pairs/box, 200 pairs/case
MSG5865	6 ½	50 pairs/box, 200 pairs/case
MSG5870	7	50 pairs/box, 200 pairs/case
MSG5875	7 ½	50 pairs/box, 200 pairs/case
MSG5880	8	50 pairs/box, 200 pairs/case
MSG5885	8 ½	50 pairs/box, 200 pairs/case
MSG5890	9	50 pairs/box, 200 pairs/case

Recommended for use in

- General, neuro, vascular, plastic and eye surgery

Regulations and quality standards

- EN ISO 13485 compliant
- Class IIa sterile medical device; compliant with Medical Device Regulation (EU) 2017/745
- Category III personal protective equipment; compliant with PPE Regulation (EU) 2016/425



Medline International France SAS
5, rue Charles Lindbergh
44110 Châteaubriant, France

This glove is a class IIa sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML574_EN 01/2021

Learn more. Contact your Medline account manager or visit: uk.medline.eu