

Signature™ Latex OR

Interventional latex surgical glove



- Made from natural rubber latex
 - Contains 50 ug/g or less water-extractable protein
 - **Caution:** Safe use of these gloves by or on latex-sensitised individuals has not been established
 - **Caution:** This product contains natural rubber latex, which may cause allergic reactions, including anaphylactic response
- Powder-free to help prevent powder-related complications
- A versatile, sterile, powder-free latex option
- Space-saving folded packaging design

Specifications

- Freedom from holes, per EN 455-1: acceptable quality level (AQL) of 0.65
- Finger thickness: 0.22 mm
- Glove length (per EN 455-2, size 7.5): 289 mm
- Colour: cream
- Texture: microtextured
- Former (mould) design: anatomical, curved to help reduce hand fatigue
- Cuff design: tapered and beaded to reduce roll-down

Ordering information

Item code	Size	Packaging
MSG7255	5 ½	50 pairs/box, 200 pairs/case
MSG7260	6	50 pairs/box, 200 pairs/case
MSG7265	6 ½	50 pairs/box, 200 pairs/case
MSG7270	7	50 pairs/box, 200 pairs/case
MSG7275	7 ½	50 pairs/box, 200 pairs/case
MSG7280	8	50 pairs/box, 200 pairs/case
MSG7285	8 ½	50 pairs/box, 200 pairs/case
MSG7290	9	50 pairs/box, 200 pairs/case

Recommended for use in

- General surgery

Regulations and quality standards

- Factory EN ISO 13485-compliant
- Class IIa sterile medical device; compliant with Medical Device Regulation (MDR) (EU) 2017-745
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



STERILE R



EC REP

©BM Medical
Collamore House, Collamore
Road, Dalkey, Co. Dublin, Ireland

This glove is a class IIa sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1164_EN 09/2023

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk