



SIGNATURE LATEX ORTHO

STERILE LATEX POWDER-FREE SURGICAL GLOVES

- Natural rubber latex formulation
 - **Caution:** safe use of these gloves by or on latex-sensitised individuals has not been established
 - **Caution:** this product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses
- Powder-free helps to reduce powder-related skin irritation and post-operative complications associated with granulomas
- Exclusive Medline SH8 glove mould with an independent thumb design allows for an anatomical fit and more natural movement in the fingers, thumb and palm
- Thicker than standard latex surgical gloves for additional protection during high-risk, trauma and orthopedic procedures
- Inner polymer coating for effortless donning glove properties

Item Number	Size	Packaging
MSG5660	6	50 Pairs/Box, 200 Pairs/Case
MSG5665	6 ½	50 Pairs/Box, 200 Pairs/Case
MSG5670	7	50 Pairs/Box, 200 Pairs/Case
MSG5675	7 ½	50 Pairs/Box, 200 Pairs/Case
MSG5680	8	50 Pairs/Box, 200 Pairs/Case
MSG5685	8 ½	50 Pairs/Box, 200 Pairs/Case
MSG5690	9	50 Pairs/Box, 200 Pairs/Case

Recommended for use in:
Orthopedic surgery

Regulations and quality standards

- Certified to EN ISO 13485
- Medical Device Directive 93/42/EEC - Class IIA
- Personal Protective Equipment Directive 89/686/EEC - Category III

SPECIFICATIONS:

- Freedom from Holes per EN 455: AQL 0.65
- Synthetic polymer coating
- Former (mould) design: anatomical
- Cuff design: tapered and beaded to reduce roll-down
- No residual of accelerators
- Finger tip thickness: 0.31 mm
- Glove length: conforms to EN 455-2
- Colour: brown
- Grip: smooth



To learn more about our products, please contact your Medline representative or visit us at www.medline.eu/uk