



STERILE EXAM GLOVES

Effective Protection in Nitrile and Latex

When Sterility is Essential

Medline sterile exam gloves help you deliver better patient care when sterile technique is appropriate. Our sterile exam gloves are ambidextrous and can be used for: labor and delivery, emergency departments, dermatology, wound care and minor non-invasive procedures.

Clinical Sterility, Medline Dependability

Certified Sterility

Sterile exam gloves are sterilized with gamma irradiation and they are packaged to assure sterility. QA testing on every lot ensures we exceed market standards.

Convenience

Small packets—pocket sized—hold a pair of ambidextrous gloves. Compare that to a bulky surgical glove package! Medline sterile exam gloves make it easy to stay compliant with sterility protocols.

Cost-Effectiveness

At about 40% the cost of surgical gloves, sterile exam gloves are a smart alternative for controlling costs while maintaining the highest level of patient care.

Brands	Description	Item Numbers
NITRILE Medline® Procedure Nitrile, Pairs	296mm Nitrile Sterile Exam Glove Pairs powder free, blue, finger-tip textured, S, M, L, XL	MDS2294 to MDS2297 50 pr/bx, 6 bx/cs, 300 pr/cs
LATEX Medline® Procedure Latex, Pairs	243mm Latex Sterile Exam Glove, Pairs powder-free, cream, textured, S, M, L, XL	MDS8114 to MDS8117 50 pr/bx, 6 bx/cs, 300 pr/cs

Regulations and quality standards

Conforms to EU Directive 93/42/EEC - Class I Medical Device. Products meet the requirements of European harmonized standards EN 455-1, EN 455-2, EN 455-3 and EN 455-4. Tested and passed in accordance with ASTM F1671 for Viral Penetration. Certified to EN ISO 13485.



To learn more about our products, please contact your Medline representative or visit us at www.medline.com/uk