

A full-page photograph of a female surgeon in an operating room. She is wearing a blue bouffant cap, a white surgical mask, and blue surgical scrubs. She is looking down at a piece of blue fabric she is holding with her gloved hands. The background shows the sterile environment of an operating room with overhead lights and medical equipment.

**DESIGNED, TESTED AND TRUSTED
BY OPERATING THEATRE PROFESSIONALS.**

Surgical Gowns

Performance, Variety and Quality

Medline Surgical Gowns

Performance, variety and quality

An extensive choice


To satisfy any preference in the operating theatre, we offer five different fabric options for surgical gowns, each with its own distinct style and reinforcement.

Dependable protection

Medline's gown fabrics offer high protection. Each gown is latex-free and meets or exceeds the minimum performance requirements of the EN 13795 standard.* The Essential and Advanced materials have a high hydrostatic head, and the extra alcohol repellence treatment gives wearers the protection they need and the comfort they deserve. The unique Ultimate gown offers all-over impervious protection with exceptional breathability.



High performance with high comfort

All of our surgical gowns offer theatre staff excellent comfort and ease of use. High-tech breathable fabrics keep the wearer comfortable and dry, even during long procedures. The gowns feature sleeves with glued seams, which contribute to greater strength and comfort.



Easy to use

All Medline surgical gowns feature intuitive colour coding for easier gown selection in the operating theatre. Purple product labels and transfer tabs indicate standard performance, and blue is for high performance. In addition, we offer three different gown collar colours to quickly identify the gown style: standard, fabric reinforced or poly-reinforced.



*The EN 13795 is the European harmonised standard specifying requirements for single-use and reusable coverings, i.e. surgical gowns and surgical drapes, used as medical devices for patients, clinical staff and equipment and intended to prevent the transmission of infective agents between patients and clinical staff during invasive surgical procedures.

OPS™ UltraGard

Effective gown, affordable price

- Lightweight SMMMS (35 gsm) material with natural fluid-repellent protection
- Comfortable and breathable
- Gown available in standard and poly-reinforced protection

Item No.	Item Description	Qty/Case
E3505CEA	Surgical gown, S/M, 115 cm	36
E3506CEA/E3506DCE*	Surgical gown, L, 120 cm	36/48*
E3545CEA/E3545DCE*	Surgical gown, XL, 130 cm	36/48*
FE3545CE	Surgical gown, XL, 140 cm	32
E3516CEA/E3516DCE*	Surgical gown, XLL, 150 cm	32/42*
E3585CEA/E3585DCE*	Surgical gown, 2XLL, 160 cm	20/42*
E3500CEA	Surgical gown, poly-reinforced, S/M, 115 cm	32
E3510CEA	Surgical gown, poly-reinforced, L, 120 cm	32
E3540CEA	Surgical gown, poly-reinforced, XL, 140 cm	28
E3513CE	Surgical gown, poly-reinforced, XLL, 150 cm	24
E3581CE	Surgical gown, poly-reinforced, 2XLL, 160 cm	20



*These gowns are available in dispenser boxes as well.

UltraGard Gowns are also available in single-layer packaging without overwrap and hand towels. This configuration helps improve efficiency by reducing set-up time and the budgetary needs for gowns.

Item No.	Item Description	Qty/Case
E3606CE	Surgical gown, L, 120 cm, direct, ST	40
E3645CE	Surgical gown, XL, 130 cm, direct, ST	40
E3616CE	Surgical gown, XLL, 150 cm, direct, ST	36
E3640CE	Surgical gown, poly-reinforced, XL, 140 cm, direct, ST	32
E3613CE	Surgical gown, poly-reinforced, XLL, 150 cm, direct, ST	28
E3681CE	Surgical gown, poly-reinforced, 2XLL, 160 cm, direct, ST	24



OPS™ Essential

Enhanced protection

- Lightweight material SMMMS (35 gsm) with extra treatment for enhanced repellence
- Comfortable: offers exceptional drapability and breathability
- Essential poly-reinforced gowns comply with the EN 13795, EN 14605 and EN 14126 standards.

Item No.	Item Description	Qty/Case
S3505CEC	Surgical gown, S/M, 115 cm	36
S3506CEC	Surgical gown, L, 120 cm	36
S3545CEC	Surgical gown, XL, 130 cm	36
FS3545CEA	Surgical gown, XL, 140 cm	32
S3507CEC	Surgical gown, XL, 140 cm, finger loops	32
S3516CEC	Surgical gown, XLL, 150 cm	32
S3585CEC	Surgical gown, 2XLL, 170 cm	32
S3504CEC	Surgical gown, fabric reinforced, S/M, 115 cm	32
S3508CEC	Surgical gown, fabric reinforced, L, 120 cm	28
S3509CEC	Surgical gown, fabric reinforced, XL, 130 cm	28
S3512CEC	Surgical gown, fabric reinforced, XLL, 150 cm	28
S3511CEC	Surgical gown, fabric reinforced, 2XLL, 170 cm	24
S3500CEC/S3500P*	Surgical gown, poly-reinforced, S/M, 115 cm	32
S3510CEC/S3510P*	Surgical gown, poly-reinforced, L, 120 cm	32
S3540CEC/S3540P*	Surgical gown, poly-reinforced, XL, 140 cm	28
S3513CEB/S3513P*	Surgical gown, poly-reinforced, XLL, 150 cm	24
S3581CEB/S3581P*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20

*Compliant with EN 14605 and EN 14126, personal protective equipment category III



Essential Gowns are also available in single-layer packaging.

Item No.	Item Description	Qty/Case
S3606CE	Surgical gown, L, 120 cm, direct	40
S3645CE	Surgical gown, XL, 130 cm, direct	40
S3616CE	Surgical gown, XLL, 150 cm, direct	36
S3685CE	Surgical gown, 2XLL, 170 cm, direct	32
S3609CE	Surgical gown, fabric reinforced, XL, 130 cm, direct	30
S3612CE	Surgical gown, fabric reinforced, XLL, 150 cm, direct	28
S3611CE	Surgical gown, fabric reinforced, 2XLL, 170 cm, direct	24



OPS™ Advanced

No compromise between comfort and protection

- Breathable: quickly transfers heat away to keep the wearer cool and comfortable
- Higher-basis-weight soft SMMMS (41 gsm) for enhanced protection and comfort
- High barrier performance
- Optional raglan cut provides more roominess and easier movement.

Item No.	Item Description	Qty/Case
OP9505CE	Surgical gown, S/M, 120 cm	28
OP9506CEA	Surgical gown, L, 130 cm	32
OP9507CEA	Surgical gown, XL, 140 cm	32
OP9516CEB	Surgical gown, XLL, 150 cm	32
OP39516CE	Surgical gown, XLL, 150 cm, raglan sleeves	24
OP9575CE	Surgical gown, XXL, 150 cm	24
OP9585CE	Surgical gown, 2XLL, 170 cm	20
OP9590CE	Surgical gown, 3XL, 135 cm	24
OP9508CEA/OP9508DCE*	Surgical gown, fabric reinforced, L, 130 cm	24/40*
OP9809CE	Surgical gown, fabric reinforced, XL, 130 cm, mask and finger loops	28
OP9509CEA/OP9509DCE*	Surgical gown, fabric reinforced, XL, 140 cm	28/36*
OP9512CE/OP9512DCE*	Surgical gown, fabric reinforced, XLL, 150 cm	24/32*
OP39512CE	Surgical gown, fabric reinforced, XLL, 150 cm, raglan sleeves	20
OP9812CE	Surgical gown, fabric reinforced, XLL, 150 cm, mask and finger loops	24
OP9583CE/OP9583DCE*	Surgical gown, fabric reinforced, 2XLL, 170 cm	20/28*
OP9500CE	Surgical gown, poly-reinforced, S/M, 120 cm	28
OP9510CE/OP9510DCE*	Surgical gown, poly-reinforced, L, 130 cm	28/40*
OP9540CE/OP9540DCE*	Surgical gown, poly-reinforced, XL, 140 cm	24/36*
OP9513CE/OP9513DCE*	Surgical gown, poly-reinforced, XLL, 150 cm	20/36*
OP9581CE/OP9581DCE*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20/32*

*These gowns are available in dispenser boxes as well.



OPS™ Ultimate

The intelligent gown with breathable membrane technology

- High performance across the entire gown and fully impervious to blood-borne pathogens (AAMI* PB70, ASTM F1671 tested, EN 13795)
- Made of 60-gsm trilaminate material
- Breathable and cool: reacts to increasing body temperature by increasing water vapour transmission
- Recommended for all procedures where breathable all-over protection is needed, such as long cardiological, orthopaedic and transplant surgeries

Item No.	Item Description	Qty/Case
LW89045CEB	Breathable impervious surgical gown, XL, 130 cm	20
LW32474CEB	Breathable impervious surgical gown, XLL, 150 cm, raglan sleeves	20
LW39079CEB	Breathable impervious surgical gown, 2XLL, 150 cm, raglan sleeves	20



We also offer surgical gowns featuring special fabrics tested against chemotherapy and cytotoxic drugs.

Item No.	Item Description	Qty/Case
DYNJPE2301P	Impervious OR gown, prevention plus, L, 114 cm	24
DYNJPE2302P	Impervious OR gown, prevention plus, XL, 142 cm	24

*AAMI: Association for the Advancement of Medical Instrumentation. The AAMI PB70 is the US reference standard establishing a system of classification for protective apparel and drapes based on their liquid barrier performance and it specifies related labelling requirements and standardised testing methods.

Eclipse®

A soft and comfortable gown

- Outstanding comfort and breathability
- Spunlace (68 gsm) with a textile-like appearance
- Suitable protection

Item No.	Item Description	Qty/Case
9505CEA	Surgical gown, S/M, 115 cm	28
9506CEA/9506DCE*	Surgical gown, L, 120 cm	28/40*
9415CEA	Surgical gown, L, 120 cm, raglan sleeves	26
9507CEA/9507DCE*	Surgical gown, XL, 130 cm	24/40*
9707CE	Surgical gown, XL, 130 cm, finger loops	24
9516CEA/9516DCE*	Surgical gown, XLL, 150 cm	24/36*
9716CE	Surgical gown, XLL, 150 cm, finger loops	24
9575CEA/9575DCE*	Surgical gown, XXL, 130 cm	24/36*
9585CE/9585DCE*	Surgical gown, 2XLL, 170 cm	24/32*
9508CEA	Surgical gown, fabric reinforced, L, 120 cm	24
9509CEA/9509DCE*	Surgical gown, fabric reinforced, XL, 130 cm	22/32*
9539CEA	Surgical gown, fabric reinforced, XL, 130 cm	7
9709CE	Surgical gown, fabric reinforced, XL, 130 cm, finger loops	22
9809CE	Surgical gown, fabric reinforced, XL, 130 cm, mask and finger loops	22
9571CEA/9571DCE*	Surgical gown, fabric reinforced, XXL, 130 cm	20/28*
9512CEA/9512DCE	Surgical gown, fabric reinforced, XLL, 150 cm	22/32*
9451CEA	Surgical gown, fabric reinforced, XLL, 150 cm	11
9712CE	Surgical gown, fabric reinforced, XLL, 150 cm, finger loops	22
9812CEA	Surgical gown, fabric reinforced, XLL, 150 cm, mask and finger loops	20
9581CEA/9581DCE*	Surgical gown, fabric reinforced, 2XLL, 170 cm	19/28*
9000CEA	Surgical gown, poly-reinforced, S/M, 115 cm	28
9010CEA/9010DCE*	Surgical gown, poly-reinforced, L, 120 cm	28/32*
9040CEA	Surgical gown, poly-reinforced, XL, 130 cm	24
9513CEA	Surgical gown, poly-reinforced, XLL, 150 cm	24
9580CE/9580DCE*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20/28*

*This item is also available packed in a dispenser box.



Surgical Gown Selection Guide

Choose the right gown with the right performance for every surgery, based on the type and length of the procedure, amount of fluids and required comfort.

Range	Eclipse			UltraGard™		Essential			Advanced			Ultimate	
	S	F	P	S	P	S	F	P	S	F	P	S	
Standard/Fabric/Poly-reinforced													
Ophthalmic	✓			✓		✓			✓				
Tracheostomy	✓			✓		✓			✓				
Ear surgery	✓			✓		✓			✓				
Cochlea implants	✓			✓		✓			✓				
Mouth	✓			✓		✓			✓				
Tonsillectomy	✓			✓		✓			✓				
Nose surgery	✓			✓		✓			✓				
Varicose veins	✓	✓		✓		✓	✓		✓	✓			
Proctology	✓	✓		✓		✓	✓		✓	✓			
Obstetrics		✓	✓		✓		✓	✓		✓	✓		✓
Maxillofacial	✓	✓		✓		✓	✓		✓	✓			
Endoabdominal	✓			✓		✓			✓				
Endocrine	✓			✓		✓			✓				
Pacemakers	✓			✓		✓			✓				
General orthopaedics	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓
Breast procedures	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Burns, plastics and grafting		✓	✓		✓		✓	✓		✓	✓		✓
Angiography/angioplasty	✓	✓		✓		✓	✓		✓	✓	✓		
Gynaecological	✓	✓		✓		✓	✓		✓	✓			
General abdominal surgery		✓	✓		✓		✓	✓		✓	✓		✓
Urogenital		✓	✓		✓		✓	✓	✓	✓	✓		✓
Orthopaedics, hips		✓	✓		✓		✓	✓		✓	✓		✓
Kidney, adrenals and bladder	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Orthopaedics, knees		✓	✓		✓		✓	✓		✓	✓		✓
Thoracic		✓	✓		✓		✓	✓		✓	✓		✓
Cardiovascular surgery		✓	✓		✓		✓	✓		✓	✓		✓
Cranial		✓	✓		✓		✓	✓		✓	✓		✓
Spinal procedures		✓	✓		✓		✓	✓		✓	✓		✓

DRG Surgical Procedures

For more information about our products, please visit www.medline.eu/uk or speak with your local account manager.



ALWAYS ON.

Medline Industries Ltd
 3rd Floor
 Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom
 Tel.: +44 844 334 5237
 Fax: +44 844 334 5238

www.medline.eu/uk
uk-customerservice@medline.com

