

Surgical Gowns

Performance, quality and variety



**Designed, tested
and trusted
by operating
theatre
professionals**

Medline Surgical Gowns

Performance, quality and variety


An extensive choice

To satisfy preferences in the operating theatre, we offer five different fabric options for surgical gowns, each with a distinct style and reinforcement.

Dependable protection

Medline's gown fabrics offer high-level protection. Each gown is latex-free and meets or exceeds the minimum performance requirements of the EN 13795-1 standard.* The Essential and Advanced materials have a good hydrostatic head, and the extra repellent treatment gives patients the protection they need. The unique Ultimate gown offers all-over impervious protection with exceptional breathability.


*The EN 13795-1 is the European standard, specifying requirements for single-use and reusable coverings—surgical gowns and drapes used as medical devices for patients, clinical staff and equipment, intended to prevent the transmission of infective agents during invasive surgical procedures.



High performance with high comfort

All our surgical gowns offer operating theatre staff excellent comfort and ease of use. High-tech, breathable fabrics keep the wearer comfortable and dry, even during long procedures. The gowns feature sleeves with glued seams, which contribute to greater strength and comfort.

Average sleeve seam tensile strength (N)





Easy to use

All Medline surgical gowns feature intuitive colour coding for easier gown selection in the operating theatre. Purple product labels and transfer tabs indicate standard performance, and blue is for high performance. In addition, we offer three different gown collar colours to quickly identify the gown style—standard, fabric-reinforced or poly-reinforced.

Purple collar:
standard gowns

White collar:
reinforced gowns

Blue collar:
poly-reinforced
gowns



Prioritised safety

The standard ranges of Medline surgical gowns offer an SMS wrap solution, preferred for its low-linting and low-memory effect, therefore, decreasing the risk of contamination.



OPS™ UltraGard

Effective gown, affordable price

- Lightweight SMMMS (35 gsm) material with natural, fluid-repellent protection
- Comfortable and breathable
- Gown available with either standard or poly-reinforced protection

Item number	Description	Qty/case
E3505CEA	Surgical gown, S/M, 115 cm	36
E3506CEA	Surgical gown, L, 120 cm	36/48*
E3545CEA/E3545DCE*	Surgical gown, XL, 130 cm	36/48*
FE3545CE	Surgical gown, XL, 140 cm	32
E3516CEA	Surgical gown, XLL, 150 cm	32/42*
E3585CEA/E3585DCE*	Surgical gown, 2XLL, 160 cm	20/42*
E3500CEA	Surgical gown, poly-reinforced, S/M, 115 cm	32
E3510CEA	Surgical gown, poly-reinforced, L, 120 cm	32
E3540CEA	Surgical gown, poly-reinforced, XL, 140 cm	28
E3513CE	Surgical gown, poly-reinforced, XLL, 150 cm	24
E3581CE	Surgical gown, poly-reinforced, 2XLL, 160 cm	20

*These gowns are also available in dispenser boxes.



OPS™ Essential

Enhanced protection

- Lightweight SMMMS (35 gsm) material with extra treatment for enhanced repellency
- Comfortable; offers exceptional drapeability and breathability
- Essential poly-reinforced gowns comply with EN 13795-1, EN 14605 and EN 14126 standards

Item number	Description	Qty/case
S3505CEC	Surgical gown, S/M, 115 cm	36
S3506CEC	Surgical gown, L, 120 cm	36
S3545CEC	Surgical gown, XL, 130 cm	36
FS3545CEA	Surgical gown, XL, 140 cm	32
S3516CEC	Surgical gown, XLL, 150 cm	32
S3585CEC	Surgical gown, 2XLL, 170 cm	32
S3504CEC	Surgical gown, fabric-reinforced, S/M, 115 cm	32
S3508CEC	Surgical gown, fabric-reinforced, L, 120 cm	28
S3509CEC	Surgical gown, fabric-reinforced, XL, 130 cm	28
S3512CEC	Surgical gown, fabric-reinforced, XLL, 150 cm	28
S3511CEC	Surgical gown, fabric-reinforced, 2XLL, 170 cm	24
S3500CEC	Surgical gown, poly-reinforced, S/M, 115 cm	32
S3510CEC	Surgical gown, poly-reinforced, L, 120 cm	32
S3540CEC	Surgical gown, poly-reinforced, XL, 140 cm	28
S3513CEB	Surgical gown, poly-reinforced, XLL, 150 cm	24
S3581CEB	Surgical gown, poly-reinforced, 2XLL, 170 cm	20
S3591CE	Surgical gown, poly-reinforced, 3XL, 140 cm	24
S3587CE	Surgical gown, poly-reinforced, 3XLL, 170 cm	20



OPS™ Essential gowns are also available with finger loops and/or masks.

Item number	Description	Qty/case
S3507CEC	Surgical gown, XL, 140 cm, finger loops	32
S3740CE	Surgical gown, poly-reinforced, XL, 140 cm, finger loops	28
S3713CE	Surgical gown, poly-reinforced, XLL, 150 cm, finger loops	24
S3840CE	Surgical gown, poly-reinforced, XL, 140 cm, mask and finger loops	28
S3813CE	Surgical gown, poly-reinforced, XLL, 150 cm, mask and finger loops	24
S3881CE	Surgical gown, poly-reinforced, 2XLL, 170 cm, mask and finger loops	20

OPS™ Advanced

No compromise between comfort and protection

- Breathable; quickly transfers heat away to keep the wearer cool and comfortable
- Higher basis weight, soft SMMMS (41 gsm) for enhanced protection and comfort
- High barrier performance
- Optional raglan cut increases roominess and provides easier movement

Item number	Description	Qty/case
OP9505CE	Surgical gown, S/M, 120 cm	28
OP9506CEA	Surgical gown, L, 130 cm	32
OP9507CEA	Surgical gown, XL, 140 cm	32
OP9516CEB	Surgical gown, XLL, 150 cm	32
OP39516CE	Surgical gown, XLL, 150 cm, raglan sleeves	24
OP9575CE	Surgical gown, XXL, 150 cm	24
OP9585CE	Surgical gown, 2XLL, 170 cm	20
OP9590CE	Surgical gown, 3XL, 135 cm	24
OP9508CEA/OP9508DCE*	Surgical gown, fabric-reinforced, L, 130 cm	24/40*
OP9509CEA	Surgical gown, fabric-reinforced, XL, 140 cm	28
OP9512CE/OP9512DCE*	Surgical gown, fabric-reinforced, XLL, 150 cm	24/32*
OP39512CE	Surgical gown, fabric-reinforced, XLL, 150 cm, raglan sleeves	20
OP9583CE/OP9583DCE*	Surgical gown, fabric-reinforced, 2XLL, 170 cm	20/28*
OP9500CE	Surgical gown, poly-reinforced, S/M, 120 cm	28
OP9510CE	Surgical gown, poly-reinforced, L, 130 cm	28/40*
OP9540CE	Surgical gown, poly-reinforced, XL, 140 cm	24/36*
OP9513CE	Surgical gown, poly-reinforced, XLL, 150 cm	20
OP9581CE	Surgical gown, poly-reinforced, 2XLL, 170 cm	20



*These gowns are also available in dispenser boxes.

OPS™ Advanced gowns are also available with finger loops and/or masks.

Item number	Description	Qty/case
OP9809CE	Surgical gown, fabric-reinforced, XL, 130 cm, mask and finger loops	28
OP9812CE	Surgical gown, fabric-reinforced, XLL, 150 cm, mask and finger loops	24

A soft and comfortable gown

- Outstanding comfort and breathability
- Spunlace (68 gsm) with a textile-like appearance
- Suitable protection

Item number	Description	Qty/case
9505CEA	Surgical gown, S/M, 115 cm	28
9506CEA/9506DCE*	Surgical gown, L, 120 cm	28/40*
9415CEA	Surgical gown, L, 120 cm, raglan sleeves	26
9507CEA/9507DCE*	Surgical gown, XL, 130 cm	24/40*
9516CEA/9516DCE*	Surgical gown, XXL, 150 cm	24/36*
9575CEA	Surgical gown, XXL, 130 cm	24
9585CE/9585DCE*	Surgical gown, 2XLL, 170 cm	24/32*
9508CEA	Surgical gown, fabric-reinforced, L, 120 cm	24
9509CEA/9509DCE*	Surgical gown, fabric-reinforced, XL, 130 cm	22/32*
9539CEA	Surgical gown, fabric-reinforced, XL, 130 cm	7
9571CEA	Surgical gown, fabric-reinforced, XXL, 130 cm	20
9512CEA/9512DCE	Surgical gown, fabric-reinforced, XLL, 150 cm	22/32*
9451CEA	Surgical gown, fabric-reinforced, XLL, 150 cm	11
9581CEA/9581DCE*	Surgical gown, fabric-reinforced, 2XLL, 170 cm	19/28*
9000CEA	Surgical gown, poly-reinforced, S/M, 115 cm	28
9010CEA/9010DCE*	Surgical gown, poly-reinforced, L, 120 cm	28/32*
9040CEA	Surgical gown, poly-reinforced, XL, 130 cm	24
9513CEA	Surgical gown, poly-reinforced, XLL, 150 cm	24
9580CE/9580DCE*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20/28*

*This item is also available in dispenser boxes.



Eclipse® gowns are also available with finger loops and/or masks.

Item number	Description	Qty/case
9707CE	Surgical gown, XL, 130 cm, finger loops	24
9716CE	Surgical gown, XLL, 150 cm, finger loops	24
9709CE	Surgical gown, fabric-reinforced, XL, 130 cm, finger loops	22
9809CE	Surgical gown, fabric-reinforced, XL, 130 cm, mask and finger loops	22
9712CE	Surgical gown, fabric-reinforced, XLL, 150 cm, finger loops	22
9812CEA	Surgical gown, fabric-reinforced, XLL, 150 cm, mask and finger loops	20

OPS™ Ultimate

The intelligent gown with breathable membrane technology

- High performance across the entire gown (EN 13795-1 compliant)
- Made of 60-gsm trilaminate material
- Breathable and cool
- Recommended for procedures where breathable protection all over is needed, such as long cardiological, orthopaedic and transplant surgeries

Item number	Description	Qty/case
LW89045CEB	Surgical gown, breathable, impervious, XL, 130 cm	20
LW32474CEB	Surgical gown, breathable, impervious, XLL, 150 cm, raglan sleeves	20
LW39079CEB	Surgical gown, breathable, impervious, 2XLL, 150 cm, raglan sleeves	20



Urology Gowns

Optimised coverage ideal for urology procedures

Item number	Description	Qty/case
DYNJPE2304	Surgical gown, impervious, size, with pleat, XL	30
DYNJPE2305	Surgical gown, impervious, size, with pleat, XLL	30



Dual Protection

Advanced viral and chemical protection

Enhanced clinician protection for your peace of mind. We offer two ranges of surgical gowns with further testing against chemicals and infective agents. This solution is classified as Medical Device class I (MDD) and Personal Protective Equipment category III (PPE).



OPS™ Essential

Essential dual-certified gowns comply with EN 14605 and EN 14126.

Item number	Description	Qty/case
S3500P	Surgical gown, poly-reinforced, S/M, 115 cm	32
S3510P	Surgical gown, poly-reinforced, L, 120 cm	32
S3540P	Surgical gown, poly-reinforced, XL, 140 cm	28
S3513P	Surgical gown, poly-reinforced, XLL, 150 cm	24
S3581P	Surgical gown, poly-reinforced, 2XLL, 170 cm	20

Compliant with EN 14605 and EN 14126, personal protective equipment category III

Prevention™ Plus

As the end user, you have enhanced protection against fluids and cytotoxic drugs (chemotherapy). This is facilitated by an innovative design with a cooling and impervious three-layer material—ideal for protecting clinicians when performing HIPEC (hyperthermic intraperitoneal chemotherapy) treatments.

Item number	Description	Qty/case
DYNJPE2301P	Surgical gown, impervious, L, 114 cm	24
DYNJPE2302P	Surgical gown, impervious, XL, 142 cm	24

Prevention Plus gowns comply with EN 14605 and PPE Regulation (EU) 2016/425.



‘Peel-and-Go’ Solution

In addition to simplifying packaging for staff and eliminating costly non-essential overwrap and hand towels, the ‘peel-and-go’ packaging reduces your hospital’s waste. With a single, peelable sterilisation pouch, the ‘peel-and-go’ is the efficient, economic, and environmentally-friendly solution to your facility’s needs.

OPS™ Ultragard gowns in single-layer packaging are available in standard and poly-reinforced ranges.

Item number	Description	Qty/case
E3606CE	Surgical gown, L, 120 cm, direct, ST	40
E3645CE	Surgical gown, XL, 130 cm, direct, ST	40
E3616CE	Surgical gown, XLL, 150 cm, direct, ST	36
E3640CE	Surgical gown, poly-reinforced, XL, 140 cm, direct, ST	32
E3613CE	Surgical gown, poly-reinforced, XLL, 150 cm, direct, ST	28
E3681CE	Surgical gown, poly-reinforced, 2XLL, 160 cm, direct, ST	24



OPS™ Essential gowns in single-layer packaging are available in standard and fabric-reinforced ranges.

Item number	Description	Qty/case
S3606CE	Surgical gown, L, 120 cm, direct	40
S3645CE	Surgical gown, XL, 130 cm, direct	40
S3616CE	Surgical gown, XLL, 150 cm, direct	36
S3685CE	Surgical gown, 2XLL, 170 cm, direct	32
S3609CE	Surgical gown, fabric-reinforced, XL, 130 cm, direct	30
S3612CE	Surgical gown, fabric-reinforced, XLL, 150 cm, direct	28
S3611CE	Surgical gown, fabric-reinforced, 2XLL, 170 cm, direct	24



Surgical Gown Selection Guide

Choose the gown with the right performance for every surgery, based on the type and length of the procedure, amount of fluids and required comfort.

	Range Standard/Fabric-reinforced/ Poly-reinforced	Eclipse			UltraGard		Essential			Advanced			Ultimate	
		S	F	P	S	P	S	F	P	S	F	P	S	
DRG Surgical Procedures	Ophthalmic	✓			✓		✓			✓				
	Tracheostomy	✓			✓		✓			✓				
	Ear surgery	✓			✓		✓			✓				
	Cochlea implants	✓			✓		✓			✓				
	Mouth	✓			✓		✓			✓				
	Tonsillectomy	✓			✓		✓			✓				
	Nose surgery	✓			✓		✓			✓				
	Varicose veins	✓	✓		✓		✓	✓		✓	✓			
	Proctology	✓	✓		✓		✓	✓		✓	✓			
	Obstetrics		✓	✓		✓		✓	✓		✓	✓		✓
	Maxillofacial	✓	✓		✓		✓	✓		✓	✓			
	Endoabdominal	✓			✓		✓			✓				
	Endocrine	✓			✓		✓			✓				
	Pacemakers	✓			✓		✓			✓				
	General orthopaedics	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓
	Breast procedures	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	Burns, plastics and grafting		✓	✓		✓		✓	✓		✓	✓		✓
	Angiography/angioplasty	✓	✓		✓		✓	✓		✓	✓	✓		
	Gynaecological	✓	✓		✓		✓	✓		✓	✓			
	General abdominal surgery		✓	✓		✓		✓	✓		✓	✓		✓
	Urogenital		✓	✓		✓		✓	✓	✓	✓	✓		✓
	Orthopaedics, hips		✓	✓		✓		✓	✓		✓	✓		✓
	Kidney, adrenals and bladder	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	Orthopaedics, knees		✓	✓		✓		✓	✓		✓	✓		✓
	Thoracic		✓	✓		✓		✓	✓		✓	✓		✓
	Cardiovascular surgery		✓	✓		✓		✓	✓		✓	✓		✓
	Cranial		✓	✓		✓		✓	✓		✓	✓		✓
	Spinal procedures		✓	✓		✓		✓	✓		✓	✓		✓

For more information on our products, please contact
your account manager or visit us at www.medline.eu/uk



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk
uk-customerservice@medline.com

FOLLOW US  

These products are class I sterile medical devices intended to be used by healthcare professionals.
Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML281_EN 08/2021

 Medline International France SAS
5 Rue Charles Lindbergh,
44110 Châteaubriant, France

 2797 (BSI)