

# Surgical Gowns

Performance, quality and variety



**Designed, tested  
and trusted  
by operating  
theatre  
professionals**

# Medline Surgical Gowns

Performance, quality and variety


## An extensive choice

To satisfy preferences in the operating theatre, we offer five different fabric options for surgical gowns, each with a distinct style and reinforcement.

## Dependable protection

Medline's gown fabrics offer high-level protection. Each gown is latex-free and meets or exceeds the minimum performance requirements of the EN 13795-1 standard.\* The Essential and Advanced materials have a good hydrostatic head and an extra repellent treatment gives patients the protection they need. The unique Ultimate gown offers all-over impervious protection with exceptional breathability.


\*The EN 13795-1 is the European standard, specifying requirements for single-use and reusable coverings—surgical gowns and drapes used as medical devices for patients, clinical staff and equipment, intended to prevent the transmission of infective agents during invasive surgical procedures.



### High performance with high comfort

All our surgical gowns offer operating theatre staff excellent comfort and ease of use. High-tech, breathable fabrics keep the wearer comfortable and dry, even during long procedures. The gowns feature sleeves with glued seams, which contribute to greater strength and comfort.

Average sleeve seam tensile strength (N)





### Easy to use

All Medline surgical gowns feature intuitive colour coding for easier gown selection in the operating theatre. Purple product labels and transfer tabs indicate standard performance, and blue is for high performance. In addition, we offer three different gown collar colours to quickly identify the gown style—standard, fabric-reinforced or poly-reinforced.

Purple collar:  
standard gowns

White collar:  
reinforced gowns

Blue collar:  
poly-reinforced  
gowns



### Prioritised safety

The standard ranges of Medline surgical gowns offer an SMS wrap solution, preferred for its low-linting and low-memory effect; therefore, decreasing the risk of contamination.



# Organising your operating theatre

## Reusable gown dispenser

Require operating theatre organisation at your facility but don't want to resort to using disposable cardboard dispensers? We understand, and we agree. To minimise waste, without impacting your system, consider pairing your Medline gown mix with our reusable dispenser solution. These dispensers, with clear fronts to view remaining stock levels, can easily be replenished while maintaining FIFO principles (first in, first out).



## OPS™ UltraGard

### Effective gown, affordable price

- Lightweight SMS (35 gsm) material with natural, fluid-repellent protection
- Comfortable and breathable
- Gown available with either standard or poly-reinforced protection

Item number	Description	Qty/case
E3505CEA	Surgical gown, S/M, 115 cm	36
E3506CEA	Surgical gown, L, 120 cm	36
E3545CEA	Surgical gown, XL, 130 cm	36
FE3545CE	Surgical gown, XL, 140 cm	32
E3516CEA	Surgical gown, XLL, 150 cm	32
E3585CEA	Surgical gown, 2XLL, 160 cm	20
E3500CEA	Surgical gown, poly-reinforced, S/M, 115 cm	32
E3510CEA	Surgical gown, poly-reinforced, L, 120 cm	32
E3540CEA	Surgical gown, poly-reinforced, XL, 140 cm	28
E3513CE	Surgical gown, poly-reinforced, XLL, 150 cm	24
E3581CE	Surgical gown, poly-reinforced, 2XLL, 160 cm	20





# OPS™ Essential

## Enhanced protection

- Lightweight SMS (35 gsm) material with extra treatment for enhanced repellency
- Comfortable; offers exceptional drapeability and breathability
- Available with standard, fabric-reinforced or poly-reinforced protection

Item number	Description	Qty/case
S3505CEC	Surgical gown, S/M, 115 cm	36
S3506CEC	Surgical gown, L, 120 cm	36
S3545CEC	Surgical gown, XL, 130 cm	36
FS3545CEA	Surgical gown, XL, 140 cm	32
S3516CEC	Surgical gown, XLL, 150 cm	32
S3585CEC	Surgical gown, 2XLL, 170 cm	32
S3504CEC	Surgical gown, fabric-reinforced, S/M, 115 cm	32
S3508CEC	Surgical gown, fabric-reinforced, L, 120 cm	28
S3509CEC	Surgical gown, fabric-reinforced, XL, 130 cm	28
S3512CEC	Surgical gown, fabric-reinforced, XLL, 150 cm	28
S3511CEC	Surgical gown, fabric-reinforced, 2XLL, 170 cm	24
S3500CEC	Surgical gown, poly-reinforced, S/M, 115 cm	32
S3510CEC	Surgical gown, poly-reinforced, L, 120 cm	32
S3540CEC	Surgical gown, poly-reinforced, XL, 140 cm	28
S3513CEB	Surgical gown, poly-reinforced, XLL, 150 cm	24
S3581CEB	Surgical gown, poly-reinforced, 2XLL, 170 cm	20
S3591CE	Surgical gown, poly-reinforced, 3XL, 140 cm	24
S3587CE	Surgical gown, poly-reinforced, 3XL, 170 cm	20



OPS™ Essential gowns are also available with finger loops and/or masks.

Item number	Description	Qty/case
S3507CEC	Surgical gown, XL, 140 cm, finger loops	32
S3740CE	Surgical gown, poly-reinforced, XL, 140 cm, finger loops	28
S3713CE	Surgical gown, poly-reinforced, XLL, 150 cm, finger loops	24
S3840CE	Surgical gown, poly-reinforced, XL, 140 cm, mask and finger loops	28
S3813CE	Surgical gown, poly-reinforced, XLL, 150 cm, mask and finger loops	24
S3881CE	Surgical gown, poly-reinforced, 2XLL, 170 cm, mask and finger loops	20

# OPS™ Advanced

No compromise between comfort and protection

- Higher basis weight, SMS (41 gsm) for enhanced protection
- Breathable fabric with additional soft finish for optimal comfort
- Available with standard, fabric-reinforced or poly-reinforced protection
- Optional raglan cut increases roominess and provides easier movement

Item number	Description	Qty/case
OP9505CE	Surgical gown, S/M, 120 cm	28
OP9506CEA	Surgical gown, L, 130 cm	32
OP9507CEA	Surgical gown, XL, 140 cm	32
OP9516CEB	Surgical gown, XLL, 150 cm	32
OP39516CE	Surgical gown, XLL, 150 cm, raglan sleeves	24
OP9575CE	Surgical gown, XXL, 150 cm	24
OP9585CE	Surgical gown, 2XLL, 170 cm	20
OP9590CE	Surgical gown, 3XL, 135 cm	24
OP9508CEA	Surgical gown, fabric-reinforced, L, 130 cm	24
OP9509CEA	Surgical gown, fabric-reinforced, XL, 140 cm	28
OP9512CE	Surgical gown, fabric-reinforced, XLL, 150 cm	24
OP39512CE	Surgical gown, fabric-reinforced, XLL, 150 cm, raglan sleeves	20
OP9583CE	Surgical gown, fabric-reinforced, 2XLL, 170 cm	20
OP9500CE	Surgical gown, poly-reinforced, S/M, 120 cm	28
OP9510CE	Surgical gown, poly-reinforced, L, 130 cm	28
OP9540CE	Surgical gown, poly-reinforced, XL, 140 cm	24
OP9513CE	Surgical gown, poly-reinforced, XLL, 150 cm	20
OP9581CE	Surgical gown, poly-reinforced, 2XLL, 170 cm	20



OPS™ Advanced gowns are also available with finger loops and/or masks.

Item number	Description	Qty/case
OP9809CE	Surgical gown, fabric-reinforced, XL, 130 cm, mask and finger loops	28
OP9812CE	Surgical gown, fabric-reinforced, XLL, 150 cm, mask and finger loops	24

## A soft and comfortable gown

- Optimal comfort and breathability
- High tear-resistance
- Spunlace material with a textile-like appearance and feel
- Available with standard, fabric-reinforced or poly-reinforced protection

Item number	Description	Qty/case
9505CEA	Surgical gown, S/M, 115 cm	28
9506CEA/9506DCE*	Surgical gown, L, 120 cm	28/40*
9415CEA	Surgical gown, L, 120 cm, raglan sleeves	26
9507CEA/9507DCE*	Surgical gown, XL, 130 cm	24/40*
9516CEA/9516DCE*	Surgical gown, XLL, 150 cm	24/36*
9575CEA	Surgical gown, XXL, 130 cm	24
9585CE/9585DCE*	Surgical gown, 2XLL, 170 cm	24/32*
9508CEA	Surgical gown, fabric-reinforced, L, 120 cm	24
9509CEA/9509DCE*	Surgical gown, fabric-reinforced, XL, 130 cm	22/32*
9539CEA	Surgical gown, fabric-reinforced, XL, 130 cm	7
9571CEA	Surgical gown, fabric-reinforced, XXL, 130 cm	20
9512CEA	Surgical gown, fabric-reinforced, XLL, 150 cm	22
9451CEA	Surgical gown, fabric-reinforced, XLL, 150 cm	11
9581CEA/9581DCE*	Surgical gown, fabric-reinforced, 2XLL, 170 cm	19/28*
9000CEA	Surgical gown, poly-reinforced, S/M, 115 cm	28
9010CEA	Surgical gown, poly-reinforced, L, 120 cm	28
9040CEA	Surgical gown, poly-reinforced, XL, 130 cm	24
9513CEA	Surgical gown, poly-reinforced, XLL, 150 cm	24
9580CE/9580DCE*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20/28*

\*This item is also available in dispenser boxes.



Eclipse® gowns are also available with finger loops and/or masks.

Item number	Description	Qty/case
9707CE	Surgical gown, XL, 130 cm, finger loops	24
9716CE	Surgical gown, XLL, 150 cm, finger loops	24
9709CE	Surgical gown, fabric-reinforced, XL, 130 cm, finger loops	22
9712CE	Surgical gown, fabric-reinforced, XLL, 150 cm, finger loops	22
9809CE	Surgical gown, fabric-reinforced, XL, 130 cm, mask and finger loops	22
9812CEA	Surgical gown, fabric-reinforced, XLL, 150 cm, mask and finger loops	20

# OPS™ Ultimate

The intelligent gown with breathable membrane technology

- High Performance (HP) throughout the entire gown
- Made of 60-gsm trilaminate material
- Breathable, microporous fabric
- Recommended for procedures where impervious protection throughout is needed, such as long cardiological, orthopaedic and transplant surgeries

Item number	Description	Qty/case
UL89045CE	Surgical gown, breathable, impervious, XL, 130 cm	20
UL32474CE	Surgical gown, breathable, impervious, XLL, 150 cm, raglan sleeves	20
UL39079CE	Surgical gown, breathable, impervious, 2XLL, 150 cm, raglan sleeves	20



## Urology Gowns

Optimised coverage ideal for fluid-heavy, seated urology procedures

Item number	Description	Qty/case
DYNJPE2304	Surgical gown, impervious, size, with pleat, L, 114 cm	30
DYNJPE2305	Surgical gown, impervious, size, with pleat, XL, 142 cm	30





# Dual Protection



## Advanced viral and chemical protection

Enhanced clinician protection for your peace of mind is available. We offer two ranges of surgical gowns with extensive testing against chemicals and infectious agents. These solutions are classified as Medical Device class I (MDD) and Personal Protective Equipment category III (PPE) EN14605 and EN14126.

### OPS™ Essential

Essential dual-certified gowns offer MDD/PPE Category III [PB]-4 protection against chemicals and infectious agents on the reinforced area. A full list of testing results is available in the related instructions for use and technical data sheets.

Item number	Description	Qty/case
S3500P	Surgical gown, poly-reinforced, S/M, 115 cm	32
S3510P	Surgical gown, poly-reinforced, L, 120 cm	32
S3540P	Surgical gown, poly-reinforced, XL, 140 cm	28
S3513P	Surgical gown, poly-reinforced, XLL, 150 cm	24
S3581P	Surgical gown, poly-reinforced, 2XLL, 170 cm	20



### Prevention™ Plus

As the end user, you have enhanced protection against fluids and cytotoxic drugs (chemotherapy) throughout the entire gown. This is facilitated by an innovative design with a microporous breathable and impervious four-layer material—ideal for protecting clinicians when performing hyperthermic intraperitoneal chemotherapy (HIPEC) treatments.

Item number	Description	Qty/case
DYNJPE2301P	Surgical gown, impervious, L, 114 cm	24
DYNJPE2302P	Surgical gown, impervious, XL, 142 cm	24



# Direct Packaging Solution

Medline offers an efficient, economic and environmentally-friendly solution to your facility's operating theatre needs. In addition to simplifying staff operations and eliminating non-essential overwrap and hand towels, the direct packaging gown solution helps reduce your facility's waste and lower the affiliated costs.

The direct packaging solution helps minimise packaging and raw material usage, which can increase supply chain efficiency, improve waste management, and most importantly, limit negative environmental impacts. Our approach is informed by insights from our product life cycle assessment (LCA).

**OPS™ Ultragard direct gowns are available in standard and poly-reinforced ranges.**

Item number	Description	Qty/case
E3605CE	Surgical gown, S/M, 115 cm, direct	40
E3606CE	Surgical gown, L, 120 cm, direct	40
E3645CE	Surgical gown, XL, 130 cm, direct	40
E3616CE	Surgical gown, XLL, 150 cm, direct	36
E3685CE	Surgical gown, 2XLL, 160 cm, direct	36
E3600CE	Surgical gown, poly-reinforced, S/M, 115 cm, direct	36
E3610CE	Surgical gown, poly-reinforced, L, 120 cm, direct	36
E3640CE	Surgical gown, poly-reinforced, XL, 140 cm, direct	32
E3613CE	Surgical gown, poly-reinforced, XLL, 150 cm, direct	28
E3681CE	Surgical gown, poly-reinforced, 2XLL, 160 cm, direct	24



**OPS™ Essential direct gowns are available in standard, fabric-reinforced, and poly-reinforced ranges.**

Item number	Description	Qty/case
S3606CE	Surgical gown, L, 120 cm, direct	40
S3645CE	Surgical gown, XL, 130 cm, direct	40
S3616CE	Surgical gown, XLL, 150 cm, direct	36
S3685CE	Surgical gown, 2XLL, 170 cm, direct	32
S3609CE	Surgical gown, fabric-reinforced, XL, 130 cm, direct	30
FS3645CE	Surgical gown, essential, standard, XL, 140 cm, direct	36
S3612CE	Surgical gown, fabric-reinforced, XLL, 150 cm, direct	28
S3611CE	Surgical gown, fabric-reinforced, 2XLL, 170 cm, direct	24
S3610CE	Surgical gown, poly-reinforced, L, 120 cm, direct	32
S3640CE	Surgical gown, poly-reinforced, XL, 140 cm, direct	28
S3613CE	Surgical gown, poly-reinforced, XLL, 150 cm, direct	24
S3681CE	Surgical gown, poly-reinforced, 2XLL, 170 cm, direct	20



# Surgical Gown Selection Guide

Choose the gown with the right performance for every surgery, based on the type and length of the procedure, amount of fluids and required comfort.

Range		Eclipse			UltraGard		Essential			Advanced			Ultimate
Standard/Fabric-reinforced/ Poly-reinforced		S	F	P	S	P	S	F	P	S	F	P	S
Surgical Procedures	Ophthalmic	✓			✓		✓			✓			
	Tracheostomy	✓			✓		✓			✓			
	Ear surgery	✓			✓		✓			✓			
	Cochlea implants	✓			✓		✓			✓			
	Mouth	✓			✓		✓			✓			
	Tonsillectomy	✓			✓		✓			✓			
	Nose surgery	✓			✓		✓			✓			
	Varicose veins	✓	✓		✓		✓	✓		✓	✓		
	Proctology	✓	✓		✓		✓	✓		✓	✓		
	Obstetrics		✓	✓		✓		✓	✓		✓	✓	✓
	Maxillofacial	✓	✓		✓		✓	✓		✓	✓		
	Endoabdominal	✓			✓		✓			✓			
	Endocrine	✓			✓		✓			✓			
	Pacemakers	✓			✓		✓			✓			
	General orthopaedics	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Breast procedures	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Burns, plastics and grafting		✓	✓		✓		✓	✓		✓	✓	✓
	Angiography/angioplasty	✓	✓		✓		✓	✓		✓	✓	✓	
	Gynaecological	✓	✓		✓		✓	✓		✓	✓		
	General abdominal surgery		✓	✓		✓		✓	✓		✓	✓	✓
	Urogenital		✓	✓		✓		✓	✓	✓	✓	✓	✓
	Orthopaedics, hips		✓	✓		✓		✓	✓		✓	✓	✓
	Kidney, adrenals and bladder	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Orthopaedics, knees		✓	✓		✓		✓	✓		✓	✓	✓
	Thoracic		✓	✓		✓		✓	✓		✓	✓	✓
	Cardiovascular surgery		✓	✓		✓		✓	✓		✓	✓	✓
	Cranial		✓	✓		✓		✓	✓		✓	✓	✓
	Spinal procedures		✓	✓		✓		✓	✓		✓	✓	✓

**For more information** on our products, please contact  
your account manager or visit us at [www.medline.eu/uk](http://www.medline.eu/uk)



**Medline Industries Ltd.**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237

[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

**Medline Ireland Ltd.**  
1 Grant's Row  
Mount Street Lower  
Dublin 2  
Ireland  
Tel.: +353 1 691 73 06

[www.medline.eu/ie](http://www.medline.eu/ie)  
[ie-customerservice@medline.com](mailto:ie-customerservice@medline.com)


FOLLOW US  

These products are class I sterile medical devices intended to be used by healthcare professionals.  
Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML281\_EN 07/2023

 Medline International France SAS  
5 Rue Charles Lindbergh,  
44110 Châteaubriant, France

 2797 (BSI)