

# Surgical Gowns

Performance, quality and variety



**Designed, tested  
and trusted  
by operating  
theatre  
professionals**

# Medline Surgical Gowns

Performance, quality and variety

### An extensive choice

To satisfy preferences in the operating theatre, we offer five different fabric options for surgical gowns, each with a distinct style and reinforcement.

### Dependable protection




Medline's gown fabrics offer high-level protection. Each gown is latex-free and meets or exceeds the minimum performance requirements of the EN 13795 standard.\* The Essential and Advanced materials have a high hydrostatic head, and the extra alcohol repellence treatment gives wearers the protection they need and comfort they deserve. The unique Ultimate gown offers all-over impervious protection with exceptional breathability.

### High performance with high comfort

All our surgical gowns offer operating theatre staff excellent comfort and ease of use. High-tech, breathable fabrics keep the wearer comfortable and dry, even during long procedures. The gowns feature sleeves with glued seams, which contribute to greater strength and comfort.

### Easy to use

All Medline surgical gowns feature intuitive colour coding for easier gown selection in the operating theatre. Purple product labels and transfer tabs indicate standard performance, and blue is for high performance. In addition, we offer three different gown collar colours to quickly identify the gown style—standard, fabric-reinforced or poly-reinforced.



\*The EN 13795 is the European standard, specifying requirements for single-use and reusable coverings—surgical gowns and drapes used as medical devices for patients, clinical staff and equipment, intended to prevent the transmission of infective agents during invasive surgical procedures.

# OPS™ UltraGard

Effective gown, affordable price

- Lightweight SMMMS (35 gsm) material with natural, fluid-repellent protection
- Comfortable and breathable
- Gown available with either standard or poly-reinforced protection

Item number	Description	Qty/case
E3505CEA	Surgical gown, S/M, 115 cm	36
E3506CEA/E3506DCE*	Surgical gown, L, 120 cm	36/48*
E3545CEA/E3545DCE*	Surgical gown, XL, 130 cm	36/48*
FE3545CE	Surgical gown, XL, 140 cm	32
E3516CEA/E3516DCE*	Surgical gown, XLL, 150 cm	32/42*
E3585CEA/E3585DCE*	Surgical gown, 2XLL, 160 cm	20/42*
E3500CEA	Surgical gown, poly-reinforced, S/M, 115 cm	32
E3510CEA	Surgical gown, poly-reinforced, L, 120 cm	32
E3540CEA	Surgical gown, poly-reinforced, XL, 140 cm	28
E3513CE	Surgical gown, poly-reinforced, XLL, 150 cm	24
E3581CE	Surgical gown, poly-reinforced, 2XLL, 160 cm	20

\*These gowns are also available in dispenser boxes.



UltraGard gowns are also available in single-layer packaging, without overwrap and hand towels. This configuration helps improve efficiency by reducing set-up time and budgetary needs for gowns.

Item number	Description	Qty/case
E3606CE	Surgical gown, L, 120 cm, direct, ST	40
E3645CE	Surgical gown, XL, 130 cm, direct, ST	40
E3616CE	Surgical gown, XLL, 150 cm, direct, ST	36
E3640CE	Surgical gown, poly-reinforced, XL, 140 cm, direct, ST	32
E3613CE	Surgical gown, poly-reinforced, XLL, 150 cm, direct, ST	28
E3681CE	Surgical gown, poly-reinforced, 2XLL, 160 cm, direct, ST	24





# OPS™ Essential

## Enhanced protection

- Lightweight material SMMMS (35 gsm) with extra treatment for enhanced repellence
- Comfortable: offers exceptional drapeability and breathability
- Essential poly-reinforced gowns comply with EN 13795, EN 14605 and EN 14126 standards

Item number	Description	Qty/case
S3505CEC	Surgical gown, S/M, 115 cm	36
S3506CEC	Surgical gown, L, 120 cm	36
S3545CEC	Surgical gown, XL, 130 cm	36
FS3545CEA	Surgical gown, XL, 140 cm	32
S3507CEC	Surgical gown, XL, 140 cm, finger loops	32
S3516CEC	Surgical gown, XLL, 150 cm	32
S3585CEC	Surgical gown, 2XLL, 170 cm	32
S3504CEC	Surgical gown, fabric-reinforced, S/M, 115 cm	32
S3508CEC	Surgical gown, fabric-reinforced, L, 120 cm	28
S3509CEC	Surgical gown, fabric-reinforced, XL, 130 cm	28
S3512CEC	Surgical gown, fabric-reinforced, XLL, 150 cm	28
S3511CEC	Surgical gown, fabric-reinforced, 2XLL, 170 cm	24
S3500CEC/S3500P*	Surgical gown, poly-reinforced, S/M, 115 cm	32
S3510CEC/S3510P*	Surgical gown, poly-reinforced, L, 120 cm	32
S3540CEC/S3540P*	Surgical gown, poly-reinforced, XL, 140 cm	28
S3513CEB/S3513P*	Surgical gown, poly-reinforced, XLL, 150 cm	24
S3581CEB/S3581P*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20
S3740CE	Surgical gown, poly-reinforced, XL, 140 cm, finger loops	28
S3713CE	Surgical gown, poly-reinforced, XLL, 150 cm, finger loops	24
S3840CE	Surgical gown, poly-reinforced, XL, 140 cm, mask and finger loops	28
S3813CE	Surgical gown, poly-reinforced, XLL, 150 cm, mask and finger loops	24



\*Compliant with EN 14605 and EN 14126, personal protective equipment category III

### Essential gowns are also available in single-layer packaging.

Item number	Description	Qty/case
S3606CE	Surgical gown, L, 120 cm, direct	40
S3645CE	Surgical gown, XL, 130 cm, direct	40
S3616CE	Surgical gown, XLL, 150 cm, direct	36
S3685CE	Surgical gown, 2XLL, 170 cm, direct	32
S3609CE	Surgical gown, fabric-reinforced, XL, 130 cm, direct	30
S3612CE	Surgical gown, fabric-reinforced, XLL, 150 cm, direct	28
S3611CE	Surgical gown, fabric-reinforced, 2XLL, 170 cm, direct	24



# OPS™ Advanced

## No compromise between comfort and protection

- Breathable: quickly transfers heat away to keep the wearer cool and comfortable
- Higher basis weight, soft SMMMS (41 gsm) for enhanced protection and comfort
- High barrier performance
- Optional raglan cut increases roominess and provides easier movement

Item number	Description	Qty/case
OP9505CE	Surgical gown, S/M, 120 cm	28
OP9506CEA	Surgical gown, L, 130 cm	32
OP9507CEA	Surgical gown, XL, 140 cm	32
OP9516CEB	Surgical gown, XLL, 150 cm	32
OP39516CE	Surgical gown, XLL, 150 cm, raglan sleeves	24
OP9575CE	Surgical gown, XXL, 150 cm	24
OP9585CE	Surgical gown, 2XLL, 170 cm	20
OP9590CE	Surgical gown, 3XL, 135 cm	24
OP9508CEA/OP9508DCE*	Surgical gown, fabric-reinforced, L, 130 cm	24/40*
OP9809CE	Surgical gown, fabric-reinforced, XL, 130 cm, mask and finger loops	28
OP9509CEA	Surgical gown, fabric-reinforced, XL, 140 cm	28
OP9512CE/OP9512DCE*	Surgical gown, fabric-reinforced, XLL, 150 cm	24/32*
OP39512CE	Surgical gown, fabric-reinforced, XLL, 150 cm, raglan sleeves	20
OP9812CE	Surgical gown, fabric-reinforced, XLL, 150 cm, mask and finger loops	24
OP9583CE/OP9583DCE*	Surgical gown, fabric-reinforced, 2XLL, 170 cm	20/28*
OP9500CE	Surgical gown, poly-reinforced, S/M, 120 cm	28
OP9510CE/OP9510DCE*	Surgical gown, poly-reinforced, L, 130 cm	28/40*
OP9540CE/OP9540DCE*	Surgical gown, poly-reinforced, XL, 140 cm	24/36*
OP9513CE	Surgical gown, poly-reinforced, XLL, 150 cm	20
OP9581CE	Surgical gown, poly-reinforced, 2XLL, 170 cm	20



\*These gowns are also available in dispenser boxes.

# OPS™ Ultimate

The intelligent gown with breathable membrane technology

- High performance across the entire gown and fully impervious to blood-borne pathogens (AAMI\* PB70, ASTM F1671 tested, EN 13795 compliant)
- Made of 60 gsm trilaminate material
- Breathable and cool: reacts to increasing body temperature by increasing water vapour transmission
- Recommended for all procedures where breathable protection all over is needed, such as long cardiological, orthopaedic and transplant surgeries

Item number	Description	Qty/case
LW89045CEB	Breathable, impervious surgical gown, XL, 130 cm	20
LW32474CEB	Breathable, impervious surgical gown, XLL, 150 cm, raglan sleeves	20
LW39079CEB	Breathable, impervious surgical gown, 2XLL, 150 cm, raglan sleeves	20



We also offer surgical gowns featuring special fabrics, tested against chemotherapy and cytotoxic drugs.

Item number	Description	Qty/case
DYNJPE2301P	Impervious OR gown, prevention plus, L, 114 cm	24
DYNJPE2302P	Impervious OR gown, prevention plus, XL, 142 cm	24

Prevention Plus gowns comply with EN 14605 and PPE Regulation (EU) 2016/425.


\*AAMI: Association for the Advancement of Medical Instrumentation. The AAMI PB70 is the US reference standard, establishing a system of classification for protective apparel and drapes, based on their liquid barrier performance. It specifies related labelling requirements and standardised testing methods.

# Eclipse®

A soft and comfortable gown

- Outstanding comfort and breathability
- Spunlace (68 gsm) with a textile-like appearance
- Suitable protection

Item number	Description	Qty/case
9505CEA	Surgical gown, S/M, 115 cm	28
9506CEA/9506DCE*	Surgical gown, L, 120 cm	28/40*
9415CEA	Surgical gown, L, 120 cm, raglan sleeves	26
9507CEA/9507DCE*	Surgical gown, XL, 130 cm	24/40*
9707CE	Surgical gown, XL, 130 cm, finger loops	24
9516CEA/9516DCE*	Surgical gown, XLL, 150 cm	24/36*
9716CE	Surgical gown, XLL, 150 cm, finger loops	24
9575CEA	Surgical gown, XXL, 130 cm	24
9585CE/9585DCE*	Surgical gown, 2XLL, 170 cm	24/32*
9508CEA	Surgical gown, fabric-reinforced, L, 120 cm	24
9509CEA/9509DCE*	Surgical gown, fabric-reinforced, XL, 130 cm	22/32*
9539CEA	Surgical gown, fabric-reinforced, XL, 130 cm	7
9709CE	Surgical gown, fabric-reinforced, XL, 130 cm, finger loops	22
9809CE	Surgical gown, fabric-reinforced, XL, 130 cm, mask and finger loops	22
9571CEA	Surgical gown, fabric-reinforced, XXL, 130 cm	20
9512CEA/9512DCE	Surgical gown, fabric-reinforced, XLL, 150 cm	22/32*
9451CEA	Surgical gown, fabric-reinforced, XLL, 150 cm	11
9712CE	Surgical gown, fabric-reinforced, XLL, 150 cm, finger loops	22
9812CEA	Surgical gown, fabric-reinforced, XLL, 150 cm, mask and finger loops	20
9581CEA/9581DCE*	Surgical gown, fabric-reinforced, 2XLL, 170 cm	19/28*
9000CEA	Surgical gown, poly-reinforced, S/M, 115 cm	28
9010CEA/9010DCE*	Surgical gown, poly-reinforced, L, 120 cm	28/32*
9040CEA	Surgical gown, poly-reinforced, XL, 130 cm	24
9513CEA	Surgical gown, poly-reinforced, XLL, 150 cm	24
9580CE/9580DCE*	Surgical gown, poly-reinforced, 2XLL, 170 cm	20/28*



\*This item is also available in dispenser boxes.

# Surgical Gown Selection Guide

Choose the gown with the right performance for every surgery, based on the type and length of the procedure, amount of fluids and required comfort.

Range		Eclipse			UltraGard		Essential			Advanced			Ultimate
		S	F	P	S	P	S	F	P	S	F	P	S
DRG Surgical Procedures	Ophthalmic	✓			✓		✓			✓			
	Tracheostomy	✓			✓		✓			✓			
	Ear surgery	✓			✓		✓			✓			
	Cochlea implants	✓			✓		✓			✓			
	Mouth	✓			✓		✓			✓			
	Tonsillectomy	✓			✓		✓			✓			
	Nose surgery	✓			✓		✓			✓			
	Varicose veins	✓	✓		✓		✓	✓		✓	✓		
	Proctology	✓	✓		✓		✓	✓		✓	✓		
	Obstetrics		✓	✓		✓		✓	✓		✓	✓	✓
	Maxillofacial	✓	✓		✓		✓	✓		✓	✓		
	Endoabdominal	✓			✓		✓			✓			
	Endocrine	✓			✓		✓			✓			
	Pacemakers	✓			✓		✓			✓			
	General orthopaedics	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Breast procedures	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Burns, plastics and grafting		✓	✓		✓		✓	✓		✓	✓	✓
	Angiography/angioplasty	✓	✓		✓		✓	✓		✓	✓	✓	
	Gynaecological	✓	✓		✓		✓	✓		✓	✓		
	General abdominal sugery		✓	✓		✓		✓	✓		✓	✓	✓
	Urogenital		✓	✓		✓		✓	✓	✓	✓	✓	✓
	Orthopaedics, hips		✓	✓		✓		✓	✓		✓	✓	✓
	Kidney, adrenals and bladder	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Orthopaedics, knees		✓	✓		✓		✓	✓		✓	✓	✓
	Thoracic		✓	✓		✓		✓	✓		✓	✓	✓
	Cardiovascular surgery		✓	✓		✓		✓	✓		✓	✓	✓
	Cranial		✓	✓		✓		✓	✓		✓	✓	✓
	Spinal procedures		✓	✓		✓		✓	✓		✓	✓	✓

**For more information** on our products, please contact your account manager or visit us at [uk.medline.eu](http://uk.medline.eu)



**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
Fax: +44 844 334 5238  
[uk.medline.eu](http://uk.medline.eu) | [uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

FOLLOW US  

These products are class I sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

 Medline International France SAS  
5, rue Charles Lindbergh  
44110 Châteaubriant, France

 2797 (BSI)

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2020 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML281\_EN 09/2020