

Non-Touch Infrared Forehead Thermometer

TH5001N



**Quick, safe
and accurate
temperature
monitoring**



Non-Touch Infrared Forehead Thermometer

TH5001N

Avoid cross-contamination
via non-touch
infrared
technology

Benefits and features

- Easy-to-use one-touch operation
- Facilitates the measurement of forehead temperature from 2-cm distance
- Provides results in less than 2 seconds
- Allows for switching between temperature units (°C/°F)
- Offers flexibility of use in a range of settings (clinical or otherwise, including home)
- Gives the end user the capacity to also measure the temperature of liquids, as well as ambient temperature
- Runs on 9-volt battery (included)
- Comes with 2-year warranty



Ordering information

| Item number | Description | Packaging |
|-------------|--|-----------|
| TH5001N | Digital non-touch forehead thermometer | 1 per box |

For more information on our products, please contact your Medline account manager or visit us at: uk.medline.eu



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
uk.medline.eu | uk-customerservice@medline.com

FOLLOW US  

This product is a class IIa medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2021 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML1102_EN 01/2021

Dongshin Technology Co., Ltd.
Floor 1-2, W3 Building, Foshan Xusheng
Industrial Estate Xilixiedamang,
518108 Nanshan District, Shenzhen, China



Shanghai International Holding Corp GmbH (Europe)
Eiffelstraße 80
20537 Hamburg, Germany

