

Trumpet valve for suction and irrigation


For laparoscopic procedures



**Easily
control
the flow.**

Trumpet valve

For suction and irrigation



Product name

Code	Description	Cannula diameter	Cannula length	Packaging
OR1200S	Suction/irrigation kit, single spike	5 mm	330 mm	12 ea/cs
OR1201D	Suction/irrigation kit, double spike	5 mm	330 mm	12 ea/cs
ONG299601CST12	Suction/irrigation kit, single spike	10 mm	330 mm	5 ea/cs

Product information

- A non-conductive polycarbonate* cannula with a satin surface to prevent light reflection (available with a 5 or 10 mm diameter)
- An ergonomic handle that can be used ambidextrously with 2 buttons: the red one for suction and the blue one for irrigation
- Single or double spike to connect 1 or 2 saline bags



KEY FEATURES

- Non-conductive cannula*
- Single-use
- DEHP-free
- Latex-free



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk | uk-customerservice@medline.com

FOLLOW US  

These products are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

*OR1200S and OR1201D:  Luxus Lebenswelt GmbH Köchstr. 1 49877 Willich, Germany  0197 (TÜV) ONG299601CST12:  Rocamed France Z.I. de Signes Allée de Stockholm 83870 Signes, France  0459 (LNE-GMED)

We reserve the right to correct any errors that may occur within this brochure.

© 2019 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML650-EN 10/2019