

# Ward Supplies

Product overview catalogue



**Patients first,  
always ...**

**Medline offers a range of products, designed to help enhance patient well-being and facilitate optimal comfort. This brochure is, therefore, a snapshot of our dedication to providing healthcare professionals with the tools they need to deliver ideal care.**



# Table of contents

<b>Emesis Bag</b>	<b>4</b>
Emesis Bag	4
Emesis Dispenser	5
<b>Lube Jelly</b>	<b>6</b>
Ultrasound Gel	6
<b>Drug Management Devices</b>	<b>7</b>
Medication Dispenser	7
<b>Patient Grooming</b>	<b>8</b>
Surgical Razor	8
Denture Cup	8
Super-Soft Toothbrush	8
AloeTouch Wet Wipes	8
Surgical Clipper	9
<b>Cotton Applicator</b>	<b>10</b>
Cotton Applicator	10
<b>Diagnostic Instruments</b>	<b>11</b>
Safety Lancets	11
Tourniquets	11
ECG Electrode Spray	12
Stethoscopes	12
<b>Hot and Cold Therapy</b>	<b>13</b>
Hot/Cold Pack	13
Instant Cold Pack	13
<b>Physical Therapy</b>	<b>14</b>
Drinking Beakers with Lid	14
<b>Miscellaneous</b>	<b>15</b>
Glycerine Swabsticks	15
Feeding Bib	15
Carrying Patient Bag	15

# Emesis Bag

Easy to use and dispose of

Medline's Emesis Bag is designed to provide a safe, disposable and easy solution for the collection, containment and disposal of expelled vomitus, by way of its twist-and-seal closure. It is an ideal solution in countless departments across the hospital and continuum of care.



**1 Efficient collection**  
Designed to prevent spills and eliminate caregiver and patient contact with emesis

**2 Safe containment**  
Designed to isolate contaminants and help prevent cross-contamination on floors and surfaces

**3 Easy disposal**  
Simply twist the bag closed, hook into the notch and dispose of the bag

**Our emesis bag is designed to offers:**

- Lower risk of splashback and mess versus a basin; allows for quicker clean-up
- Lower risk of cross-contamination, fewer linen changes and floor cleanings due to an isolating plastic ring
- Improved privacy and discreetness for the patient
- More independence, facilitating self-help.



Item number	Description	Packaging
EMB01	Emesis bag, blue	144 pcs/cs
EMB02	Emesis bag, clear	144 pcs/cs

## Emesis Dispenser

Convenient and easy to install

Medline's Emesis Bag Dispenser is a wall dispenser designed to hold up to 24 units of the Emesis Bag. Additionally, it can be placed conveniently in any location and easily installed.



- +** **Removable cap**  
Designed to help prevent contamination
- +** **Compact**  
16.5 cm long
- +** **Light-weight**  
0.4kg
- +** **High capacity**  
Designed to hold 24 emesis bags

Item number	Description	Packaging
DISP01	Emesis bags dispenser	6 pcs/cs

# Ultrasound Gel\*

Safe and easy to remove

Medline's ultrasound gel is designed to be water-soluble. It is easy to wipe off, without drying, flaking or leaving residue behind. In addition, it is hypoallergenic, making it non-irritating and safe for use on sensitive skin.

## Non-Sterile Ultrasound Gel

For examining areas with intact skin

### Non-sterile ultrasound gel in bottle

Item number	Description	Packaging
782010	SONOSID® ultrasound gel, 250 ml, bottle	16 pcs/cs
782020	SONOSID ultrasound gel, 500 ml, bottle	25 pcs/cs
782030	SONOSID ultrasound gel, 1 litre, bottle	15 pcs/cs



### Non-sterile ultrasound gel in canister

Item number	Description	Packaging
782040	SONOSID ultrasound gel, 5 litres, canister	4 pcs/cs



### Non-sterile ultrasound gel in cubitainer

Item number	Description	Packaging
782050	SONOSID ultrasound gel, 5 litres, cubitainer	1 pc/cs



The 5-litre cubitainer comes with nozzles and refill bottles

\*Available in Germany only.

# Medication Dispenser

For efficient tracking of daily medication intake

Medline's one-day medication dispenser is designed to help patients take their medication at different times of the day. It also helps healthcare providers track their patients' daily dosage.

Item number	Description	Packaging
733010	Medication dispenser, small compartments, EN/DE	200 pcs/cs
733020	Medication dispenser, large compartments, EN/DE	200 pcs/cs
733010FRIT	Medication dispenser, small compartments, FR/IT	200 pcs/cs



**2 compartment sizes available**  
Designed to meet a range of dosage needs

**Time indicators**  
Designed to facilitate time management via words and icons



**Braille-enabled**  
Designed to be inclusive for use by visually-impaired or blind patients



**Plastic cover/slider**  
Designed to help ensure that every pill stays in its designated compartment



## Surgical Razor

The disposable Surgical Razor is designed for use in preoperative shaving periods, when stainless steel blades, with ergonomically angled heads, are required.

Item number	Description	Packaging
726070	Surgical single-blade razor, blue	1,000 pcs/cs
726075	Surgical double-blade razor, green	1,000 pcs/cs



## Denture Cup

The one-piece Denture Cup with snap-on lid is designed for storing dentures and braces.

Item number	Description	Packaging
T12-900-00-901	Denture cup, 200 ml, green	25 pcs/cs



## Super-Soft Toothbrush

The Super-Soft Toothbrush is designed to be well-suited to use by/on patients who have sensitive teeth or sensitive gums, thanks to its super-soft nylon bristles.

Item number	Description	Packaging
TB001	Pediatric Toothbrush	500 pcs/cs
MDS096082	Adult Toothbrush	72 pcs/cs



## AloeTouch Wet Wipes

Pre-moistened wet wipes that are ideal for continence care and personal clean-up. These gender-neutral, lightly scented wipes are dermatologically tested and contain a pH-balanced formula.

Item number	Description	Packaging
MSCE263654	Spunlace wipes; scented, 20 x 30 cm	48 wipes/pcs; 12 pcs/cs





# Surgical Clipper

Medline's Surgical Clipper is designed to allow for optimal preoperative hair removal. The clipper is also designed to be easy to handle and facilitate one-handed manoeuvrability and better control, thanks to the handle grips on its sides.

Item number	Description	Colour	Packaging
CCSB4	Close shave blade	Black	25/cs
CUB1	Universal blade	Grey	50/cs
CNB2	Sensitive area blade	Red	50/cs
CCB3	Neuro blade	Blue	50/cs
DYND70800EU	Clipper handle and charger, EU plug	White	1
DYND70800UK	Clipper handle and charger, UK plug	White	1
DYND70801EU	Clipper handle only	White	1
DYND70802EU	Charger only, EU plug	White	1
DYND70802UK	Charger only, UK plug	White	1

Close Shave Blade



729010

Universal Blade



729010

Sensitive Blade



729010

Neuro Blade



729010

# Cotton Applicator

Cotton applicators are designed to help clinicians apply substances, such as medication, to the ear/nose/throat area, or for wound care treatments.

Item number	Description	Packaging
830000A	Cotton applicator; small, non-sterile	100/pk, 10/bx, 20/cs
830010A	Cotton applicator; small, sterile	2/pk, 100/pch, 15/cs
831510A	Cotton applicator; big, sterile	2/pk, 50/pch, 15/cs
831578A	Cotton applicator; big, non-sterile	50/pk, 20/cs



## Safety Lancets

Medline's Safety Lancet Needle are designed to enable finger pricking, which is necessary when extracting microsamples of capillary blood from adults.

Item number	Description	Colour	Colour	Colour	Packaging
SLN2118	Normal flow	Needle; 21 G	1.8 mm	Green	200/cs
SLN1818	High flow	Needle; 18 G	1.8 mm	Orange	200/cs
SLN2816	Micro flow	Needle; 28 G	1.6 mm	Blue	200/cs
SLBLADE12	Pediatric	Blade	1.2 mm	Purple	200/cs

Normal flow



SLN1818

High flow



SLN2816

Micro flow



SLM2118

Pediatric



SLBLADE12

## Tourniquets

Our single-use tourniquet solution is designed for use with patients requiring venous blood sampling and intravenous cannulation.

Item number	Description	Packaging
DYNDE75025	Disposable rolled tourniquets 2.5 cm x 53.3 cm x 0.05 cm	25 pcs/box 1000 pcs/cs
DYNDE75020A	Disposable rolled tourniquets 2.5 cm x 45.7 cm x 0.05 cm	25 pcs/box 1000 pcs/cs



## ECG Electrode Spray\*

The ECG Electrode Spray is designed to enhance conductivity and improve precision and reliability during electrocardiogram assessments.



Item number	Description	Packaging
780820	ECG/EEG electrodes spray	1 e/a box; 20 c/s

## Stethoscopes\*

Color-coded stethoscopes are designed to offer improved acoustics and clinical precision.

Item number	Description	Packaging
729010	Stethoscope, blue	1 e/a box; 50 c/s
729020	Stethoscope, green	1 e/a box; 100 c/s
729030	Stethoscope, red	1 e/a box; 100 c/s
729040	Stethoscope, black	1 e/a box; 100 c/s



729010



729020



729030



729040

\*Available in Germany only.

## Hot/Cold Pack

Hot and cold packs, with soft-touch non-woven (fabric-like) lining, are designed to help prevent condensation and to be placed directly on patients' skin.

Item number	Description	Size	Packaging
NON4410E	Ice bag, clamp-close, 4 ties	35.7 x 13.3 cm	25 ea/bx; 2 bx/cs
NON4420E	Ice bag, clamp close, 4 ties	41.3 x 17.9 cm	25 ea/bx; 2 bx/cs
NON4430E	Ice bag, bilateral	35.5 x 13 cm	12 ea/bx; 2 bx/cs
NON4440E	Ice bag, *EENT style	28.5 x 19.9 cm	15 ea/bx; 2 bx/cs
NON4460E	Ice bag, clamp close, hook	41.3 x 17.9 cm	15 ea/bx; 2 bx/cs
NON4480E	Ice bag, perineal, flexible	35.7 x 11 cm	50 ea/bx; 2 bx/cs



## Instant Cold Pack

Instant cold pack, with soft-touch non-woven (fabric-like) lining, is designed to help prevent condensation and to be placed directly on patients' skin.

Item number	Description	Size	Target temperature	Packaging
CPACKP	Delux Perineal Instant Cold Pack	11.4 x 36.2 cm	-	24pcs/cs
CPACKS	Delux Instant Cold Pack	12.7 x 19 cm	0° C/32° F	24pcs/cs
CPACKL	Delux Instant Cold Pack	16 x 24 cm	0° C/32° F	24pcs/cs



## Drinking Beakers with Lid\*

The Drinking Beaker with lid is designed to provide an easy solution to giving water, or other liquids, to bedridden patients. Depending on the capability of the patient, in terms of swallowing large or small quantities, this product is available with two lid sizes: narrow and wide.

Item number	Description	Packaging
713030	Drinking beaker, 250 ml, narrow spout	50 pcs/cs
713035	Drinking beaker, 250 ml, wide spout	50 pcs/cs
312000	Drinking cup, 220 ml	20 pcs/cs
312001	Narrow spout	20 pcs/cs
312002	Wide spout	20 pcs/cs

\*Available in Germany only.



713030



713035



312000



312001



312002

# Glycerine Swabsticks

Lemon-flavoured cotton buds are designed to refresh the mouth and throat.

Item number	Description	Packaging
727060	Lemon glycerine swabsticks, 15 cm	3 pcs/pk 30 pcs/cs



# Feeding Bib

The tie-on disposable bib is designed to meet a variety of patient needs, including post-op recovery and geriatric care.

Item number	Description	Packaging
713010	Disposable feeding bib, 65 x 37 cm	1,000 pcs/cs



# Carrying Patient Bag

The patient bags provide a convenient storage solution where patients can keep their possessions during surgical ambulatory care, or while in the emergency department.

Item number	Description	Packaging
714010-PE	Carrying patient bag, *PE, 58 x 58 cm	250 pcs/cs
PBAG5060	Patient belonging bag, 50 x 60 cm	250 pcs/cs



**With our wide-ranging ward supplies, putting your patients first is that much easier.**

**For more information on our products, please contact your account manager or visit us at: [www.medline.eu/uk](http://www.medline.eu/uk)**



**Medline Industries Ltd.**  
 3rd Floor  
 Quayside Wilderspool Business Park  
 Greenalls Avenue  
 Warrington WA4 6HL  
 United Kingdom  
 Tel.: +44 844 334 5237  
[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

**Medline Ireland Ltd.**  
 1 Grant's Row  
 Mount Street Lower  
 Dublin 2  
 Ireland  
 Tel.: +353 1 691 73 06

[www.medline.eu/ie](http://www.medline.eu/ie)  
[ie-customerservice@medline.com](mailto:ie-customerservice@medline.com)

FOLLOW US

- Emesis Bags are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France
- SONOSID ultrasound gel is a class I non-sterile medical device intended to be used by healthcare professionals. GELL'D® GmbH Gettechnik  
Harmate 28  
48683 Ahaus, Germany
- Surgical clippers are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France
- Cotton applicators are class I non-sterile medical devices intended to be used by healthcare professionals. CMC Medical Devices & Drugs S.L.  
C/ Horacio Lengo No18  
CP 29006, Málaga, Spain
- Safety lancets are class IIa sterile medical devices intended to be used by healthcare professionals. Shanghai International Holding Corp.GmbH (Europe)  
Eiffelstraße 60  
20537 Hamburg, Germany 0123 (TÜV Rheinland)
- Tourniquets are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France
- ECG Electrode spray is a class I non-sterile medical device intended to be used by healthcare professionals. SAV Liquid Production GmbH  
Hochriesstr. 2  
D-83126 Flintsbach a. Inn, Germany
- Stethoscopes are class I non-sterile medical devices intended to be used by healthcare professionals. ZOUSTECH S. L. Pso. Castellana  
141, Pso. 19  
28046 Madrid, Spain
- Hot/Cold Packs are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France
- Instant Cold Packs are class IIa non-sterile medical devices intended to be used by healthcare professionals. LOTUS NL B.V.  
Koningin Julianaplein 10  
1e Verd., 2595AA, The Hague, Netherlands 1639 (SGS)
- Drinking cup, ref. 312000, is a class I non-sterile medical device intended to be used by healthcare professionals. Medical Device Management Ltd.  
Block B, The Crescent Building, Northwood, Santry,  
Dublin 9, D05 C6X8, Ireland.
- Glycerine swabstick, ref. 727060, is a class I non-sterile medical device intended to be used by healthcare professionals. WellKang Ltd/Enterprise Hub,  
NW Business Complex, 1 Beraghmore Road,  
Derry, BT48 8SE, Northern Ireland.