

INSERT/BAG LEVEL

FRONT SIDE

3.94 in  
10 cm

BACK SIDE

3.94 in  
10 cm

5.91 in  
15 cm



REF

OR169

x10

en

de

fr

nl

it

es

pt

sv

fi

da

no

cs

pl

sl

SPECIMEN TRAP/ SPECIMEN SOCK

GEWEBEKOLLEKTOR WEIT

DISPOSITIF DE PRELEVEMENT, TYPE CHAUSSETTE

WEEFSELOPVANGZAK

SACCHETTO RACCOLTA PARTICELLARE

ACCESORIO DE TOMA DE MUESTRAS

DISPOSITIVO DE RECOLHA DE AMOSTRA

PREPARATUPPSAMLINGSBURK

NÄYTTEENKERÄYSSIIVILÄ

BEHOLDER TIL PRÆPERATOPSAMLING

PREPARATOPPSAMLINGSKANNE

ZÁCHYTNÁ NÁDOBA NA VZORKY / NÁVLEK NA VZORKY

SYFON NA PRÓBKII/RĘKAW NA PRÓBKII

LOVILNIK/NOGAVICA ZA VZORCE



Latex free, Latexfrei, Sans latex, Latexvrij, Senza lattice,  
Sin látex, Sem látex, Latexfri, Lateksiton, Latexfri,  
Lateksfri,Bez Latexu, Bezlateksowe, Brez lateksa

Medical Device - Medizinprodukt - Dispositif médical - Medical device - Dispositivo Medico -  
Producto Sanitario - Dispositivo médico - Medicinteknisk Produkt - Lääkinnällinen laite -  
Medicinsk udstyr - Medisinsk utstyr - Zdravotnickým prostředkem - Wyrób medyczny -  
Medicinski pripomoček

Manufactured in China for:  
 Medline Industries, Inc., Three Lakes Drive, Northfield, IL 60093 USA



Medline International France SAS,  
5 rue Charles Lindbergh,  
44110 Châteaubriant, France

RG19SJJP

GTIN (01): 10888277433991



(01)10888277433991


 Date of Mfg: YYYY-MM-DD


 Exp (17): YYYY-MM-28

 Lot (10): XXXYYWWN




(17)YYMM28(10)XXXYYWWN





CMYK



PMS  
COLORS  
USED

2935 C

2935 C  
60% opacity

White

**Note for the printing :**  
UDI area: generate and apply at manufacture  
Date of manufacture is to be written as follow: YYYY-MM-DD  
Expiry date is to be written as follow: YYYY-MM-28  
Caution: The month should be one less than the month from the date of manufacture  
Example: MFG date: 2018-06-04 - EXP date: 2020-05-28  
LOT number format: If possible, please follow below rule.  
XXXYYWWN: Vendor code (XXX) - Year (YY) - W (Week) - WN (for week number)

8.00 in  
20.32 cm

4.00 in  
10.16 cm



REF OR169  
x100

- en** SPECIMEN TRAP/ SPECIMEN SOCK  
**de** GEWEBEKOLLEKTOR WEIT  
**fr** DISPOSITIF DE PRELEVEMENT,  
TYPE CHAUSSETTE  
**nl** WEEFSELOPVANGZAK  
**it** SACCHETTO RACCOLTA PARTICELLARE  
**es** ACCESORIO DE TOMA DE MUESTRAS  
**pt** DISPOSITIVO DE RECOLHA DE AMOSTRA



Latex free, Latexfrei, Sans latex, Latexvrij, Senza lattice,  
Sin látex, Sem látex, Latexfri, Lateksiton, Latexfri,  
Lateksfri, Bez Latexu, Bezlateksowe, Brez lateksa


Medical Device - Medizinprodukt - Dispositif médical - Medical device - Dispositivo Medico -  
Producto Sanitario - Dispositivo médico - Medicinteknisk Produkt - Lääkinnällinen laite -  
Medicinsk udstyr - Medicinsk utstyr - Zdravotnickým prostředkem - Wyrób medyczny -  
Medicinski pripomoček

Manufactured in China for:

**Medline Industries, Inc.**, Three Lakes Drive, Northfield, IL 60093 USA




Medline International France SAS,  
5 rue Charles Lindbergh,  
44110 Châteaubriant, France RG19SJP



REF OR169  
x100

- sv** PREPARATUPPSAMLINGSBURK  
**fi** NÄYTTEENKERÄYSSIIVILÄ  
**da** BEHOLDER TIL PRÆPERATOPSAMLING  
**no** PREPARATUPPSAMLINGSKANNE  
**cs** ZÁCHYTNÁ NÁDOBA NA VZORKY /  
NÁVLEK NA VZORKY  
**pl** SYFON NA PRÓBKII/RĘKAW NA PRÓBKII  
**sl** LOVILNIK/NOGAVICA ZA VZORCE



GTIN (01): 40888277433992



(01)40888277433992



Date of Mfg: YYYY-MM-DD



Exp (17): YYYY-MM-DD



Lot (10): XXXYYWWN

PO:

Gross weight:

Net weight:

PMS COLORS USED

CMYK



Black

#### Note for the printing :

UDI area: generate and apply at manufacture

Date of manufacture is to be written as follow: YYYY-MM-DD

Expiry date is to be written as follow: YYYY-MM-28

Caution: The month should be one less than the month from the date of manufacture

Example: MFG date: 2018-06-04 - EXP date: 2020-05-28

LOT number format: If possible, please follow below rule.

XXXYYWWN: Vendor code (XXX) - Year (YY) - W (Week) - WN (for week number)

Gross/Net weight and PO info are to be added somewhere on the case