



www.medline.eu/uk



Technical Data Sheet

SensiCare[®] Ice

Non sterile Nitrile – Powder Free Glove

Reference & Size

486800	486801	486802	486803	486804
XS	S	M	L	XL

Primary Material

Latex Free
Powder-Free Nitrile

Color

Blue

Surface

Finger Tip Textured

Cuff Design

Beaded cuff

Chemical additive (accelerators)

Carbamate

Powder content (per EN455-3)

Powder Content:0.58 milligrams per glove

Thickness (based on average)

Finger Tip	0.09 mm
Palm	0.07 mm
Cuff	0.05 mm

Packaging

200 gloves per box, 10 boxes per case (XL – 180 gloves per box)
Box dimensions :245mm L x 125mm W x 74mm H
Case dimensions :382mm L x 257mm W x 255mm H

Features:

- MicroStop Film
- Flat Layer Packaging
- Color coded sizing



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-48680X -
EN02

01.INT.03.30

Length & Width Per EN455-2

486800	486801	486802	486803	486804
78 mm	87 mm	98 mm	108 mm	117 mm
246 mm	250 mm	250 mm	243 mm	250 mm

Force @ Break Before Challenge Testing (per EN455-2 \geq 6 N)

6.9 N

Force @ Break After Challenge Testing (per EN455-2 \geq 6 N, 14 days 70°C in an oven)

6.4 N

Freedom from Holes (per EN 455-1)

Pass 1.0 AQL (standard requests compliance with AQL1.5)

Viral Penetration

Tested and passed, in accordance with EN374-5

Chemical Resistance

The resistance to chemicals has been assessed in accordance with EN 16523-1 and chemo drugs under ASTM D 6978.

Results and recommendations for use with chemicals and chemo drugs can be obtained on request

Expiration Date

35 Months from Date of Manufacture

Manufacture and Expiration Dates are printed on packaging (YYYY-MM format)

Regulations and Quality Standards

Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485

PPE Certification

Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5

Food Contact

Tested according to European Regulation (EU) No 10/2011

Storage Recommendations

Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.

Country of Origin

Malaysia

Legal Manufacturer

Medline International France SAS

5 rue Charles Lindbergh, 44110 Châteaubriant, France

